

BactoReal[®] Typing Kit *E. coli* v1.1

Manual



For veterinary use only



DVET001



2 x 96 wells



ingenetix GmbH
Arsenalstrasse 11
1030 Vienna, Austria
T +43 (0)1 36 198 0 198
F +43 (0)1 36 198 0 199
office@ingenetix.com
www.ingenetix.com

Index

1. Product description.....	3
2. Pathogen information	3
3. Principle of real-time PCR.....	3
4. General Precautions.....	4
5. Contents	4
6. Additionally required materials and devices.....	5
7. Preparation of real-time PCR	5
7.1. Pipetting scheme	5
7.2. Programming of the real-time PCR instrument.....	5
8. Interpretation of Melting Curve data.....	6
Tube 1&4: F4, F18	7
Tube 2&11: F5, fimA	7
Tube 3&5: F6, F41	7
Tube 6&12: STa, AIDA.....	7
Tube 7&8: STb, Stx2e.....	7
Tube 9&10: LT, EAST	7
Tube 13&14: pAA, escV	7
Tube 15: cnf1	8
Tube 16: iucD	8
9. Annex – symbols.....	8

1. Product description

BactoReal® Typing Kit *E. coli* contains PCR primer sets for the amplification and detection of porcine *E. coli* virulence-associated genes. The method is based on an intercalating dye using real-time PCR. Amplification curve and gene specific melting curve at 515 nm indicate the amplification of *E. coli* specific genes. Screening for porcine *E. coli* adhesion and toxin genes requires picking of several bacterial colonies per plate. With one BactoReal® Typing Kit *E. coli* kit 20 colonies (2 x 96 wells) can be typed.

Mix	primer sets	fimbrial adhesins
1	F4	fimbria F4 (K88)
2	F5	fimbria F5 (K99)
3	F6	fimbria F6 (987P)
4	F18	fimbria F18
5	F41	fimbria F41 (F7)
		enterotoxins
6	STa1	heat-stable enterotoxin (estA)
7	STb2	heat-stable enterotoxin (estB)
8	Stx2e	shiga toxin variant, subtype Stx2e
9	LT	heat-labile enterotoxin (eltB)
10	EAST1	enteroaggregative heat-stable toxin (astA)
		other adhesins
11	fimA	type 1 fimbriae
12	AIDA-I	autotransporter adhesin
13	pAA	porcine attaching and effacing-associated adhesin
14	escV	translocator, type III secretion system, located in pathogenicity island LEE
15	cnf1	cytotoxic necrotizing factor 1
16	iucD	aerobactin
17	papC	pyelonephritis-associated pilus (gene of the P fimbriae for EAEC)
18	pic	serine protease autotransporter

2. Pathogen information

Most types of diarrhea are caused by strains of enterotoxigenic *E. coli* (ETEC). The main virulence factors associated with ETEC in diarrhoea are enterotoxins and fimbrial adhesins. The fimbrial types F4 (K88) and F18 are commonly found in pathogenic *E. coli* isolated from weaned pigs. F18 is most commonly found in Shiga toxin-producing *E. coli* (STEC) and the Stx2e variant has been associated with oedema disease in pigs. The fimbrial types F5 (K99), F6 (987P) and F41 are more frequently associated with *E. coli* causing neonatal diarrhea. Another virulence factor of ETEC is EAST1, which has also been detected in *E. coli* of different pathogenic strains, such as ETEC, enteropathogenic *E. coli* (EPEC) and STEC from humans and animals.

References:

Wilson RA, Francis DH. Fimbriae and enterotoxins associated with *E. coli* serogroups isolated from clinical cases of porcine colibacillosis. *Am J Vet Res.* 1986;47:213–217

3. Principle of real-time PCR

A specific DNA sequence of the pathogen genome is amplified by a set of two PCR primers that flank the target region and the generated PCR-product is detected by an intercalating nucleic acid dye with subsequent melt curve analysis.

4. General Precautions

The user should always pay attention to the following:

- Always include a negative control per PCR-run (water instead of sample).
- Be careful when handling the positive control.
- Store and extract positive material (specimens, controls and amplicons) separately from all other reagents and add it to the reaction mix in a spatially separated workspace.
- Periodically decontaminate benches and devices.
- Use sterile pipette tips with filters.
- Thaw all components thoroughly at room temperature before starting an assay. When thawed, mix the components and centrifuge briefly. The 5x HR Master Mix can be stored up to 1 month at room temperature or at 4°C, routine storage is at -20°C.

5. Contents

<i>E. coli</i> Assay Mix (green cap)	Content Primer for detection of	Amount
Tube 1 and Tube 4	F4 and F18	á 1 x 50 µl
Tube 2 and Tube 11	F5 and fimA	á 1 x 50 µl
Tube 3 and Tube 5	F6 and F41	á 1 x 50 µl
Tube 6 and Tube 12	STa and AIDA	á 1 x 50 µl
Tube 7 and Tube 8	STb and Stx2e	á 1 x 50 µl
Tube 9 and Tube 10	LT and EAST	á 1 x 50 µl
Tube 13	pAA	1 x 50 µl
Tube 14	escV	1 x 50 µl
Tube 15	iucD	1 x 50 µl
Tube 16	cnf1	1 x 50 µl
Tube 17	papC	1 x 50 µl
Tube 18	pic	1 x 50 µl

Labelling (white cap)	Content	Amount
5x HR Master Mix	Master Mix with internal reference (ROX)	2 x 400 µl

Labelling (red cap)	Content	Amount
<i>E. coli</i> Positive Control	Control-DNA	1 x 50 µl

The components of BactoReal® Typing *E. coli* should be stored at 4°C to -25°C and are stable until the expiry date stated on the label.

6. Additionally required materials and devices

- Reagents and devices for DNA-extraction
- PCR-grade water
- Disposable powder-free gloves
- Sterile pipette tips with filters
- Real-time PCR instrument
- Appropriate 96 well reaction plates or reaction tubes with corresponding (optical) closing material.

For further information a group-specific typing of *E. coli* can be performed with Ingenetix's BactoReal® Kit *E. coli* Typing which includes a series of different real-time PCR kits for the detection of stx1 and stx2 (BactoReal® Kit *E. coli* Typing STEC), and of eae (intimin) or ehxA (BactoReal® Kit *E. coli* Typing EHEC, EPEC).

7. Preparation of real-time PCR

Please make sure that at least one negative control (water), as well as one positive control (red cap) are included per PCR run.

7.1. Pipetting scheme

Sample: Pick several *E. coli* colonies and isolate the DNA (e.g. with InstaGene™ Matrix, Bio-Rad), the extracted DNA should be used in a dilution (e.g. 1:1000), so that Ct / Cp values of amplification are around 25 cycles.

Positive Control: As positive control use 1 µl of the *E. coli* Positive Control.

Prepare **one** master mix for **each** *E. coli* Assay Mix

or use them in **multiplex**: *E. coli* Assay Mix Tube 1 with 4, Tube 2 with 11,
 Tube 3 with 5, Tube 6 with 12,
 Tube 7 with 8, Tube 9 with 10, optional Tube 13 with 14:

		Per sample and assay
Preparation of master mix (mix well)	H ₂ O	14.0 µl
	5x HR Master Mix	4.0 µl
	<i>E. coli</i> Assay Mix	1.0 µl +/- 1.0 µl*
	Total volume	19.0 µl
Preparation of PCR	Master mix	19.0 µl
	Sample	1.0 µl
	Total volume	20.0 µl

*You don't have to calculate the master mix for the additional µl.

When eight PCRs per colony are on a 96-well plate, prepare the respective master mixes in the rows and transfer with a multichannel pipette on the 96-well plate. Pipette the first sample into the wells A1 to H1 and so on. A total of 10 colonies per plate can be analysed (A11-H11: NTC; A12-H12 Pos. Control.).

If 11-12 different PCRs per colony are run simultaneously on a 96-well plate switch the pipetting scheme by 45° (max. 6 colonies can be analysed).

7.2. Programming of the real-time PCR instrument

Please find further information on programming the real-time PCR instrument in the respective operator's manual. Please be aware that some PCR-platforms have to be calibrated with the corresponding dye before performing multiplex-PCR.

Experiment Properties	Reagents: SYBR® Green Reagents
	Ramp speed: Standard
Select dyes	Reporter SYBR, Quencher None
Select passive reference	ROX
Sample Volume	20 µl

Temperature Profile

Holding Stage Cycles: 1 Analysis: None	Holding Stage Cycles: 1 Analysis: None	Cycling Stage Cycles: 35 Analysis: Quantification Acquisition at 60°	Melt Curve Stage Continuous
50°C 2 min	95°C 10 min	95°C 15 sec 60°C 1 min	95°C 15 sec 60°C 1 min 95°C 30 sec 60°C 15 sec

8. Interpretation of Melting Curve data


For analysis of PCR results gained with BactoReal® Typing Kit *E. coli* please verify the amplification plot and melt curves.


A specific amplification and melting curve of the sample has to be interpreted in context of the Ct/Cp and Tm [°C] values of the negative control.

	Target 1	Melting Point Tm around [°C]	Target 2	Melting Point Tm around [°C]
Tube 1 & 4	F4	83.0	F18	86.0
Tube 2 & 11	F5	84.0	fimA	89.0
Tube 3 & 5	F6	79.5	F41	83.0
Tube 6 & 12	STa	79.5	AIDA	85.0
Tube 7 & 8	STb	80.0	Stx2e	82.5
Tube 9 & 10	LT	80.5	EAST	87.0
Tube 13&14	pAA	80.5	escV	84.5
Tube 15	cnf1	80.0		
Tube 16	iucD	80.0		
Tube 17	papC	85.0		
Tube 18	pic	82.0		

Annotation: Due to sequence variation and other factors a shift of melting point Tm can be observed.

It is possible that due to a shift of Tm the two melting points can't be identified in the multiplex PCRs. Please set up the assay mix as singleplex reaction.





Melt curve analysis was performed on an Applied Biosystems® 7500 instrument (Thermo Fisher Scientific).

Annotation: Due to the multiplexing of primers a formation of primer dimers can be observed.

9. Annex – Symbols



Batch code



Catalogue number



Contains sufficient for <n> tests



Use by



Manufactured by



Store at

Permitted Use: The 5x HR Master Mix contains 5x HOT FIREPol® SolisGreen qPCR Mix which is a product of Solis BioDyne OÜ of 51014 Tartu, Estonia (www.sbd.ee). CYGREEN dye is a component of this qPCR Master Mix. CYGREEN is a U.S. registered trademark of Enzo Life Sciences, Inc. U.S. Patent Nos. 8,153,802 and 7,569,695 CYGREEN is a U.S. Registered trademark of Enzo Life Sciences, Inc. and this dye is the subject of the U.S. patents 12,456,502 and 7,569,695.