# **ENDOSTATIN**

**(EN)** ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF  
ENDOSTATIN IN HUMAN URINE, SERUM, CITRATE PLASMA, EDTA PLASMA,  
AND HEPARIN PLASMA  
CAT. NO. BI-20742 . 12 X 8 TESTS

**(DE)** ENZYM IMMUNOASSAY ZUR QUANTITATIVEN BESTIMMUNG VON  
ENDOSTATIN IN HUMANEM URIN, SERUM, CITRAT PLASMA, EDTA PLASMA  
UND HEPARIN PLASMA  
KAT. NR. BI-20742 . 12 X 8 TESTS

FOR RESEARCH USE ONLY  
NOT FOR USE IN DIAGNOSTIC PROCEDURES

rev.no. 200819 (replacing 190402)

This kit was developed and manufactured by:  
Biomedica Medizinprodukte GmbH, A-1210 Wien, Divischgasse 4  
Tel. +43/1/291 07 45, Fax +43/1/291 07 6389, E-mail [info@bmgrp.com](mailto:info@bmgrp.com)

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[www.bmgrp.com](http://www.bmgrp.com)

## 1) INTRODUCTION

Endostatin, a 20-kDa C-terminal proteolytic fragment of collagen XVIII, is an endogenous angiogenesis inhibitor localized in the vascular basement membrane in various organs (<http://www.uniprot.org/uniprot/P39060>). The biological functions of the endostatin-network involve SPARC, thrombospondin-1, glycosaminoglycans, collagens, and integrins. Endostatin is expressed during the progression of renal fibrosis in tubular cells of injured tissue. In renal micro-vascular disease, observed in late stages of patients with chronic kidney disease, increased endostatin levels are possibly the consequence of enhanced extracellular matrix degradation. Thus, endostatin may become an important marker for progressive microvascular renal disease in patients with chronic kidney disease. Endostatin levels in blood are also likely to increase in patients with other microvascular tissue injuries, including atherosclerosis, myocardial- and brain ischemia. In ischemic stroke patients, high endostatin plasma levels predict a worse long-term clinical outcome. Endostatin has recently been shown to be an independent risk factor of graft loss after kidney transplant as demonstrated in a prospective observational cohort study in 574 maintenance kidney transplant recipients.

### Areas of interest:

- Micro-vascular injury
- Atherosclerosis
- Graft loss
- Chronic kidney disease
- Ischemia

## 2) CONTENT OF THE KIT

CONT	KIT COMPONENTS	QUANTITY
PLATE	Goat polyclonal anti-Endostatin antibody, pre-coated microtiter strips in a strip holder, packed in an aluminium bag with desiccant	12 x 8 tests
WASHBUF	Wash buffer concentrate 20x, natural cap	1 x 50 ml
AB	Rabbit polyclonal anti-Endostatin antibody – biotin labelled, green cap, green dye, ready to use	1 x 7 ml
STD	Standards 1-7, (0; 25; 50; 100; 200; 400; 800 pmol/l), white caps, ready to use	7 x 400 µl
CTRL	Controls, A + B, (exact concentration on the label), yellow cap, ready to use	2 x 400 µl
ASYBUF	Assay Buffer, red cap, ready to use	1 x 100 ml
CONJ	Conjugate, (streptavidin-HRPO), amber bottle, amber cap, ready to use	1 x 22 ml
SUB	Substrate (TMB solution), amber bottle, blue cap, ready to use	1 x 22 ml
STOP	Stop solution, white cap, ready to use	1 x 7 ml

## 3) ADDITIONAL MATERIAL IN THE KIT

- 2 self-adhesive plastic films
- Protocol sheet
- Quality control protocol
- Instruction for use

## 4) MATERIAL AND EQUIPMENT REQUIRED BUT NOT SUPPLIED

- Precision pipettes calibrated to deliver 10 µl, 20 µl, 50 µl, 100 µl, 200 µl, 300 µl, 1000 µl and disposable tips
- Distilled or deionised water
- Plate washer is recommended for washing, alternative multichannel pipette or manifold dispenser
- ELISA reader capable of measuring absorbance at 450 nm (with correction wavelength at 630 nm)

## 5) REAGENTS AND SAMPLE PREPARATION

All reagents as supplied in the kit are stable at 4°C (2-8°C) until expiry date stated on the label of each reagent. Serum, plasma (EDTA, heparin, citrate), and urine are suitable for use in this assay. Do not change sample type during studies.

### Sample preparation:

Urine samples: We recommend first-morning-void urine. **Use urine samples undiluted.**

If urine samples read higher than STD7, it is recommended to dilute urine samples with assay buffer 1+3 and test again.

### Serum and plasma samples:

Collect venous blood samples by using standardized blood collection tubes. Perform plasma or serum separation by centrifugation according to supplier's instructions of the blood collection devices. Assay the acquired samples immediately or aliquot and store at -25°C or lower. Lipemic or haemolyzed samples may give erroneous results. Samples are stable for up to 4 freeze-thaw cycles. Samples should be mixed well before assaying. We recommend duplicates for all values.

**Serum and plasma samples must be diluted 1+100 with assay buffer (ASYBUF), eg. 10 µl sample + 1000 µl ASYBUF.** Diluted samples are stable at 4°C (2-8°C) overnight. Thus, dilutions can be prepared one day before analysis. Dilution factor has to be considered when calculating the final concentration of the sample.


For further information on sample stability please visit our website [www.bmgrp.com](http://www.bmgrp.com) (see Validation Data) or contact our customer service by e-mail [info@bmgrp.com](mailto:info@bmgrp.com) or by phone +43/ 1/ 29107-45.

### **Reconstitution/Handling:**

**WASHBUF (Wash buffer):** Dilute the concentrate 1:20, e.g. 50 ml WASHBUF + 950 ml distilled water. Crystals in the buffer concentrate will dissolve at room temperature. The undiluted WASHBUF is stable at 4°C (2-8°C) until expiry date stated on label. The diluted WASHBUF is stable up to one month at 4°C (2-8°C). Only use diluted WASHBUF when performing the assay.

## **6) PRINCIPLE OF THE ASSAY**

This kit is a sandwich enzyme immunoassay for the determination of Endostatin in human urine, serum or plasma samples. In a first step, STD/sample/CTRL and detection antibody (rabbit polyclonal anti-human Endostatin-Biotin) are pipetted into the wells of the microtiter strips, which are pre-coated with anti-Endostatin antibody. Endostatin present in the STD/sample/CTRL binds to the pre-coated antibody in the well and forms a sandwich with the detection antibody. In the washing step all non-specific unbound material is removed. In a second step, the conjugate (Streptavidin-HRP) is pipetted into the wells and reacts with the detection antibody. After another washing step, the substrate (TMB Tetramethylbenzidine) is pipetted into the wells. The enzyme catalysed color change of the substrate is directly proportional to the amount of Endostatin. This color change is detectable with a standard microtiter plate ELISA reader. The concentration of Endostatin in the sample is determined directly from the dose response curve.

## **7) ASSAY PROTOCOL**

All reagents and samples must be at room temperature (18-24°C) before use in the assay.

Mark position for STD/SAMPLE/CTRL (Standard/Sample/Control) on the protocol sheet.

Take microtiter strips out of the aluminium bag. Store unused strips with desiccant at 4°C (2-8°C) in the aluminium bag. Strips are stable until expiry date stated on the label.


1. Pipette 100 µl ASYBUF (Assay Buffer, red cap) into each well.
2. Add 20 µl of STD/pre-diluted\* SAMPLE/CTRL (Standard/Sample/Control) in duplicate into respective well.  
*\*Use urine samples undiluted. Use 1+100 pre-diluted serum and plasma samples, see chapter 5) Reagents and sample preparation.*
3. Add 50 µl AB (Antibody) into each well. Swirl gently.
4. **Cover tightly and incubate 3 hours at room temperature (18-24°C).**
5. Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer). After final wash, remove remaining WASHBUF by strongly tapping plate against paper towel.

6. Add 200 µl CONJ (Conjugate, amber cap) into each well.
7. **Cover tightly and incubate for 1 hour at room temperature (18-24°C).**
8. Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer). After final wash, remove remaining WASHBUF by strongly tapping plate against paper towel.
9. Add 200 µl SUB (Substrate, blue cap) into each well.
10. **Incubate for 30 min at room temperature (18-24°C) in the dark.**
11. Add 50 µl STOP (Stop solution, white cap) into each well. Swirl gently.
12. Measure absorbance immediately at 450 nm with reference 630 nm, if available.

## 8) CALCULATION OF RESULTS

Read the optical density (OD) of all wells with a plate reader using 450 nm wavelength (correction wavelength 630 nm). Construct the standard curve from the OD values of the STD. Use commercially available software or graph paper. Obtain sample concentration from this standard curve. The assay was evaluated with 4PL algorithm. Different curve fitting methods need to be evaluated by the user. Respective dilution factors, e.g. 1:101 for serum and plasma samples, have to be considered when calculating the final concentration of the sample.

### Example typical STD-curve:

The quality control (QC) protocol supplied with the kit shows the results of the final release QC for each kit at production date. Data for OD obtained by customers may differ due to various influences and/or due to the normal decrease of signal intensity during shelf life. However, this does not affect validity of results as long as an OD of 1.50 or more is obtained for the STD with the highest concentration and the value of the CTRL is in range (target range see label).

## 9) ASSAY CHARACTERISTICS

Method:	Sandwich ELISA, HRP/TMB, 12x8-well strips	
Sample type:	Urine, serum, EDTA plasma, citrate plasma, and heparin plasma	
Standard range:	0 to 800 pmol/l (7 standards and 2 controls)	
Conversion factor:	Endostatin: 1 pg/ml = 0.05 pmol/l (MW: 20 kDa)	
Sample volume:	20 µl / well (urine: undiluted; serum and plasma: diluted 1+100 in ASYBUF)	
Incubation time:	3 h / 1 h / 30 min	
Sensitivity:	LOD: (0 pmol/l + 3 SD); 4 pmol/l; LLOQ: 3 pmol/l	
Specificity:	This assay recognizes endogenous and recombinant human Endostatin. There is no cross reactivity with human recombinant Collagen Type XV (COL15).	
Precision:	Intra-assay (n=5) ≤ 6% Inter-assay (n=16) ≤ 5%	
	Urine (n=8): 87%	
	Serum (n=16): 120%	EDTA plasma (n=12): 106%

Spike/Recovery (average recovery spiked with 100 pmol/l rec. Endostatin):	Citrate plasma (n=8): 85%	Heparin plasma (n=4): 96%	
Dilution linearity (average recovery of expected Endostatin after dilution):	<u>Dilution:</u>	<u>1+1</u>	<u>1+3</u>
	Urine (n=6)	103%	123%
	Serum (n=4)	116%	125%
	EDTA plasma (n=4)	120%	134%
	Citrate plasma (n=4)	101%	112%
Values from apparently healthy individuals:	Heparin plasma (n=4) 111% 130%		
	Median urine (n=789): 20 pmol/l Median serum (n=59): 5151 pmol/l (5.2 nmol/l) Median citrate plasma (n=30): 4747 pmol/l (4.7 nmol/l)  Each laboratory should establish its own reference range for the samples under investigation. Do not change sample type during the study.		

For further information on assay characteristics please visit our website [www.bmgrp.com](http://www.bmgrp.com) (see Validation Data) or contact our customer service by e-mail [info@bmgrp.com](mailto:info@bmgrp.com) or by phone +43/ 1/ 29107-45.

### 10) Precision

Intra-assay: 2 samples of known concentrations were tested 5 times within 1 kit lot by 1 operator.

Inter-assay: 2 samples of known concentrations were tested 16 times within 2 different kit lots and by 3 different operators.

Intra-assay (n=5)	Sample 1	Sample 2	Inter-assay (n=16)	Sample 1	Sample 2
Mean (pmol/l)	52	416	Mean (pmol/l)	50	405
SD (pmol/l)	0.9	23.6	SD (pmol/l)	2.1	21.6
CV (%)	2	6	CV (%)	4	5

### 11) TECHNICAL HINTS

- Do not mix or substitute reagents with those from other lots or sources.
- Do not mix stoppers and caps from different reagents or use reagents between lots.
- Do not use reagents beyond expiration date.
- Protect reagents from direct sunlight.
- Substrate solution should remain colorless until added to the plate.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.

### 12) PRECAUTIONS

All test components of human source were tested against HIV-Ab, HCV-Ab and HBsAg and were found to be negative. Nevertheless, they should be handled and disposed as if they were infectious.

Liquid reagents contain ≤0.1% Proclin 950 as preservative. Avoid contact with skin and mucous membrane. Proclin 950 is not toxic in concentrations used in this kit. It may cause allergic skin reactions – avoid contact with skin or eyes.

- Do not pipette by mouth.
- Do not eat, drink, smoke or apply cosmetics where reagents are used.
- Wear gloves, glasses and lab jacket while performing this assay.
- Sulfuric acid is irritating to the eyes and skin. Avoid contact with skin and mucous. Irritations are possible. – Flush with water if contact occurs!!

### 13) LITERATURE

1. Endostatin Is an Independent Risk Factor of Graft Loss after Kidney Transplant. *Chu C et al., 2020, Am J Nephrol 51:373–380.*
2. Validation of an enzyme-linked immunosorbent assay (ELISA) for quantification of endostatin levels in mice as a biomarker of developing glomerulonephritis. *Wallwitz J et al., 2019, PLoS One, 14(8): e0220935.*
3. Higher Parathyroid Hormone Level Is Associated with Increased Arterial Stiffness in Type 1 Diabetes *Zobel EH et al., 2017, Diabetes Care 40(3):dc162428.*
4. Endostatin in Chronic Kidney Disease: Associations with Inflammation, Vascular Abnormalities, Cardiovascular Events and Survival. *Kanbay M et al., 2016 European Journal of Internal Medicine : 81–87.*
5. Elevated plasma levels of endostatin are associated with chronic kidney disease. *Chen et al., 2012, Am J Nephrol 35(4): 335-340.*
6. Early-onset coronary artery disease after pediatric kidney transplantation: implicating the angiogenesis inhibitor, endostatin. *Igbal CW et al., 2011, Am Surg 77(6): 731-735.*
7. A defective angiogenesis in chronic kidney disease. *Futrakul N et al., 2008, Ren Fail 30(2): 215-217.*
8. Excretion of anti-angiogenic proteins in patients with chronic allograft dysfunction. *Moskowitz-Kassai E et al., 2012, Nephrol Dial Transplant 27(2): 494-497.*
9. Endostatin and angiostatin are increased in diabetic patients with coronary artery disease and associated with impaired coronary collateral formation. *Sodha NR et al., 2009, Am J Physiol Heart Circ Physiol, 296: H428-H434.*
10. A large screening of angiogenesis biomarkers and their association with neurological outcome after ischemic stroke. *Navarro-Sobrinho M et al., 2011, Atherosclerosis, 216(1): 205-211.*

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## **ASSAY PROTOCOL AND CHECKLIST**

### **PREPARATION OF REAGENTS:**

- Bring all reagents to room temperature (18-24°C).
- Prepare reagents and samples as instructed.
- Take microtiter strips out of the aluminium bag and mark positions on the protocol sheet.
- Bring unused and prepared components to the storage temperature mentioned in the package insert.

### **TEST PROCEDURE:**

- Step 1) Pipette 100 µl ASYBUF (red cap) into respective well.
- Step 2) Add 20 µl of **STD/pre-diluted\* SAMPLE/CTRL** (Standard/Sample/Control) in duplicate into respective well.
- Step 3) Add 50 µl AB (biotinylated anti-Endostatin antibody, green cap, green dye) into each well, swirl gently.
- Step 4) Cover tightly and incubate 3 hours at room temperature (18-24°C).**
- Step 5) Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer). After final wash, remove remaining WASHBUF by strongly tapping plate against paper towel.
- Step 6) Add 200 µl CONJ (Conjugate, amber cap) into each well.
- Step 7) Cover tightly and incubate for 1 hour at room temperature (18-24°C).**
- Step 8) Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer). After final wash, remove remaining WASHBUF by strongly tapping plate against paper towel.
- Step 9) Add 200 µl SUB (Substrate, blue cap) into each well.
- Step 10) Incubate for 30 min at room temperature (18-24°C) in the dark.**
- Step 11) Add 50 µl STOP (Stop solution, white cap) into each well. Swirl gently.
- Step 12) Read Optical Density

*\*For sample dilution see chapter 5) Reagent and Sample Preparation.*