ENDOSTATIN MOUSE/RAT

(EN) ELISA FOR THE QUANTITATIVE DETERMINATION OF MOUSE/RAT
ENDOSTATIN IN SERUM AND PLASMA
Cat. No. BI-20742MR . 12 x 8 TESTS

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FOR RESEARCH USE ONLY
NOT FOR USE IN DIAGNOSTIC PROCEDURES

rev.no. 190111 (replacing 170424)

This kit was developed and manufactured by:
Biomedica Medizinprodukte GmbH, A-1210 Wien, Divischgasse 4
Tel. +43/1/291 07 45, Fax +43/1/291 07 6389, E-mail info@bmgrp.com

CONTENT / INHALT

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Detailed information on the Endostatin mouse/rat ELISA, e.g. assay validation data, sample matrix comparisons, and stability data is available on our website.

www.bmgrp.com

1) INTRODUCTION

Endostatin, a 20-kDa C-terminal proteolytic fragment of collagen XVIII, is an endogenous angiogenesis inhibitor localized in the vascular basement membrane in various organs (<http://www.uniprot.org/uniprot/P39060>). The biological functions of the endostatin-network involve SPARC, thrombospondin-1, glycosaminoglycans, collagens, and integrins.

In animal studies, renal Endostatin expression preceded deteriorating kidney function and induced renal fibrosis in aging mice. In humans, Endostatin is expressed during the progression of renal fibrosis in tubular cells of injured tissue. In renal micro-vascular disease, observed in late stages of patients with chronic kidney disease, increased endostatin levels are possibly the consequence of enhanced extracellular matrix degradation. Thus endostatin may become an important marker for progressive microvascular renal disease in patients with chronic kidney disease. Endostatin levels in blood are also likely to increase in patients with other microvascular tissue injuries, including atherosclerosis, myocardial- and brain ischemia. In ischemic stroke patients, high endostatin plasma levels predict a worse long-term clinical outcome. In a cohort of critically ill patients, plasma endostatin improved AKI prediction based on clinical risk factors. Endostatin has evolved as a molecular target and is currently under investigation in clinical trials.

Areas of interest:

- Micro-vascular injury
- Atherosclerosis
- Sepsis
- Chronic kidney disease
- Ischemia
- Preeclampsia

2) CONTENTS OF THE KIT

CONTENTS	KIT COMPONENTS	QUANTITY
PLATE	Anti mouse endostatin pre-coated microtiter strips in a strip holder, packed in an aluminium bag with desiccant	12 x 8 tests
WASHBUF	Wash buffer concentrate 20x, natural cap	1 x 50 ml
ASYBUF	Assay buffer, red cap, ready to use	1 x 30 ml
STOCK STD	Stock standard (32 nmol/l), recombinant mouse endostatin, red cap, lyophilized	1 vial
CTRL	Control, yellow cap, lyophilized (exact concentration after reconstitution see label)	1 vial
CONJ	Conjugate, (goat anti-endostatin-HRPO), amber bottle, amber cap, ready to use	1 x 7 ml
SUB	Substrate (TMB solution), amber bottle, blue cap, ready to use	1 x 13 ml
STOP	Stop solution, white cap, ready to use	1 x 7 ml

3) ADDITIONAL MATERIAL IN THE KIT

- 1 self-adhesive plastic film
- Quality control protocol
- Protocol sheet
- Instruction for use

4) MATERIAL AND EQUIPMENT REQUIRED BUT NOT SUPPLIED

- Precision pipettes calibrated to deliver 5 µl, 10 µl, 50 µl, 150 µl, 250 µl, 500 µl and disposable tips
- Distilled or deionised water
- Plate washer is recommended for washing, alternative multichannel pipette or manifold dispenser
- Refrigerator with 4°C (2-8°C)
- ELISA reader capable of measuring absorbance at 450 nm (with correction wavelength at 630 nm)
- Graph paper or software for calculation of results

5) REAGENTS AND SAMPLE PREPARATION

All reagents as supplied in the kit are stable at 4°C (2-8°C) until expiry date stated on the label of each reagent. Bring all reagents to room temperature before use.

Reconstituted STOCK STD and CTRL are stable at -25°C or lower until expiry date on label. STOCK STD and CTRL can undergo 3 freeze-thaw cycles.

Serum and plasma are suitable for use in this assay. Do not change sample type during studies.

A. Sample preparation:

Collect venous blood samples by using standardized blood collection tubes for serum or plasma. Perform serum and plasma separation by centrifugation according to supplier's instructions of the blood collection devices and measure the acquired serum or plasma samples as soon as possible. For longer storage aliquot and store at -25°C or lower. Samples are stable for 3 freeze-thaw cycles. Thawed samples should be assayed as soon as possible. Lipemic or haemolysed samples may give erroneous results. Samples should be mixed well before assaying. We recommend duplicates for all values.

Samples must be diluted 1+50 with ASYBUF (assay buffer) prior to the assay, e.g. 5 µl sample + 250 µl ASYBUF (or 10 µl sample + 500 µl ASYBUF). *Note: 50 µl pre-diluted sample is required / well.*

B. Preparation of CTRL (control):

Reconstitute the CTRL (control) in 200 µl ASYBUF (assay buffer), leave at room temperature (18-26°C) for 15 min and mix well prior to making dilution. The exact concentration is stated on the label.

CTRL must be diluted 1+50 with ASYBUF prior to the assay, e.g. 5 µl CTRL + 250 µl ASYBUF (or 10 µl CTRL + 500 µl ASYBUF). *Note: 50 µl pre-diluted CTRL is required / well.*

C. Preparation of STOCK STD (stock standard):

Reconstitute the Endostatin STOCK STD (stock standard) in 200 µl ASYBUF (assay buffer). Leave at room temperature (18-26°C) for 15 min and mix well prior to making dilutions.

D. Preparation of the standard curve:


Use polypropylene tubes

- Mark tubes as shown below (Graph1).
- Pipette 500 µl of ASYBUF (assay buffer) into tube marked as STD7.
- Pipette 150 µl of ASYBUF each into tubes marked as STD6 to STD1.
- Pipette 10 µl of the reconstituted STOCK STD into tube marked as STD7. Mix thoroughly.
- Prepare a two-fold serial dilution to obtain STD6 to STD2 (Graph1).
- ASYBUF serves as the zero standard (=STD1, 0 nmol/l).

Note: 50 µl pre-diluted STD is required / well.

Always mix each tube thoroughly before the next step!

Graph1: preparation of the standard curve

Standards diluted for sample measurement

Established standard curve already includes the 1+50 assay pre-dilution, thus the sample results can be read directly from the standard curve (also see chapter calculation of results).

Pre-diluted STDs, CTRLs and samples should be measured as soon as possible (do not store or freeze).


Preparation of WASHBUF (wash buffer):

Dilute the concentrate 1:20 (e.g. 50 ml WASHBUF + 950 ml distilled water). Crystals in the buffer concentrate will dissolve at room temperature (18-26°C). The diluted wash buffer is stable at 4°C (2-8°C) for one month. Use only diluted wash buffer for the assay performance.

For further information on sample stability please visit our website www.bmgrp.com (see Validation Data) or contact our customer service by e-mail info@bmgrp.com or by phone +43/ 1/ 29107-45.

6) PRINCIPLE OF THE ASSAY

This kit is a sandwich enzyme immunoassay for the direct determination of endostatin in mouse/rat serum and plasma samples. In a first step, 1+50 pre-diluted STD/sample/CTRL and conjugate (goat anti-endostatin-HRPO) are pipetted into the wells of the microtiter strips, which are pre-coated with anti mouse endostatin antibody. Endostatin present in the STD/sample/CTRL binds to the pre-coated antibody in the well and forms a sandwich with the detection antibody in the conjugate. In the washing step all non-specific unbound material is removed. In a second step, the substrate (TMB Tetramethylbenzidine) is pipetted into the wells. The enzyme catalysed colour change of the substrate is directly proportional to the amount of endostatin. This colour change is detectable with a standard microtiter plate ELISA reader. A dose response curve of the absorbance (optical density, OD at 450 nm) vs. standard concentration is generated, using the values obtained from the standard. The concentration of mouse/rat endostatin in the sample is determined directly from the dose response curve.

7) ASSAY PROTOCOL

All reagents and samples must be at room temperature (18-26°C) before use in the assay.

Mark position for STD/SAMPLE/CTRL (Standard/Sample/Control) on the protocol sheet.

Take microtiter strips out of the aluminium bag. Store unused strips with desiccant at 4°C (2-8°C) in the aluminium bag. Strips are stable until expiry date stated on the label.

1. Pipette 50 µl of pre-diluted (1+50) STD/SAMPLE/CTRL (Standard/Sample/Control) (see reagents and sample preparation) in duplicate into respective well.

2. Add 50 µl CONJ (Conjugate, amber cap) into each well, swirl gently.

3. **Cover tightly and incubate for 2 hours at room temperature (18-26°C).**
4. Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer, natural cap). Remove remaining WASHBUF by strongly tapping plate against paper towel after the last wash.
5. Add 100 µl SUB (Substrate, blue cap) into each well, swirl gently.
6. **Incubate for 30 minutes at room temperature (18-26°C) in the dark.**
7. Add 50 µl STOP (Stop solution, white cap) into each well, swirl gently.
8. Measure absorbance immediately at 450 nm (with reference 630 nm, if available).


8) CALCULATION OF RESULTS

Read the optical density (OD) of all wells on a plate reader using 450 nm wavelength (correction wavelength 630 nm). Construct the standard curve from the OD values of the STD. Use commercially available software or graph paper. Obtain sample concentration from this standard curve. The assay was evaluated with logit-log and 4PL algorithm curve fitting. Different curve fitting methods need to be evaluated by the user.

Samples, control, and standards are all diluted 1+50 prior to the assay, so there is no need to take this dilution factor into account. Sample results can be read directly from the standard curve.

Note: Sample dilutions above 1+50 have to be considered when calculating the final sample concentration.

Example typical STD-curve:

The quality control (QC) protocol supplied with the kit shows the results of the final release QC for each kit lot. Data for OD obtained by customers may differ due to various influences and/or due to the normal decrease of signal intensity during shelf life. However, this does not affect validity of results as long as an OD of 1.00 or more is obtained for the STD with the highest concentration and the value of the CTRL is in range (target range see labels).

9) ASSAY CHARACTERISTICS

Method:	Sandwich ELISA, HRP/TMB, 12x8-well strips
Sample type:	Mouse or rat serum, plasma
Standard range:	0-32 nmol/l (0 / 1 / 2 / 4 / 8 / 16 / 32)
Conversion factor:	1 ng/ml = 0.049 nmol/l or 1 nmol/l=20.376 ng/ml (MW: 20.4 kDa)
Sample volume:	5 µl / sample
Incubation time:	2 h / 30 min – room temperature
Sensitivity:	LOD (0 nmol/l + 3 SD): 0.24 nmol/l; LLOQ: 0.5 nmol/l
Specificity:	This assay detects recombinant and endogenous mouse and rat Endostatin. Mouse and rat Endostatin share a 95.7% homology.

Precision:	Intra-assay (n=5) ≤ 9%, Inter-assay (n=15) ≤ 10%			
Spike/Recovery (average recovery spiked with 25 nmol/l recombinant mouse Endostatin):	Mouse serum (n=7): 95%	Rat serum (n=4): 97%		
	Mouse plasma (n=5): 91%	Rat plasma: n.a.		
Dilution linearity of recombinant and endogenous Endostatin (average recovery of expected) Endostatin values after a 1+1; 1+3; 1+7 dilution in ASYBUF):	Recovery (%):	Endostatin recombinant / endogenous		
	Dilution:	1+1	1+3	1+7
	Mouse serum (n=6)	108 / 101	117 / 88	114 / *
	Mouse plasma (n=5)	101 / 107	104 / 109	103 / *
	Rat serum (n=7)	99 / 96	n.a. / 89	n.a. / *
	Rat plasma (n=4)	n.a. / 89	n.a. / 80	n.a. / *
Values from various mouse and rat samples:	Mouse sera C57BL6.J0laHsd, 12 weeks (n=11): 6.7 ± 0.8 nmol/l Wildtype normal mice sera, 12 weeks, male (n=10): 5.4 ± 1.2 nmol/l Wildtype normal rat sera, 12 weeks, male (n=8): 2.5 ± 0.4 nmol/l <i>Each laboratory should establish its own reference range for the samples under investigation.</i>			

*not detectable, n.a.: not analysed

For further information on assay characteristics please visit our website www.bmgrp.com (see Validation Data) or contact our customer service by e-mail info@bmgrp.com or by phone +43/ 1/ 29107-45.

10) PRECISION

Intra-assay: 2 samples of known concentration were tested 5 times in 1 kit lot by 1 operator.

Inter-assay: 2 samples of known concentration were tested 15 times in 3 different assays in 3 days by 2 different operators.

Intra-assay (n= 5)	Sample 1	Sample 2	Inter-assay (n= 15)	Sample 1	Sample 2
Mean (nmol/l)	0.98	31.98	Mean (npmol/l)	1.01	31.99
SD (nmol/l)	0.09	0.55	SD (nmol/l)	0.10	0.48
CV (%)	9	2	CV (%)	10	2

Detailed information on the mouse/rat Endostatin ELISA, e.g. assay performance characteristics, matrix comparisons, and stability data is available on our website www.bmgrp.com (see Validation Data).

11) TECHNICAL HINTS

- Do not mix or substitute reagents with those from other lots or sources.
- Do not mix stoppers and caps from different reagents or use reagents between lots.
- Do not use reagents beyond expiration date.
- Protect reagents from direct sunlight.
- Substrate solution should remain colourless until added to the plate.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.

12) PRECAUTIONS

All liquid reagents contain ≤ 0.1% Proclin 950 as preservative. Avoid contact with skin and mucous membrane. Proclin 950 is not toxic in concentrations used in this kit. It may cause allergic skin reactions – avoid contact with skin or eyes.

- Do not pipette by mouth.
- Do not eat, drink, smoke, or apply cosmetics where reagents are used.
- Wear gloves, glasses, and lab coat while performing this assay.
- Sulfuric acid is irritating to the eyes and skin. Avoid contact with skin and mucous. Irritations are possible. Flush with water if contact occurs!

13) LITERATUR

1. Endostatin expression in the murine model of ischaemia/reperfusion-induced acute renal failure. *Bellini MH et al., Nephrology 2007; 12(5):459–65.*
2. Endostatin and kidney fibrosis in aging: a case for antagonistic pleiotropy? *Lin CHS et al., Am J Physiol Heart Circ Physiol 2014; 306:H1692–9.*
3. Plasma endostatin may improve acute kidney injury risk prediction in critically ill patients. *Mårtensson J et al., Ann Intensive Care 2016; 6:6.*
4. Endostatin in chronic kidney disease: Associations with inflammation, vascular abnormalities, cardiovascular events and survival. *Kanbay et al., European Journal of Internal Medicine 2016; 63: 81-87.*
5. The association between endostatin and kidney disease and mortality in patients with type 2 diabetes. *Carlsson AC et al., Diabetes & Metabolism 2016 ; 42(5) 351-357.*
6. Elevated plasma levels of endostatin are associated with chronic kidney disease. *Chen et al., Am J Nephrol 2012; 35(4): 335-340.*
7. Early-onset coronary artery disease after pediatric kidney transplantation: implicating the angiogenesis inhibitor, endostatin. *Igbal CW et al., Am Surg 2011; 77(6): 731-735.*
8. A large screening of angiogenesis biomarkers and their association with neurological outcome after ischemic stroke. *Navarro-Sobrinho M et al., Atherosclerosis 2011; 16(1): 205-211.*

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BI-20742MR ENDOSTATIN MOUSE/RAT

PREPARATION OF REAGENTS:

- Bring all reagents to room temperature (18-26°C).
- Prepare reagents and samples as instructed.
- Bring unused and prepared components to the storage temperature mentioned in the package insert.
- Take microtiter strips out of the aluminium bag and mark positions on the protocol sheet.

TEST PROCEDURE:

- Step 1) Pipette 50 µl of pre-diluted (1+50) STD/SAMPLE/CTRL (Standard/Sample/Control) (see reagents and sample preparation) in duplicate into respective well.
- Step 2) Add 50 µl CONJ (Conjugate, amber cap) into each well, swirl gently.
- Step 3) Cover tightly and incubate for 2 hours at room temperature (18-26°C).**
- Step 4) Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer, natural cap). Remove remaining WASHBUF by strongly tapping plate against paper towel after the last wash.
- Step 5) Add 100 µl SUB (Substrate, blue cap) into each well, swirl gently.
- Step 6) Incubate for 30 minutes at room temperature (18-26°C) in the dark**
- Step 7) Add 50 µl STOP (Stop solution, white cap) into each well, swirl gently.
- Step 8) Measure absorbance immediately at 450 nm with reference 630 nm, if available.