NovaLisa[®]

Chikungunya Virus IgM μ-capture

ELISA

CE

Only for in-vitro diagnostic use

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ENGLISH

1. INTENDED USE

The Chikungunya Virus IgM μ -capture ELISA is intended for the qualitative determination of IgM class antibodies against Chikungunya Virus in human serum or plasma (citrate, heparin).

2. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific IgM-class antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) μ -capture technique.

Microtiterplates are coated with anti-human IgM antibodies to bind the corresponding antibodies of the sample. After washing the wells to remove all unbound sample material, antigen is added. This antigen binds to the captured specific IgM antibodies. After a further washing step biotinylated antibody is pipetted into the wells. After washing horseradish peroxidase (HRP) labelled streptavidin is added that binds to the captured specific immune complex. After a further washing step, the immune complexes are visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is proportional to the amount of specific IgM antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

3. MATERIALS

3.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with anti-human IgM; in resealable aluminium foil.
- **DIL:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 \pm 0.2; coloured yellow; ready to use; white cap; \leq 0.0015% (v/v) CMIT/ MIT (3:1).
- **SOLN STOP:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **WASH BUF 20x:** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 \pm 0.2, for washing the wells; white cap; 0.2% (w/v) 5-Bromo-5-nitro-1,3-dioxane.
- **Antigen, lyophilized:** 6 bottles containing lyophilized Chikungunya virus Antigen Solution; red cap; \leq 0.02% (v/v) MIT.
- **Antibody Solution:** 1 bottle containing 6 mL of biotinylated Chikungunya virus antibody, ready to use; coloured blue; white cap; \leq 0.02% (v/v) MIT.
- **Conjugate:** 1 bottle containing 6 mL Streptavidin conjugated with peroxidase, ready to use; coloured red; black cap.
- **SUB TMB:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 1.5 mL control; coloured yellow; ready to use; red cap; \leq 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; green cap; \leq 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 1.5 mL control; coloured yellow; ready to use; blue cap; \leq 0.0015% (v/v) CMIT/ MIT (3:1).

For hazard and precautionary statements see 11.1

3.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

3.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000 μ L
- Vortex tube mixer
- Distilled water
- Disposable tubes

4. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

5. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

5.1. Microtiterplate

The break-apart snap-off strips are coated with anti-human IgM. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

5.2. Antigen, lyophilized

The bottles contain Antigen, lyophilized. The content of each vial has to be resolved in 1 mL **WASH BUF 1x** by turning it slowly (no vortex) and 15 min incubation at room temperature (20...25 °C). The reconstituted solution is stable for 1 day at 2...8 °C.

5.3. **WASH BUF 20x**

Dilute **WASH BUF 20x** 1 + 19; e. g. 10 mL **WASH BUF 20x** + 190 mL distilled water. The diluted buffer (**WASH BUF 1x**) is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

5.4. **SUB TMB**

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. **SUB TMB** should be colourless or could have a slight blue tinge. If **SUB TMB** turns into blue, it may have become contaminated and should be thrown away.

6. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise, they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

6.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with **DIL**. Dispense 10 µL sample and 1 mL **DIL** into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

7. ASSAY PROCEDURE

Please read the instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of **WASH BUF 1x** from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 11. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 50 µL standards/controls and diluted samples into their respective wells. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of **WASH BUF 1x**. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results
5. Dispense 50 µL reconstituted Chikungunya virus Antigen into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).**
7. Repeat step 4.
8. Dispense 50 µL Chikungunya virus Antibody Solution into all wells except for the Blank well A1.
9. **Incubate for 30 min at room temperature (20...25 °C).**
10. Repeat step 4.
11. Dispense 50 µL Streptavidin peroxidase conjugate into all wells except for the Blank well A1.
12. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
13. Repeat step 4.
14. Dispense 100 µL **SUB TMB** into all wells.
15. **Incubate for exact 15 min. at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
16. Dispense 100 µL **SOLN STOP** into all wells in the same order and at the same rate as for **SUB TMB**, thereby a colour change from blue to yellow occurs.
17. Measure the absorbance at 450/620 nm within 30 min after addition of **SOLN STOP**.

7.1. Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

8. RESULTS

8.1. Run Validation Criteria

In order for an assay run to be considered valid, these instructions for use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value < **0.100**
- **Negative Control:** Absorbance value < **Cut-off**
- **Cut-off Control:** Absorbance value **0.150 – 1.300**
- **Positive Control:** Absorbance value > **Cut-off**

If these criteria are not met, the test is not valid and must be repeated.

8.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43

$$\text{Cut-off} = 0.43$$

8.2.1. Results in Units [NTU]

$\frac{\text{Sample (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec-Units} = \text{NTU}]$

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ NTU}$

8.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks.
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.		

8.3.1. Antibody Isotypes and State of Infection

Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection

9. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

9.1. Precision

Intraassay	n	Mean (OD)	CV (%)
#1	24	0.287	8.53
#2	24	0.769	5.82
#3	24	0.618	5.50
Interassay	n	Mean (NTU)	CV (%)
#1	12	32.42	4.88
#2	12	26.96	5.82
#3	12	5.95	12.84

9.1. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100% (95% confidence interval: 95.01% - 100%).

9.2. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100% (95% confidence interval: 96.19% - 100%).

9.3. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

9.4. Cross Reactivity

Cross reactivity with antibodies against Borrelia, CMV and Toxoplasma cannot be excluded. Interference with polyclonal stimulation of EBV infections is likely. In the presence of infectious Mononucleosis (Pfeiffer's Disease, EBV infection) polyclonal stimulation of B lymphocytes can occur. This may result in non-specific reactions in the detection of antibodies of the IgM class. Therefore, it is recommended to exclude an EBV infection by differential diagnosis. Cross reactivity with antibodies against other alpha viruses cannot be excluded.

10. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

11. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- **The Chikungunya virus antigens are inactivated. All materials should still be regarded and handled as potentially infectious. Wear gloves while performing the test. We recommend using the antigen under BSL2 cabinet (clean bench).**
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

11.1. Safety note for reagents containing hazardous substances

Reagents may contain CMIT/MIT (3:1) or MIT (refer to 3.1).
Therefore, the following hazard and precautionary statements apply.

	Warning	H317	May cause an allergic skin reaction.
		P261	Avoid breathing spray.
		P280	Wear protective gloves/ protective clothing.
		P302+P352	IF ON SKIN: Wash with plenty of soap and water.
		P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
		P362+P364	Take off contaminated clothing and wash it before reuse

Reagents may contain 5-Bromo-5-nitro-1,3-dioxane (refer to 3.1)
Therefore, the following hazard and precautionary statements apply.

	Warning	H315	Causes skin irritation.
		H319	Causes serious eye irritation
		P280	Wear protective gloves/ protective clothing.
		P302+P352	IF ON SKIN: Wash with plenty of soap and water.
		P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
		P337+P313	If eye irritation persists: Get medical advice/attention.

Further information can be found in the safety data sheet.

11.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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Chikungunya Virus IgM μ -capture

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all samples and standards/controls.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate Blank (A1)	Negative control	Cut-off control	Positive control	Sample (diluted 1+100)
Negative control	-	50 μ L	-	-	-
Cut-off control	-	-	50 μ L	-	-
Positive control	-	-	-	50 μ L	-
Sample (diluted 1+100)	-	-	-	-	50 μ L
Cover wells with foil supplied in the kit					
Incubate for 1 h at 37 \pm 1 $^{\circ}$C					
Wash each well three times with 300 μ L of WASH BUF 1x					
Reconstituted Antigen	-	50 μ L	50 μ L	50 μ L	50 μ L
Incubate for 30 min at room temperature (20...25 $^{\circ}$C)					
Wash each well three times with 300 μ L of WASH BUF 1x					
Antibody Solution	-	50 μ L	50 μ L	50 μ L	50 μ L
Incubate for 30 min at room temperature (20...25 $^{\circ}$C)					
Wash each well three times with 300 μ L of WASH BUF 1x					
Streptavidin conjugate	-	50 μ L	50 μ L	50 μ L	50 μ L
Incubate for 30 min at room temperature (20...25 $^{\circ}$C)					
Do not expose to direct sunlight					
Wash each well three times with 300 μ L of WASH BUF 1x					
SUB TMB	100 μ L	100 μ L	100 μ L	100 μ L	100 μ L
Incubate for exactly 15 min at room temperature (20...25 $^{\circ}$C) in the dark					
SOLN STOP	100 μ L	100 μ L	100 μ L	100 μ L	100 μ L
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

CHIM0590_IFU_rev01_fromLot_112N