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FD Darm check professionell

Literatur

- Knittel, T.; Rudi, V.; Lutz, C.; Abels, L.; Frost, G.; High prevalence of elevated fecal calprotectin levels in stool samples from patients with diverticular disease: (2017) UEG Week; *Gastroenterologie*, May 2018, Volume 154, Issue 6, Supplement 1, Page S664
- Knittel, Th., Lüthgens, K.; Higher sensitivity for the detection of adenomas by a Hb/Hp bed site test using a novel sampling device (2013)
- Tao, S., Brenner, H.; Well adjusted qualitative immunochemical faecal occult blood tests could be a promising alternative for inexpensive, high-quality colorectal cancer screening; *European Journal of Cancer Prevention* 2013, 22:305-310 Sieg A., Thoms C., Lüthgens K., John M.R., Schmid Gayk H.; Detection of colorectal neoplasm by the highly sensitive hemoglobin-haptoglobin complex in feces; *Int J Colorectal Dis* (1999) 14: 267-271
- Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; *Archive of Internal Medicine*; (1995) 155:2389-2402
- Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; *British Medical Journal*; (1988) 296:1092-1094
- Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? *Cancer Supplement*; (1994) 74(7) : 2023-2027
- Miller, A.B.; An Epidemiological Perspective on Cancer Screening; *Clinical Biochemistry* (1995) 28(1) : 41-48
- Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
- Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; *Morbidity and Mortality Weekly Report*; (1995) 45 (5) : 107- 110
- St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; *Gastroenterology*; (1993) 104:1661-1668
- Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: *Hepatogastroenterology*; 2000 Mar-Apr; 47(32):396-9.
- and *Mortality Weekly Report*; (1995) 45 (5) : 107-110
- Hundt S, Haug U, Brenner H; *American College of Physicians* 2009; 150(3):162-169
- Devlin SM and Panaccione R: Evolving inflammatory bowel disease treatment paradigms: top-down versus step-up. *Med Clin North Am*. 2010;94:1-18
- Pineton de Chambrun G et al.: Clinical implications of mucosal healing for the management of IBD. *Nat Rev Gastroenterol Hepatol* 2010; 7(1):15-29
- Lichtenstein GR and Rutgeerts P: Importance of mucosal healing in ulcerative colitis. *Inflamm Bowel Dis*. 2010;16:338-346
- Smith MA et al.: Pharmacogenomics in the treatment of inflammatory bowel disease. *Pharmacogenetics*, 2010;11(3):421-437
- Lin MV et al.: What is the optimal therapy for Crohn's disease: step-up or top-down? *Expert Rev Gastroenterol Hepatol*. 2010;4(2):167-180
- Strauch U and Schölmerich J.: Emerging drugs to treat Crohn's disease. *Expert Opin Emerg Drugs*, 2010;15(2):309-322
- Isaacs KL: How rapidly should remission be achieved? *Dig Dis* 2010;28(3):548-555
- Schwartz M and Regueiro M: Prevention and treatment of postoperative Crohn's disease recurrence: an update for a new decade. *Curr Gastroenterol Rep*. 2011 Feb;13(1):95-100
- Ha C and Kornbluth A: Mucosal healing in inflammatory bowel disease: where do we stand? *Curr Gastroenterol Rep*. 2010;12(6):471-478.
- Fagerberg UL et al.: Fecal calprotectin: a quantitative marker of colonic inflammation in children with inflammatory bowel disease. *J Pediatr Gastroenterol Nutr*. 2007;45(4):414-420
- Rutgeerts P et al.: Biological therapies for inflammatory bowel diseases. *Gastroenterology*, 2009;136(5):1182-1197
- Jalocha L et al.: Mucosal healing in Crohn disease. *Pol Merkur Lekarski*. 2009;26(155):554-555;
- Baert F et al.: Mucosal healing predicts sustained clinical remission in patients with early-stage Crohn's disease. *Gastroenterology*, 2010;138(2):463-468
- Allez M and Lémann M: Role of endoscopy in predicting the disease course in inflammatory bowel disease. *World J Gastroenterol*. 2010;16:2626-2632
- Lasson A: Calprotectin in feces a well-documented marker of gastrointestinal inflammation. Indicates disease intensity-- normalization of values predict mucosal healing. *Läkartidningen*, 2010;107(143):2645-2649
- Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. *Clinica Chimica Acta* 381 Feb 2007, p. 63-68
- Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. *Inflamm Bowel Dis*, Vol 15 No 5; May 2009, p. 756-759
- Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. *Q J Med* 2005, Vol 98, May 2005, p. 435-441.
- Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol*. 2008; Vol 103, p. 162-169.
- Gisbert, J.P. et al. Role of Biological Markers in Inflammatory Bowel Disease. *Gastroenterol. Hepatol*. May 2007, 30(3), p. 117-29.

FD Darm check professionell

Intended Use

The **FD Darm check professional** rapid test is based on a visual, immunochromatographic method for the qualitative detection of human **hemoglobin**, **hemoglobin/haptoglobin complex** and **calprotectin** in fecal samples. The test is intended for professional in vitro diagnostic use in hospitals, medical laboratories and medical practices.

Medical Importance

The **FD Darm check professional** combines the detection of hemoglobin, hemoglobin/haptoglobin complex with the detection of calprotectin. Blood in the stool does not necessarily have to originate from polyps or adenomas or carcinomas, but can also indicate other pathological and/or inflammatory changes in the intestine.

The hemoglobin and hemoglobin/haptoglobin complex rapid test exclusively detect human hemoglobin or the human hemo/hapto complex by using antibodies. Hemo/haptoglobin complex plays an important role in the recovery of hemoglobin from lysed erythrocytes and is relatively stable against acids and proteolytic degradation. This means that the hb/hp complex is detectable even after longer intestinal passage. Thus, blood accumulations of polyps and carcinomas located higher in the intestine can also be detected. In general, the combined use of hemoglobin and hemoglobin/haptoglobin complex appears to show increased sensitivity for colorectal adenomas, depending on the location, size, and degree of dysplasia.

Calprotectin is a leucocyte marker and in the stool sample a specific marker for a bowel inflammation. Various organic disorders of the gastrointestinal tract may damage the intestinal mucosa. These damage can result in increased permeability of the mucous membrane, inflammation or formation of ulcers. The intestinal content contains many bacteria and other microorganisms, which can release toxic or chemotactic substances. These substances stimulate e.g. Leukocytes, in particular polymorphonuclear granulocytes (PMN), to the intestinal lumen to secrete their contents, among other antimicrobial substances such as calprotectin. As a sensitive marker calprotectin specifically exhibits a daring inflammation and can be remarkable increased in chronic inflammatory diseases (CED), such as ulcerative colitis or Crohn's disease.

FD Darm check professionell

Principle of the Test

The immunological **FD Darm check** professional rapid test is based on an immunochromatographic method in which specific antibodies specifically detect hemoglobin, hemoglobin/haptoglobin complexes and calprotectin. The respective parameter to be detected reacts with the specific antibodies which are labelled with gold particles. This complex is distributed over the membrane and reaches the test band (T). On this band, an antibody against this complex is applied. If this complex is detected and the result is positive, a color change takes place.

In the case of a positive result, the corresponding molecules (hemoglobin hemoglobin/haptoglobin complex respectively calprotectin) from the stool sample bind to the gold labelled antibodies on test band (T) whereby a pink-red coloration of the band becomes visible. If the result is negative, there are no molecules that can bind to the test band (T) and no coloration is shown. The coloration of the control band (C) guarantees that sample application and migration is done correctly and that the test is valid.

Our buffer allows performing **all three tests with the same tube.**

Batch specificity

Standard in-house validations make sure that it is possible to combine stool collecting tubes and test cassettes of different batches. This combination has no erroneous impact on test results. Batches are compatible and interchangeable with each other.

Contents

- 20 individually packed hb and hb/hp double test cassettes
- 20 individually packed calprotectin test cassettes
- 20 sets fecal sampling sets with sample-taking tubes, instructions on taking stool samples and stool collectors
- test instructions

Additional Required Materials

- Watch
- gloves

FD Darm check professionell

Warnings

- Valid for in vitro diagnostic and professional use only.
- Do not use the test after expiration date.
- Only open the test packaging if the test is actually going to be performed.
- All patient samples should be handled as contaminated.
- The extraction solution contains a small amount of sodium acid.
- Patients should follow the instructions for a sample collection exactly and should not take any sample during menstruation or in case of bleeding hemorrhoids. It is recommended to use enclosed stool container in order to avoid contact of stool with blood from urine samples or interfering components from the toilet water.


Storage and stability of reagents

The test can be stored at 2° to 30°C in the sealed bag until expiration date. The expiry date is mentioned on both the packaging of the rapid test cassette and on the label of the tube.

FD Darm check professionell

FD Dripping stool sample tube

Contents:

- (1) Ejecting piston
- (2) Sampling rod
- (3) Intermediate section with wiper
- (4) Dripping Tube
- (5) Dripping tip
- (6) Closing cap


Dosing tip collects 40 mg stool sample

FD Darm check professionell


Specimen Collection and Preparation


1. **Collect** the stool using the stool collector enclosed.


2. Use sample-taking tube and **detach sample rod (2)** (with corrugated cap) from tube by turning the cap.

3. **Stick** white sampling rod (2) into the **stool** and withdrawing it with a turning. Take special care that opening of sampling rod is completely filled with stool.

4. **Put** white sampling rod (2) completely in the intermediate section (3) and close it with a **quarter turn in clockwise direction**. Take care that the tube is completely sealed.


- If the sample is taken at patient's home, the sample should be stored protected from light and in cool place.

FD Darm check professionell


Procedure of the Test


1. First **turn** the yellow ejection piston (1) **clockwise a quarter rotation** and then push it simultaneously down to transfer stool sample from the sampling rod (2) into the extraction buffer.

2. **Shake** stool-collecting tube (ideally on vortexer) until stool specimen has completely **dissolved** into extraction buffer.

3. **Remove** grey closing cap (6).

4. **Discard** the first drop of the sample and add
 - **three drops** (hb and hb/hp cassette) and
 - **three drops** (calprotectin cassette)
 of solution via dripping tip (5) into the round sample window on respective test cassette. First results can be obtained after **8 minutes**. Strongly positive samples may be readable earlier. **Viscous stool samples** may lead to extended running time (up to **15 minutes**).




- The stool sample in buffer solution should be stored cool at 2° - 8°C if possible. The test should be evaluated within 5 days.

FD Darm check professionell


Interpretation of results for Hb and Hb/Hp complex



Positive: 2 pink-colored lines appear; one in the control region (C) and one in the test region (T). A similar intensity of both lines means that the test is positive. A weak staining of (T) indicates a weak positive result. (T) can be strongly colored than (C), by highly positives samples. If the Test is positive, a further examination with different clinical methods is recommended.

Negative: Just two pink-colored lines appear in the control region (C). In the test region (T) no line appears.

Invalid: If no line appears in control region (C), it is a sign that the test is not running properly or that the test materials so not work in the correct way. Too many solid particles in the stool sample can also lead to an invalid result. In this case repeat the test with a new test cassette or contact the manufacturer for technical support.


FD Darm check professionell

Interpretation of results for calprotectin


Positive: A red line appears in control region (C) and a red one in test region (T). A weak staining of (T) indicates a weak positive result. If the Test is positive, a further examination with different clinical methods is recommended.

Negative: Only the red line in control region (C) appears. In test region (T) no line appears.

Invalid: If no line appears in control region (C), it is a sign that the test is not running properly or that the test materials do not work in the correct way. Too many solid particles in the stool sample can also lead to an invalid result. In this case repeat the test with a new test cassette or contact the manufacturer for technical support.

FD Darm check professionell

Restrictions

The test is used for qualitative determination of human hemoglobin and hemoglobin-haptoglobin complex and Calprotectin in stool samples.

Stool samples should not be taken during or 3 days before or after menstruation, in the case of bleeding caused by constipation, bleeding hemorrhoids, or in the case of medicine being administered rectally. False positive results could be obtained.

Blood-thinning medications, such as, acetylsalicylic acid or coumarin, as well as iron preparations can lead to bleeding which is not caused by tumors.

No special diet is necessary before the performance of the test.

A positive test occurs if human hemoglobin or hemoglobin-haptoglobin complex is present in the sample. In addition to intestinal bleeding, blood in stool may also have different origins, such as hemorrhoids, menstruation bleeding or the mixing of the stool sample with urine containing blood. Not all intestinal bleedings are the result of benign or malignant adenoids. The results obtained by performing this test should be verified by using other methods.

Since carcinomas and adenoids can bleed to different extents and also intermittently it is advisable to check several stool samples of the same patient.

Urine and excessive dilution of samples with water from the toilet bowl can lead to false test results. It is therefore recommended to use the stool container included in the test.

Stool samples from patients treated with non-steroidal anti-inflammatory drugs may lead to false-positive results in the calprotectin test.

Negative results cannot exclude an inflammation. Some diseases, such as celiac disease/sprue and microscopic colitis polyps a result mainly of mononuclear inflammation, can react negatively in calprotectin test.

With a positive result of calprotectin the quantitative determination in the laboratory is recommended to get an initial value for a possible following therapy monitoring.

The **FD professionell Darm check** should be evaluated as positive as soon as one of the three parameters (hemoglobin, hemoglobin/haptoglobin complex or calprotectin) is positive. In this case, even if just one parameter is positive, further diagnostic investigations should be initiated, because the health of the intestine may be impaired.

FD Darm check professionell

Performance Characteristics

Hb/Hp complex rapid test cassette

Sensitivity:

The limit of detection of the test are 60ng for free and 25ng for complex hemoglobin per milliliter of buffer solution.

Prozone effect:

The Hb/Hp complex rapid test cassette also works reliably at extremely high hemoglobin values (>500,000ng/mL).

Specificity:

The Hb/Hp complex rapid test cassette is specific for human hemoglobin or the hemoglobin-haptoglobin complex and shows no cross-reactivity with hemoglobin from cattle, pigs, rabbits, horses, or sheep at a concentration of up to 0.5 mg/mL in the extraction buffer. The hemoglobin-haptoglobin complex also displays no cross-reactivity with bilirubin, vitamin C, and horseradish peroxidase.

Calprotectin rapid test cassette

Sensitivity:

Samples having a calprotectin concentration $\geq 50\mu\text{g} / \text{g}$ stool lead to a positive result in caürptectin rapid test.

The test showed in comparison to other commercial immunoassay, for the detection of human calprotectin, a diagnostic sensitivity > 94%.

Specificity:

The calprotectin rapid test showed in comparison to other commercial immunoassay, for the detection of human calprotectin, a diagnostic specificity of 93%.

The test is specific to human fecal calprotectin.

FD Darm check professionell

Literature

- Knittel, T.; Rudi, V.; Lutz, C.; Abels, L.; Frost, G.; High prevalence of elevated fecal calprotectin levels in stool samples from patients with diverticular disease: (2017) *UEG week; Gastroenterologie, May 2018, Volume 154, Issue 6, Supplement 1, Page S664*
- Knittel, Th., Lüthgens, K.; Higher sensitivity for the detection of adenomas by a Hb/Hp bed site test using a novel sampling device (2013)
- Tao, S., Brenner, H.; Well adjusted qualitative immunochemical faecal occult blood tests could be a promising alternative for inexpensive, high-quality colorectal cancer screening; *European Journal of Cancer Prevention* 2013, 22:305-310 Sieg A., Thoms C., Lüthgens K., John M.R., Schmid Gayk H.; Detection of colorectal neoplasm by the highly sensitive hemoglobin-haptoglobin complex in feces; *Int J Colorectal Dis* (1999) 14: 267-271
- Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; *Archive of Internal Medicine*; (1995) 155:2389-2402
- Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; *British Medical Journal*; (1988) 296:1092-1094
- Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? *Cancer Supplement*; (1994) 74(7) : 2023-2027
- Miller, A.B.; An Epidemiological Perspective on Cancer Screening; *Clinical Biochemistry* (1995) 28(1) : 41-48
- Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
- Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; *Morbidity and Mortality Weekly Report*; (1995) 45 (5) : 107- 110
- St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; *Gastroenterology*; (1993) 104:1661-1668
- Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: *Hepatogastroenterology*; 2000 Mar-Apr; 47(32):396-9.
- and *Mortality Weekly Report*; (1995) 45 (5) : 107-110
- Hundt S, Haug U, Brenner H; *American College of Physicians* 2009; 150(3):162-169
- Devlin SM and Panaccione R: Evolving inflammatory bowel disease treatment paradigms: top-down versus step-up. *Med Clin North Am.* 2010;94:1-18
- Pineton de Chambrun G et al.: Clinical implications of mucosal healing for the management of IBD. *Nat Rev Gastroenterol Hepatol* 2010; 7(1):15-29
- Lichtenstein GR and Rutgeerts P: Importance of mucosal healing in ulcerative colitis. *Inflamm Bowel Dis.* 2010;16:338-346
- Smith MA et al.: Pharmacogenomics in the treatment of inflammatory bowel disease. *Pharmacogenetics*, 2010;11(3):421-437
- Lin MV et al.: What is the optimal therapy for Crohn's disease: step-up or top-down? *Expert Rev Gastroenterol Hepatol.* 2010;4(2):167-180
- Strauch U and Schölmerich J.: Emerging drugs to treat Crohn's disease. *Expert Opin Emerg Drugs*, 2010;15(2):309-322
- Isaacs KL: How rapidly should remission be achieved? *Dig Dis* 2010;28(3):548-555
- Schwartz M and Regueiro M: Prevention and treatment of postoperative Crohn's disease recurrence: an update for a new decade. *Curr Gastroenterol Rep.* 2011 Feb;13(1):95-100
- Ha C and Kornbluth A: Mucosal healing in inflammatory bowel disease: where do we stand? *Curr Gastroenterol Rep.* 2010;12(6):471-478.
- Fagerberg UL et al.: Fecal calprotectin: a quantitative marker of colonic inflammation in children with inflammatory bowel disease. *J Pediatr Gastroenterol Nutr.* 2007;45(4):414-420
- Rutgeerts P et al.: Biological therapies for inflammatory bowel diseases. *Gastroenterology*, 2009;136(5):1182-1197
- Jalocha L et al.: Mucosal healing in Crohn disease. *Pol Merkur Lekarski.* 2009;26(155):554-555;
- Baert F et al.: Mucosal healing predicts sustained clinical remission in patients with early-stage Crohn's disease. *Gastroenterology*, 2010;138(2):463-468
- Allez M and Lémann M: Role of endoscopy in predicting the disease course in inflammatory bowel disease. *World J Gastroenterol.* 2010;16:2626-2632
- Lasson A: Calprotectin in feces a well-documented marker of gastrointestinal inflammation. Indicates disease intensity-- normalization of values predict mucosal healing. *Läkartidningen*, 2010;107(143):2645-2649
- Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. *Clinica Chimica Acta* 381 Feb 2007, p. 63-68
- Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. *Inflamm Bowel Dis*, Vol 15 No 5; May 2009, p. 756-759
- Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. *Q J Med* 2005, Vol 98, May 2005, p. 435-441.
- Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol.* 2008; Vol 103, p. 162-169.
- Gisbert, J.P. et al. Role of Biological Markers in Inflammatory Bowel Disease. *Gastroenterol. Hepatol.* May 2007, 30(3), p. 117-29.

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SYMBOL EXPLANATIONS

	Only for in-vitro diagnostic use		Only for single use
	Content		Expiry date
	Lot number		Storage temperature

FROST DIAGNOSTIKA GmbH

Speyerer Straße 74

D-67166 Otterstadt

Telefon: +49 (0) 6232 600487 0

Telefax: +49 (0) 6232 600487 60

www.frostdiagnostika.de

eMail: info@frostdiagnostika.de