

**DIAGNOSTIC KIT
FOR DETERMINATION
OF FIBRINOGEN CONCENTRATION**

Cat. No. GLAU964



AU Series – FIBRINOGEN

INTRODUCTION

Fibrinogen is the plasma protein precursor of fibrin, which when cross-linked becomes the principal component of fibrin clot.

Fibrinogen is a sensitive acute phase protein whose concentration rises several fold during inflammation. It can be increased as a response to inflammatory states, with infections, during pregnancy and after trauma.

Evidence has shown that plasma levels above the reference range constitute a significant independent risk factor for both coronary artery and cerebrovascular diseases.

Fibrinogen can be deficient in congenital afibrinogenemia. Levels may also fall significantly as the result of sequestration in extravascular spaces such as the body cavities and sites of recent trauma.

METHOD PRINCIPLE

Fibrinogen presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to fibrinogen concentration in the sample.

REAGENTS

Package	
1-Reagent	2 x 20.0 ml
2-Reagent	2 x 5.0 ml

Buffer (1-Reagent) stored at 2-8°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package. Store closed and avoid contamination.

Concentrations in the test

TRIS (pH 8.0); PEG; sodium chloride; anti human fibrinogen antiserum; HEPES buffer (pH 7.4); sodium azide (< 1 g/l); stabilizers.

Warnings and notes

- Products for in vitro diagnostic use only.
- The reagents must be used only for the intended purpose, by suitably qualified laboratory personnel, under appropriate laboratory conditions.
- Products from human source have been tested for the HIV antibody, HbsAg and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
- Products contain sodium azide (< 1 g/l) as a preservative. Avoid contact with skin and mucous membranes.

SPECIMEN

Plasma (sodium citrate).
Before analysis calibrator and samples should be diluted 20-fold with 0.9% NaCl and mixed gently.
It is recommended to perform the assay with freshly collected samples!

PROCEDURE

These reagents may be used in automatic analysers Olympus AU400/AU640.
1-Reagent and 2-Reagent are ready to use.
For reagent blank 0.9% NaCl is recommended.

APPLICATION

Reagent ID: 505

General		LIH	ISE	Range
Test name:	FIBR			
Type:	Serum			
Operation:	Yes			
Sample: Volume	5 μ L	Dilution	0 μ L	Pre-Dilution Rate: 1
Reagents: R1 Volume	250 μ L	Dilution	0 μ L	Min OD
R2 Volume	50 μ L	Dilution	0 μ L	Max OD
			L	-2.0000
			H	2.5000
Wavelength:	Pri. 340	Sec. 700	Reagent OD Limit:	
Method:	END		First L	-2.0000
Reaction Slope:	+		First H	2.5000
Measuring Point 1: First	0	Last 27	Last L	-2.0000
Measuring Point 2: First	0	Last 10	Last H	2.5000
Linearity:			Dynamic Range:	
No-Lag-Time:			Correlation Factor:	
			A	1.000
			B	0.000
			On-board Stability Period:	

General		LIH	ISE	Range			
Test name:	FIBR						
Type:	Serum						
Value/Flag:	#	Level L:	#	Level H:	#		
Normal Ranges:							
	Sex	Age L Year	Month	Age H Year	Month	L	H
1.	#	#	#	#	#	#	#
2.	#	#	#	#	#	#	#
3.	#	#	#	#	#	#	#
4.	#	#	#	#	#	#	#
5.	#	#	#	#	#	#	#
6.	#	#	#	#	#	#	#
7.	None Selected					#	#
8.	Out of Range					#	#
Panic Value:		L	#	H	#	Unit:	mg/dl
						Decimal Places:	1

General		ISE		
Test name:	FIBR			
Type:	Serum			
Calibration Type:	2AB	Formula: Polygonal		
Counts:	3	Process: CONC		
Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H
Point 1:	#	**	-2.0000	2.5000
Point 2:	#	*	-2.0000	2.5000
Point 3:				
Point 4:				
Point 5:				
Point 6:				
Point 7:				
1-Point Cal.Point:		with CONC=0	Slope Check: None	Advanced Calibration: #
MB Type Factor:		Calibration Stability Period:		

- # User defined
- * Calibrator value
- ** Saline should be used as calibrator 1

REFERENCE VALUES⁴

adults	200 – 400 mg/dl
children	125 – 300 mg/dl

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the control plasma with each batch of samples.
For the calibration of automatic analysers systems the FIBRINOGEN CALIBRATOR (Cat. No GL9964) is recommended.
The calibration curve should be prepared every 4 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analyser AU640. Results may vary if a different instrument is used.

▪ **Analytical range:** 20 mg/dl to 1000 mg/dl.

▪ **Interferences:**

Hemoglobin up to 0.32 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 2000 mg/dl, sodium citrate up to 5 g/l do not interfere with the test.

▪ **Precision**

Repeatability (run to run) n = 10	Mean [mg/dl]	SD	CV [%]
level 1	274.4	3.9	1.4
level 2	146.3	3.1	2.1
level 3	104.4	5.2	5.0

Reproducibility (day to day) n = 10	Mean [mg/dl]	SD	CV [%]
level 1	269.7	7.5	2.8
level 2	141.9	5.8	4.1
level 3	98.1	5.0	5.1

▪ **Method comparison**

A comparison between the Glenbio reagent (y) and commercially available assay (x) using 17 samples gave following results:

$$y = 0.86x + 86.9 \text{ mg/dl};$$

$$R = 0.9556 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Bergstrom, K. & Lefvert, A.K. Scand.J.clin.Lab.Invest. 40, (1980), 637.
2. Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, (1994).
3. Roitt, I., Essential Immunology, Blackwell, Oxford, (1991).
4. Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: Moss D. W., Henderson A. R., 1813, (1999).



GLENBIO LTD
10 Kilbegs Road, Antrim, Co. Antrim, BT41 4NN
Tel/Fax: +44(0)2879659842
Email: info@glenbio.com
Web: www.glenbio.com



GLENBIO IRELAND LTD
17b Fota Business Park, Carrigtwohill, Co. Cork,
T45 PK77, Ireland

OLYMPUS AU SERIES ARE REGISTERED
TRADEMARKS OF OLYMPUS CORPORATION,
LAKE SUCCESS, NEW YORK, USA

