NovaLisa[®]

Free T4

ELISA

CE

Only for in-vitro diagnostic use

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ENGLISH

1. INTENDED USE

The Free T4 ELISA is intended for the quantitative determination of free T4 (Thyroxine) in human serum or plasma (citrate, heparin).

2. PRINCIPLE OF THE ASSAY

The quantitative immunoenzymatic determination of free T4 is based on the competitive ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific anti-T4 antibodies. Free T4 in the sample competes with added horseradish peroxidase (HRP) conjugated T4 for antibody binding. Subsequently, the wells are washed to remove all unbound material. The immune complexes are visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is inversely proportional to the amount of free T4 in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

3. MATERIALS

3.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with anti-T4 antibodies; in resealable aluminium foil.
- **SOLN STOP:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **WASH BUF 20x:** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap; 0.2% (w/v) 5-Bromo-5-nitro-1,3-dioxane.
- **Conjugate:** 1 bottle containing 15 mL of peroxidase labelled T4; coloured red; ready to use; black cap; ≤ 0.02% (v/v) MIT.
- **SUB TMB:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Control:** 1 vial containing 1 mL; brownish colour; ready to use; white cap; ≤ 0.01% (v/v) CMIT/ MIT (3:1).
- **Standards:** 6 vials, each containing 1 mL standard; brownish colour; ready to use; yellow cap; ≤ 0.01% (v/v) CMIT/ MIT (3:1).

Standard A:	0 ng/L
Standard B:	4 ng/L
Standard C:	10 ng/L
Standard D:	20 ng/L
Standard E:	35 ng/L
Standard F:	75 ng/L

For hazard and precautionary statements see 11.1.

3.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

3.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes of 50 and 100 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

4. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C. The standards can be kept for a maximum of 6 months after opening.

5. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

5.1. Microtiterplate

The break-apart snap-off strips are coated with anti-T4 antibodies. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

5.2. **WASH BUF 20x**

Dilute **WASH BUF 20x** 1 + 19; e. g. 10 mL **WASH BUF 20x** + 190 mL distilled water. The diluted buffer (**WASH BUF 1x**) is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g., in a water bath. Mix well before dilution.

5.3. Standards

The standards are ready to use. They are stable until the expiry date stated on the label. Once opened, they can be kept for 6 months if stored at 2...8 °C; but no longer than the expiry date stated on the label.

5.4. **SUB TMB**

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. **SUB TMB** should be colourless or could have a slight blue tinge. If **SUB TMB** turns into blue, it may have become contaminated and should be thrown away.

6. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise, they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

7. ASSAY PROCEDURE

Please read the instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the volume of **WASH BUF 1x** from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 11. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 50 µL standards/controls and samples into their respective wells.
2. Dispense 100 µL conjugate into all wells.
3. Cover wells with the foil supplied in the kit.
4. **Incubate for 1 hour ± 5 min at 37 ± 1 °C.** Do not expose to direct sunlight.
5. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well five times with 300 µL of **WASH BUF 1x**. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!

Note: Washing is important! Insufficient washing results in poor precision and false results.

6. Dispense 100 µL **SUB TMB** into all wells.
7. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
8. Dispense 100 µL **SOLN STOP** into all wells in the same order and at the same rate as for the **SUB TMB**, thereby a colour change from blue to yellow occurs.
9. Measure the absorbance at 450/620 nm within 30 min after addition of the **SOLN STOP**.

7.1. Measurement

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

8. RESULTS

8.1. Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Standard A:** Absorbance value > **1.000** and > **Standard B**
- **Standard B:** Absorbance value > **Standard C**
- **Standard C:** Absorbance value > **Standard D**
- **Standard D:** Absorbance value > **Standard E**
- **Standard E:** Absorbance value > **Standard F**
- **Standard F:** Absorbance value < **0.200**
- **Control:** Result in ng/L within the range indicated on the label.


Standard A > Standard B > Standard C > Standard D > Standard E > Standard F

If these criteria are not met, the test is not valid and must be repeated.

8.2. Calculation of Results

In order to obtain **quantitative results in ng/L** plot the (mean) absorbance values of the Standards A - F (linear y-axis) against their corresponding concentrations (logarithmic x-axis) and draw the best-fit curve (e.g., Four Parameter Logistic) through the plotted points. Read results from this standard curve employing the (mean) absorbance values of each sample and control.

8.3. Typical Standard Curve

8.4. Interpretation of Results

Normal value ranges for this ELISA should be established by each laboratory based on its own patient populations in the geographical areas serviced.

Test results should be verified concerning the entire clinical status of the patient. Also, every decision for therapy should be taken individually.

Diagnosis of a disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as laboratory findings.

8.4.1. Normal Range

A total of 173 samples from apparently healthy adults aged 18 - 68 years were assayed to establish the reference range interval for the Free T4 ELISA according to approved guideline CLSI EP28-A3c. The reference interval for the Free T4 ELISA was found to be **7 – 22 ng/L** based on the central 95 % of the frequency distribution. The mean value, standard deviation and expected range are:

	Mean (ng/L)	SD	Range (ng/L)
Adults (n=173)	11.10	5.27	7-22

8.4.2. Samples from pregnant women

Samples from apparently healthy pregnant women were used to determine expected free T4 values for the Free T4 ELISA during pregnancy. The mean value, standard deviation and expected range are:

	Mean (ng/L)	SD	Range (ng/L)
Pregnancy (n=30)	10.03	2.19	6-16

9. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

9.1. Precision

<u>Intraassay</u>	<u>n</u>	<u>Mean (E)</u>	<u>CV (%)</u>
# 1	24	2.086	2.69
# 2	24	1.337	5.47
# 3	24	0.224	6.28
<u>Interassay</u>	<u>n</u>	<u>Mean (ng/L)</u>	<u>CV (%)</u>
# 1	12	11.18	6.68
# 2	12	24.29	4.55
# 3	12	36.96	3.99

9.2. Analytical Sensitivity

The Limit of Detection (LoD) for the Free T4 ELISA is 1.06 ng/L, determined consistent with the guidelines in CLSI document EP17-A2 based on the proportions of false positives (α) less than 5 % and false negatives (β) less than 5 %; using 240 determinations, with 120 analyte-free samples and 120 low level samples; and a Limit of Blank (LoB) of 0.00 ng/L.

9.3. Analytical Specificity

9.3.1. Cross-Reactivity

The following substances, when tested in samples at free T4 concentrations of approx. 10 ng/L and 25 ng/L according to CLSI protocol EP07-A2, caused no cross-reactivity up to the concentrations indicated in the table below:

Substance Tested	Test Concentration [ng/L]	Cross-reactivity [%]
3,3',5-Triiodo-L-thyronine (T3)	10 000	< 0.060
3,3',5'-Triiodo-L-thyronine (Reverse T3)	100 000	< 0.001
3,5-Diiodo-L-thyronine (3,5-T2)	50 000	< 0.001
3,3',5-Triiodothyroacetic Acid	10 000	< 0.010
3,3',5,5'-Tetraiodothyroacetic Acid	100 000	< 0.001

9.3.2. Interferences

The following substances, when tested in samples at free T4 concentrations of approx. 10 ng/L and 25 ng/L according to CLSI protocol EP07-A2, were found not to interfere at the concentrations indicated. A bias less than 10 %* is not considered a significant interference.

Substance Tested	Test Concentration
Bilirubin	200 µg/mL
Hemoglobin	2 mg/mL
Triglycerides	35 mg/mL
Albumin	63 mg/mL

* upper limit of 95 % confidence interval

Investigation of a rheumatoid factor positive specimen panel did not reveal significant evidence of falsely elevated or decreased free T4 levels due to interference. However, in free-hormone immunoassays interference caused by rheumatoid factors should be considered.

9.4. Correlation

The Free T4 ELISA was compared with another immunoassay. 130 specimens were used (the values ranged from approx. 1 ng/L to approx. 40 ng/L). The linear regression curve was calculated.

$$y = 0.973 + 0.145x$$

$$R^2 = 0.917$$

9.5. Measurement Range

The measurement range is 1.06 ng/L – 75 ng/L.

10. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

11. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

11.1. Safety note for reagents containing hazardous substances

Reagents may contain CMIT/MIT (3:1) or MIT (refer to 3.1).
Therefore, the following hazard and precautionary statements apply.

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing spray.
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362+P364	Take off contaminated clothing and wash it before reuse.

Reagents may contain 5-Bromo-5-nitro-1,3-dioxane (refer to 3.1).
Therefore, the following hazard and precautionary statements apply.

Warning

H315	Causes skin irritation.
H319	Causes serious eye irritation
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice/attention.

Further information can be found in the safety data sheet.

11.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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Free T4

Test Preparation

Prepare reagents and samples as described.
 Establish the distribution and identification plan for all samples and standards/controls.
 Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Standard A - F	Control	Sample
Standard A - F	50 µL	-	-
Control	-	50 µL	-
Sample	-	-	50 µL
Conjugate	100 µL	100 µL	100 µL
Cover wells with foil supplied in the kit Incubate for 1 h at 37 ± 1 °C Do not expose to direct sunlight Wash each well five times with 300 µL of WASH BUF 1x			
SUB TMB	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature (20...25 °C) in the dark			
SOLN STOP	100 µL	100 µL	100 µL
Photometric measurement at 450 nm (reference wavelength: 620 nm)			

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

FT41050_IFU_rev01_fromLot_031N