

- de** Gebrauchsanweisung
- en** User's Manual
- fr** Mode d'emploi
- it** Istruzioni per l'uso
- es** Instrucciones de empleo
- pt** Manual de operação
- ru** Руководство по применению

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QM-1-034J

<b>de</b>	Kehlkopfspiegel .....	4
<b>en</b>	Laryngeal mirrors .....	8
<b>fr</b>	Miroirs laryngiens .....	12
<b>it</b>	Specchi laringei .....	16
<b>es</b>	Espéculos laríngeos .....	20
<b>pt</b>	Espelhos laríngeos .....	24
<b>ru</b>	Ларингеальное зеркало .....	28

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Dear Customer, thank you for choosing a KaWe product. Our products are known for their high quality and longevity. This KaWe product meets EC Standards 93/42/EEC (standards for medical products).

**Please read these instructions thoroughly and carefully before attempting to use this product and heed the given care instructions.**

**Use:** These laryngeal mirrors are only to be used by authorized, trained personnel. Laryngeal mirrors are instruments used to inspect the larynx (indirect laryngoscopy).

**Use of the instrument:** During a laryngoscopy procedure, the tongue is carefully brought forward and a flat mirror is then inserted up to the uvula. For this purpose the laryngeal mirror must be inserted into the provided opening on the instrument handle and tightened with the fastening screw. It is important to ensure that the direction of the reflected image is correct.

**Unsuited use/contraindication of the product:** Any use other than that described here is not in accordance with the intended use of the instrument. The manufacturer is not liable for any resulting damages. The user alone bears the risk.

**Further information, maintenance, storage:** When used and stored properly, this product will serve you reliably for many years.

**Guarantee:** When used under normal circumstances and as directed in our User's Manual, we guarantee this product for two years after the date of purchase. Should you need further information or should your instrument require repair, please contact your dealer!

**Warning:** Ensure that the laryngoscope has been firmly screwed to the handle in order to avoid its unintended detachment! When cleaning reusable instruments special attention is required. Do not exceed 134°C. Before the first sterilisation of the laryngoscope, any residue that may be present (such as oil or grease) must first be removed.

**Limitation of reuse preparation:** Frequent re-preparation has only a minimal effect on this product. The end of the product's useful life is normally due to wear and tear and damage by use. It is recommended that instruments be re-prepared as soon as possible after their use.

### Instructions:

**Storage:** Do not store in humid rooms (such as bathrooms) or transport along with wet objects. Protect from prolonged exposure to the sun during storage and transport. Store sterilised laryngoscopes and handles in a dry, dust-free location that is protected from contaminants.

**Cleaning preparation:** Clean and disinfect the instrument immediately after using it. A mechanical cleaning method (instrument washing machine) should always be used. Prepare the cleaning and disinfection treatment as described in the chapter "Instrument Preparation". Any suitable cleaning or disinfecting agent can be used that is listed by DGHM.

**Manual cleaning:** In accordance with the recommendations given by the Robert Koch Institute (RKI), a machine is to be used for further preparation of the instrument. Manual preparation is not recommended.

### Automatic cleaning: equipment:

CM 310 (Maquet), cleaning agent (neodisher®FA forte 0.4%/neodisher®Z 0.2%)

G7828 (Miele), cleaning agent (Mucapur®XL 0.4 %/Mucapur®Z 0.15 %)

WD 390 (Belimed), cleaning agent (Mucapur®AF 0.5 %/Mucapur®Z 0.1 %)

**1.** Rinse the instruments thoroughly with running water immediately before placing them in the washing machine in order to prevent residual cleaning/disinfecting agents from entering the machine. **2.** The instruments are to be placed in a suitable instrument stand. **3.** Position the instrument stand in the washing machine so that the instruments are directly sprayed by the jet. **4.** Put the cleaning agent into the machine as directed. **5.** Start des Vario TD program with the thermal disinfection function. The thermal disinfection process is dependent on the AO value and national regulations (EN/ISO 15883). **6.** After the program has finished, remove the instruments from the RDG and dry them (preferably with pressurised air as recommended by the RKI). Ensure that difficult-to-reach areas are dried when using an instrument stand. **7.** Check for damage and cleanliness with a magnifying glass. (Experience has shown

that a magnifying glass with 8-times magnification is optimal for proper visual inspection.) If residual contaminants are visible on the instrument after being machine washed, repeat the cleaning and disinfecting process until no more contaminants can be detected. **8. Caution!** If the instruments are only mechanically cleaned (without demonstrable disinfection), a final thermal disinfection in a steam sterilizer in suitable racks or trays is required.

**Disinfection:** The effectiveness of the used cleaning and disinfecting method must be recognised (such as in the list published by Robert Koch Institute/DGHM of tested and approved disinfection agents and methods) and approved. For agents such as Helipur®, please follow the directions given by the manufacturer regarding the concentrations and reaction times.

**Maintenance:** Remove and replace damaged parts. Surgical grease/lubricant can be used on all movable parts such as the fastening screw.

**Inspection and function test:** Check moving parts such as the fastening screw for ease of movement. Inspect the mirror for cracks or breakage. All instruments: Visually inspect for damage and wear.

**Packaging:** individual/non-sterile: A standard polyethylene bag can be used. The bag used must be large enough to accommodate the instrument such that the bag closure is not strained. Sterile: The sterile packaging must comply with DIN EN 868/ISO 11607 standards and be suitable for steam sterilisation.

**Sterilisation:** The laryngoscope and handles can be effectively sterilized with the following sterilising parameters: fractionated pre-vacuum (3-fold), sterilisation temperature: min. 132 °C, max. 137 °C, retention time: at least 3 min (complete cycle), drying time: 10 min.

**Validated reprocessing procedure:** The effectiveness of this procedure was proven using the disinfectant Desinfektor GM 310, which is made by Maquet, and the cleaning agent Neodisher® FA forte with Neodisher® Z 0.2% as a neutralizer, which is made by the company Chemische Fabrik Dr. Weigert, Hamburg.

**Disposal:** This product is to be **sterilised** and then disposed of in accordance with local regulations at a separate waste collection site/recycling centre.

**Additional information:** When sterilizing several instruments at the same time in the autoclave, ensure that the maximum capacity specified by the manufacturer of the sterilizer is not exceeded.

**Common accessories:** Further information about these articles can be found on our homepage: [www.kawemed.de](http://www.kawemed.de).

Manufacturer: **KaWe**

**Contacting the manufacturer:** Write to or phone your dealer or dial: +49-7141-68188-0.

#### Symbol key:

	Manufacturer		Caution!
	Date of manufacture		Complies with relevant EU guidelines
	Heed the User's Manual		Lot code
	GOST-R certification of exports to Russia		

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**KIRCHNER & WILHELM** GmbH + Co. KG