

**GenVInSet**  
**MTHFR A1298C**

# Instructions for Use

Kit for detection of MTHFR A1298C  
polymorphism

For In Vitro diagnostic use

Product code GVS-A1298C-24 (24 tests)  
GVS-A1298C-48 (48 tests)

Store from  $-18$  to  $-30^{\circ}\text{C}$

Rev 01 / 2020-11-23



Blackhills Diagnostic Resources, S.L.U.  
Camino del Pilón 86, Casa 7 Local. 50011- Zaragoza - Spain  
[www.bdrdiagnostics.com](http://www.bdrdiagnostics.com)

# Index

Intended use	3
Summary and explanation	4
Procedure principles	5
Kit contents	6
Kit storage	7
Materials required but not supplied	8
Sample collection and preparation	9
Usage procedures	10
A) PCR setup	
B) Thermal cycler configuration	
Results	13
Quality control	15
Specific operation data	16
Procedure limitations	17
Troubleshooting guide	18
References	20
Notice to purchaser	22
Explanation of symbols used on the labels	23

**GenVinSet**

## **MTHFR A1298C**

## **Intended use**

GENVINSET® MTHFR A1298C is a kit for the A->C transition in the 1298 position, located in the exon 4 of the MTHFR gen, using Real Time PCR technology with specific TaqMan® probes.

3



Camino del Pilón 86, Casa 7 Local  
50011. Zaragoza - Spain  
[www.bdrdiagnostics.com](http://www.bdrdiagnostics.com)

## Summary and explanation

Methylenetetrahydrofolate reductase (MTHFR) is a key enzyme in the folate metabolism. It catalyzes the intracellular conversion of 5,10-methylenetetrahydrofolate (5,10-methylene THF) to 5-methyltetrahydrofolate (5-methyl THF), the predominant circulatory form of folate and primary methyl donor for the conversion of homocysteine to methionine. Although the MTHFR enzyme does not participate in the clotting cascade as other proteins do, such as Factor II (FII) and Factor V (FV) (which belong to the called "coagulation factor" group), its involvement in the folate metabolic pathway may induce the development of thrombophilia, as well as increase the risk of Alzheimer's disease (1-4).

The mechanism through the MTHFR protein is involved in thrombophilia development is not completely understood. It has been suggested that the inhibition of the folate cycle due to the lack of precursors leads to the accumulation of homocysteine in blood.

The accumulation of this metabolite can cause a damage in the endothelial cells, which constitute the inner layer of blood vessels, thus triggering the formation of clots. These clots might then travel through circulation to other organs, which could lead to severe complications, such as a heart attack or a pulmonary embolism (5,6).

The MTHFR gene is located on chromosome 1p36.3 and is composed of 2.0 kb, containing 11 exons. The product of the MTHFR gene is a 77 kDa protein (EC 1.5.1.20) (1,7). It have been described approximately 60 polymorphisms, as well as 41 rare yet deleterious mutations in the MTHFR gene (8). To date, the most studied and most common functional variant is the MTHFR C677T polymorphism (rs1801133) but, also, the MTHFR A1298C mutation has been studied.

The hyperhomocysteinemia caused by the MTHFR A1298C polymorphism in the MTHFR enzyme can be related to thrombotic events, although there is not a consensus regarding its role in the pathogenesis of thrombosis yet. The frequency of the MTHFR A1298C variant differs among the different ethnicities. The prevalence in Asian populations is relatively low (18%) but it can be higher in populations with a high level of miscegenation as European population (30%) or Latin American one (20%) (14-15).

## Procedure principles

The detection method used by GENVINSET® MTHFR A1298C is based on the Real Time PCR technology, using TaqMan® probes that specifically anneal to position 1298 of the exon 4 of the MTHFR gene, monitoring the presence of A and/or C nucleotides.

This technique provides high resolution, high sensitivity, specificity and reproducibility.

(\*) See Section "Procedure Limitations".

## MTHFR A1298C

# Kit contents

### GVS-A1298C-24 (24 tests)

- GVS-A1298C-PM: 1 vial x 110 uL Primer Mix (PM)
- GVS-A1298C-MM: 1 vial x 138 uL Master Mix (MM)
- GVS-A1298C-C1: 1 vial x 5 uL wt/wt control (C1)
- GVS-A1298C-C2: 1 vial x 5 uL mut/mut control (C2)
- GVS-RB: 1 vial x 100 uL Reaction Blank (RB)

### GVS-A1298C-48 (48 tests)

- GVS-A1298C-PM: 2 vials x 110 uL Primer Mix (PM)
- GVS-A1298C-MM: 2 vials x 138 uL Master Mix (MM)
- GVS-A1298C-C1: 1 vial x 5 uL wt/wt control (C1)
- GVS-A1298C-C2: 1 vial x 5 uL mut/mut control (C2)
- GVS-RB: 1 vial x 100 uL Reaction Blank (RB)

## **MTHFR A1298C**

### **Kit storage**

In order to ensure a proper performance, reagents should be stored from -18°C to -30°C until their expiration date, indicated on the label of the vial. Do not perform more than 3 freeze/thaw cycles to the Primer Mix (GVS-A1298C-PM) and Master Mix (GVS-A1298C-MM) vials, as this could decrease the sensitivity of the assay and impair results.

Due to the photosensitivity nature of the reagents, avoid continuous exposure to light.

# Materials required but not supplied

## General

- Gloves
- Lab coat

## Consumables

- Filter tips (P1000, P200 & P10)
- 1.5 ml autoclaved tubes
- q-PCR instrument specific reagents (in the case of using RotorGene Q, only 0.1 ml tubes are allowed)

## Equipment

- q-PCR instrument. The following devices have been validated:
  - StepOne™, Applied Biosystems™
  - 7500 Real-Time PCR System, Applied Biosystems™
  - LightCycler® 96 System, Roche
  - Mic Real Time PCR Cycler, Biomolecular Systems
  - DTLite Real-Time PCR System, DNA-Technology
- Vortex mixer
- Pipettes (P1000, P200 & P10)



## Sample collection and preparation

The present test should only be performed with complete blood samples treated with EDTA anti coagulation agents or citrate. Heparin can interfere with the PCR process and should not be used in this procedure.

The technique is compatible with several DNA extraction methods. Before delivering results with a diagnostic purpose, a validation assay with such extraction method should be done.

### **i** Caution

All biological and blood samples should be treated as possibly infectious. When manipulating them, the corresponding basic (universal) precautions should be taken. Samples should always be handled wearing the appropriate personal protection equipment.

# Usage procedures

## A) PCR preparation

### **i** Precautions

- Thaw all of the kit components before starting the assay, mix and centrifuge them.
- Work on ice or over a cool block.
- The PCR should be setup in the pre-PCR area.
- Use only filtered tips and 1.5 ml autoclaved tubes.
- Use gloves and lab coat at all times.
- In each test performed it is recommended to use one control sample for each A/C variant.

10

**1.** Thaw the samples. Prepare a mix with the Master Mix and the Primer Mix for n+1 samples:

	Vol. per sample (µL)
Master Mix	5
Primer Mix	4

**2.** Pipette 9 µL of this mix into the PCR tubes and add 1µL of DNA or Reaction Blank, in case of the contamination control well.

**3.** Seal the plate with convenient sealer, and centrifuge 1 min at 360 xg to ensure that all the volume settles to the bottom of the tube.

**4.** Place the plate in the thermal cycler and start the following program.

B) Thermal cycler configuration

1. Set up the following amplification program:

	Cycle number	Temperature (°C)	Time (mm:ss)	Ramp (%)	Analysis
Denaturation	1	95	05:00	100	X
Cycles	50	95	00:15	100	X
		62	01:00	100	Single
Cooling	1	15	∞	100	X

2. Set up the reading channels.

The emitted fluorescence must be read in FAM (495-520 nm) and HEX (535-554 nm) channels.

**NOTE – Special settings for Rotor Gene Q:**

- a. Open the Rotor-Gene Q – Pure Detection software. Select the tab “Advanced” in the window “New Run”, and click “New”.
- b. Select the type of rotor used (only 0.1 ml tubes accepted, see section 6). Select the “Locking Ring Attached” box and continue by clicking “Next”.
- c. Set the “Reaction Volume” as 10 µl, and identify the operator and the samples.
- d. Click on “Edit Profile” and set up the amplification program (see section 8.C.1). Select the step 60 sec at 62 °C, and click on “Acquiring to Cycling A”. Set “Green” and “Yellow” as the fluorescence acquisition channels. Press “OK”. Click on “OK” to accept and close the “Edit Profile” window.
- e. Click “Gain Optimisation” in the “Run New Wizard” dialog box to open the “Auto-Gain Optimization Setup” window. Scroll down the “Channel Settings” menu and select “Acquiring Channels”. Then click on “Add”. In the window “Auto-Gain Optimisation Channel Settings”, set the following parameters for each channel (“Green” and “Yellow”):

# MTHFR A1298C

Tube position = 1

Target Sample Range: 5 FI up to 10 FI

Acceptable Gain Range: -10 to 10

- f. Check the box "Perform Optimisation Before 1st Acquisition", and click on "Close".
- g. Select "Next" and then "Start Run" in the "New Run Wizard" window.

# Results

GenVInSet® MTHFR A1298C kit constitutes a qualitative technique to detect the presence of an A and/or a C on the 1298 position of the MTHFR gene (exon 7).

It is not necessary to select any passive reference.

Using the present kit, the following results can be obtained:

## Detection of A at 1298 position

On the VIC/HEX channel, the following amplification plot can be observed:




Figure 1. Positive and negative samples in the VIC/HEX channel.

Those samples that report an amplification curve can be considered as positive for A allele and they are identified by a numeric value called "Crossing Point" (Cp). This value corresponds to the cycle in which fluorescence is detected and, thereof, the amplification can be considered as positive.

## Detection of C at 1298 position

On the FAM channel, the following amplification plot can be observed:





Figure 2. Positive and negative samples in the FAM channel.

# MTHFR A1298C

Those samples that report an amplification curve can be considered as positive for C allele.

## Genotyping analysis

In “Genotyping” or “Allelic Discrimination” analysis types, select FAM channel (mutation) on the Y axis, and HEX channel (wildtype) on the X axis. Results will appear similarly as in Figure 2, in which each dot consists of an x and y component. The different kind of samples, wt/wt, wt/mut and mut/mut, will be distributed into 3 groups within the plot. The Reaction Blank will be placed near the origin of coordinates (0,0).



14

Figure 3. Plot showing one heterozygous sample (wt/mut), one normal sample (wt/wt), one homozygous mutated sample (mut/mut) for A1298C mutation and the reaction blank (RB) using Genvinset® MTHFR A1298C kit.

## Quality control

Due to the qualitative nature of this test, it will not be necessary to perform a calibration.

The following criteria have to be taken into account in order to validate an assay:

- The Reaction Blank should report negative results in both FAM and VIC/HEX channels. An amplification curve with a  $C_p > 35$  value should be considered as negative. A  $C_p < 35$  value is associated with a contamination in the Reaction Blank and, therefore, all results from the experiment should be discarded.
- An heterozygous sample (A/C) should report positive results in both FAM and VIC/HEX channels.

The assay must be performed according to the kit recommendations, as well as other quality control procedures that comply with local, federal and/or certifying agencies specifications.

# Specific operation data

## 1. Analytical specificity

The probes align specifically on the 1298 position at the exon 7 of the MTHFR gene. No unspecific alignments have been detected. No cross-reaction phenomena with genomic DNA have been reported.

## 2. Analytical sensitivity

A dilution assay was performed, using 1:4 serial dilutions of a wild type sample for A1298C polymorphism (no mutated alleles), an heterozygous sample (one mutated allele) and another homozygous one (two mutated allele), obtained by a conventional extraction system, at an initial concentration of 66.2, 114.7 and 55.1 ng/μL, respectively. The following analytical sensitivity results were obtained:

- Detection Limit of wild type and mutated allele = 6 ng/μl (\*)

(\*) Cp < 35

## 3. Diagnostic sensitivity and specificity

In a human genomic DNA study, 47 samples obtained from a laboratory were analysed. They were previously genotyped by another commercial kit.

All the tested samples were validated. The following results were obtained:

GENVINSET® MTHFR A1298C		A/A	A/C	C/C
Previous method	A/A	22	0	0
	A/C	0	17	0
	C/C	0	0	8

There is a 100% match in the results obtained with GENVINSET® MTHFR A1298C and the genotyping previously obtained with another commercial kit.



## Procedure limitations

- The method detects the A/C SNP at position 1298 of the exon 4 of the MTHFR gene.
- Mutations or polymorphisms at annealing primer/probe sites are possible and may result of the lack of allele definition. Other technologies could be necessary to resolve the typing.
- All the aforementioned conditions for the setup of the PCR should be carefully controlled. Any performances that do not meet such indications, can lead to poor results.
- All GENVINSET® components manipulation must be done according to general lab best practices and be adjusted to local regulations.
- The q-PCR thermal cycler must be calibrated according to the manufacturers' recommendations and should be used in accordingly to the manufacturer's instructions.
- Do not mix components from other kits or lot numbers.
- Do not use the kit after its expiration date.
- Do not use the kit if there are any suspicions of possible loss of reactivity, contamination, external box deterioration or any other incidence that might affect the kits performance.
- Data and result interpretation should be revised by qualified personnel.
- Discard expired reagents according to applicable regulations.

# Troubleshooting guide

## Problem

- Probable cause(s)
  - Suggested corrective measure(s)

## Reaction Blank (H<sub>2</sub>O) is positive

- **Primer Mix/Master Mix/Reaction Blank contamination**
  - Repeat the experiment with new Primer Mix/Master Mix/Reaction Blank aliquots
  - Handle the kit components always according to accepted lab practices in order to avoid contamination.
  - Verify manipulation and storage conditions
  - Discard contaminated reagents
- **Pre-PCR area is contaminated**
  - Confirm that all necessary precautions in the pre-PCR area have been followed
  - Check for possible contamination problems in other PCR techniques
  - Confirm suitability of the used reagents (1.5 ml tubes, pipette tips)
- **Pipetting error**
  - Check that the sample added corresponds to the one indicated on the worksheet

## Low or no signal in all samples. Control samples are OK

- **Bad quality of DNA**
  - Repeat the sample extraction verifying each step (Hemoglobin can interfere with the PCR)
- **Blood processed without previous frozen step**
  - Repeat extraction with a new blood aliquot previously frozen
- **Samples with very low DNA concentration**
  - Check DNA concentration
- **DNA samples with high concentration**
  - Perform the assay using diluted samples

**Fluorescence intensity too low**

- **Kit degradation**
  - Check that the storage of the kit is correct, reviewing both proper temperature conditions and light exposure (which should be avoided)
  - Avoid more than 3 freeze/thaw cycles of the Primer Mix vial
  - Aliquote the reagents if necessary
  - Repeat the test with new reagents

**Control C1 is positive**

- **Cross contamination**
  - Always handle the kit components following all necessary practices to avoid contamination
- **Pipetting error**
  - Check that the sample added corresponds to the one indicated on the worksheet

**Control C2 is negative**

- **Pipetting error**
  - Check that the sample added corresponds to the one indicated on the worksheet.

**Fluorescence intensity varies**

- **The dirtiness on the outside of the tube walls interferes with the signal**
  - Handle all consumables wearing gloves
- **Volume is not settled to the bottom of the well or there are bubbles**
  - Perform a brief centrifugation to ensure that the volume settles to the bottom of the well and to remove all bubbles
- **Pipetting error**
  - Check that the correct volume has been added

**There is no fluorescence signal**

- **Incorrect reading channels selected**
  - Set the appropriate reading channels
- **Pipetting error or reagent absence**
  - Control the pipetting and the reaction setup
  - Repeat the PCR
- **No reading channel was selected in the thermal cycler program.**
  - Check and modify the thermal cycler program

## References

1. Fodinger M, Horl WH, Sunder-Plassmann G. Molecular biology of 5,10-methylenetetrahydrofolate reductase. *J Nephrol.* 2000; 13:20–33.
2. Larsson SC, Giovannucci E, Wolk A. Folate intake, MTHFR polymorphisms, and risk of esophageal, gastric, and pancreatic cancer: a meta-analysis. *Gastroenterology.* 2006; 131:1271–1283.
3. Frosst P, Blom HJ, Milos R, Goyette P, Sheppard CA, Matthews RG, Boers GJH, den Heijer M, Kluijtmans LAJ, van den Heuvel LP, Rozen R. A candidate genetic risk factor for vascular disease: a common mutation in methylenetetrahydrofolate reductase. *Nat Genet.* 1995 May;10:111–113.
4. Friso S, Choi SW. Gene-nutrient interactions in one-carbon metabolism. *Curr Drug Metab.* 2005; 6:37–46.
5. Goldhaber S, Morrison RB. Pulmonary embolism and deep vein thrombosis. *Circulation.* 2002; 106: 1436–1438.
6. Goldhaber S, Fanikos J. Prevention of deep vein thrombosis and pulmonary embolism. *Circulation.* 2004; 110: e445–e447.
7. Goyette P, Sumner JS, Milos R, Duncan AMV, Rosenblatt DS, Matthews RG, Rozen R. Human methylenetetrahydrofolate reductase: isolation of cDNA, mapping and mutation identification. *Nat Genet.* 1994 Jun;7:195–200.
8. Leclerc D, Rozen R. Molecular genetics of MTHFR: polymorphisms are not all benign. *Med Sci (Paris).* 2007 Mar; 23(3):297–302.
9. Bailey LB, Gregory JF Jr. Polymorphisms of methylenetetrahydrofolate reductase and other enzymes: metabolic significance, risks and impact on folate requirement. *J Nutr.* 1999; 129:919–922.
10. Jacques PF, Bostom AG, Williams RR, Ellison RC, Eckfeldt JH, Rosenberg IH, Selhub J, Rozen R. Relation between folate status, a common mutation in methylenetetrahydrofolate reductase, and plasma homocysteine concentrations. *Circulation.* 1996; 93:7–9.
11. Ashfield-Watt PA, Pullin CH, Whiting JM, Clark ZE, Moat SJ, Newcombe RG, Burr ML, Lewis MJ, Powers HJ, McDowell IF. Methylenetetrahydrofolate reductase 677 C>T genotype modulates homocysteine responses to a folate-rich diet or low dose folic acid supplement: a randomized controlled trial. *Am J Clin Nutr.* 2002; 76:180–186.

12. Ulvik A, Ueland PM, Fredriksen A, Meyer K, Vollset SE, Hoff G, Schneede J. Functional inference of the methylenetetrahydrofolate reductase 677 C>T and 1298 A>C polymorphisms from a large-scale epidemiological study. *Hum Genet.* 2007; 121:57–64.

13. Botto LD, Yang Q. 5,10-Methylenetetrahydrofolate Reductase Gene Variants and Congenital Anomalies: A HuGE Review. *Am J Epidemiol.* 2000; 151(9):862–877.

14. Rengifo Ramos L, Gaviria Arias D. Genotypic analysis of c677t and a1298c polymorphisms in the methylene tetrahydrofolate reductase gene and a66g polymorphism in the methionine sintase reductase gene in down syndrome BAG. *J Basic Appl Genet* 2014;25:21-30.

15. González-Mercado M, Rivas F, Salazar-Páramo M, et al. Homocisteína, polimorfismos MTHFR C677T, A1298C y variables clínico-bioquímicas en población mexicana. *Acta Bioquímica Clínica Latinoamericana* 2014;48:23-31

16. Martínez-Echevarría MT, Casanueva-Calero K, Hernández-Acea B, Torres-Yibar W, García-Menéndez G. Polimorfismo MTHFR A1298C en pacientes cubanos con trombofilia. *Rev Hematol Mex.* 2019;20(1):5-10.

## Notice to purchaser

- This product has been developed for in vitro diagnostics purposes.
- BLACKHILLS DIAGNOSTIC RESOURCES, S.L.U. products should not be re-sold, modified for reselling or be used to manufacture other commercial products without written consent of BLACKHILLS DIAGNOSTIC RESOURCES, S.L.U.
- All information contained in this document can suffer modifications without prior notice. BLACKHILLS DIAGNOSTIC RESOURCES, S.L.U. does not assume any responsibility for possible errors in the document. This document is considered complete and accurate at the time of its publication. In any case will BLACKHILLS DIAGNOSTIC RESOURCES, S.L.U. be responsible for accidental, special, multiple or derived damages from the use of this document.
- The purchase of this product grants rights to the purchaser under certain Roche patents, only used to provide in vitro diagnostic services. It does not grant any generic patent or any other patents aimed at any other usage apart from the one specified.
- FAM™ and HEX™ are trademarks of Life Technologies Corporation.
- FAM™, HEX™ may be covered by one or more patents owned by Applied Biosystems, LLC. The purchase price of this product includes limited, nontransferable rights.
- TaqMan® is a registered trademark of Roche Molecular Systems, Inc.
- GENVINSET® is a trade mark of BLACKHILLS DIAGNOSTIC RESOURCES, S.L.U.

## Explanation of symbols used on the labels



For in vitro diagnostic use



This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical device



Catalogue number



Lot number



Expiration date



Contents sufficient for <n> tests



Manufactured by



Store at



Keep away from sunlight



Positive control