TSH

ELISA

Only for in-vitro diagnostic use

Instructions for use / Gebrauchsanweisung / Instrucciones de uso

English	2
Deutsch	7
Español	12
Abbreviations / Abkürzungen / Abreviaciones	18
Packaging materials / Verpackungsmaterialien / Materiales de embalaje	18
Symbols Key / Symbolschlüssel / Símbolos	19
Summary of Test Procedure / Kurzanleitung Testdurchführung / Resumen de la técnica	20

REF

TSH1030BA (96 Determinations)

ENGLISH

1. INTENDED USE

The TSH ELISA is intended for the quantitative determination of TSH (thyroid-stimulating hormone) in human serum or plasma (citrate, heparin).

2. PRINCIPLE OF THE ASSAY

The quantitative immunoenzymatic determination of TSH is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

The assay system utilizes two antibodies specific for different antigenic determinants of the TSH beta subunit. The first antibody is immobilized on the surface of the Microtiterplate. The second antibody is conjugated to horseradish peroxidase (HRP). The test sample is allowed to react simultaneously with the two antibodies, resulting in TSH molecules being sandwiched between the solid phase and the enzyme-linked antibodies. Subsequently, the wells are washed to remove all unbound material. The immune complexes are visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is proportional to the amount of TSH in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

3. MATERIALS

3.1. Reagents supplied

- **Microtiterplate:** 12 break apart 8-well snap-off strips coated with anti-TSH antibodies; in resealable aluminium foil.
- **SOLN STOP:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **WASH BUF 20x:** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap; 0.2% (w/v) 5-Bromo-5-nitro-1,3-dioxane.
- **Conjugate:** 1 bottle containing 15 mL of peroxidase labelled antibodies to TSH in phosphate buffer (10 mM); coloured red; ready to use; black cap.
- **SUB TMB:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Control:** 1 bottle containing 1 mL of a lot-specific control; coloured yellow; ready to use; white cap; ≤ 0.02% (v/v) MIT.
- **Standards:** 7 vials, each containing 1 mL standard; coloured yellow; ready to use; yellow cap; ≤ 0.02% (v/v) MIT.

Standard A:	0.0	µIU/mL
Standard B:	0.2	µIU/mL
Standard C:	0.5	µIU/mL
Standard D:	2.5	µIU/mL
Standard E:	5.0	µIU/mL
Standard F:	10.0	µIU/mL
Standard G:	20.0	µIU/mL

The standards are calibrated in accordance with the "WHO International Standard; Thyroid Stimulating, Human for Immunoassay; NIBSC 81/565".

For hazard and precautionary statements see 11.1.

3.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

3.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplate
- Pipettes to deliver volumes of 50 and 100 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

4. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

5. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

5.1. Microtiterplate

The break-apart snap-off strips are coated with anti-TSH antibodies. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

5.2. WASH BUF 20x

Dilute WASH BUF 20x 1 + 19; e.g. 10 mL WASH BUF 20x + 190 mL distilled water. The diluted buffer (WASH BUF 1x) is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g., in a water bath. Mix well before dilution.

5.3. **SUB TMB**

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. **SUB TMB** should be colourless or could have a slight blue tinge. If **SUB TMB** turns into blue, it may have become contaminated and should be thrown away.

6. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise, they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

7. ASSAY PROCEDURE

Please read the instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the volume of (**WASH BUF 1x**) from 300 µL up to 350 µL to avoid washing effects. Pay attention to chapter 11. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 50 µL standards/controls and samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Dispense 100 µL conjugate into all wells except for the blank well A1.
3. Cover wells with the foil supplied in the kit.
4. **Incubate for 60 min ± 5 min at 37 ± 1 °C.** Do not expose to direct sunlight.
5. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well five times with 300 µL of **WASH BUF 1x**. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
6. Dispense 100 µL **SUB TMB** into all wells except for the Substrate Blank well A1.
7. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
8. Dispense 100 µL **SOLN STOP** into all wells in the same order and at the same rate as for the **SUB TMB**, thereby a colour change from blue to yellow occurs.
9. Measure the absorbance at 450/620 nm within 30 min after addition of the **SOLN STOP**.

7.1. Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

8. RESULTS

8.1. Run Validation Criteria

In order for an assay run to be considered valid, these instructions for use have to be strictly followed and the following criteria must be met:

- **Substrate blank:** Absorbance value < **0.100**
- **Standard A:** Absorbance value < **0.100**
- **Standard B:** Absorbance value > **Standard A**
- **Standard C:** Absorbance value > **Standard B**
- **Standard D:** Absorbance value > **Standard C**
- **Standard E:** Absorbance value > **Standard D**
- **Standard F:** Absorbance value > **Standard E**
- **Standard G:** Absorbance value > **1.000**
- **Control:** Result in µIU/mL within the range indicated on the label.


Standard A < Standard B < Standard C < Standard D < Standard E < Standard F < Standard G

If these criteria are not met, the test is not valid and must be repeated.

8.2. Calculation of Results

In order to obtain **quantitative results in $\mu\text{IU/mL}$** plot the (mean) absorbance values of the Standards A – G (y-axis) on graph paper in a system of coordinates against their corresponding concentrations (x-axis) and draw a standard curve. Read results from this standard curve employing the (mean) absorbance values of each patient sample and control.

8.3. Typical Standard Curve

8.4. Interpretation of Results

Normal value ranges for this ELISA should be established by each laboratory based on its own patient populations in the geographical areas serviced.

The following values should be considered as a guideline:

normal: 0.3 – 4.5 $\mu\text{IU/mL}$

Test results should be verified concerning the entire clinical status of the patient. Also, every decision for therapy should be taken individually.

Diagnosis of a disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.

9. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

9.1. Precision

Intraassay	n	Mean (E)	CV (%)
#1	24	0.064	10.5
#2	24	0.519	6.7
#3	24	0.207	5.4
Interassay	n	Mean ($\mu\text{IU/mL}$)	CV (%)
#1	12	5.426	5.2
#2	12	2.156	10.2
#3	12	2.448	4.0
#4	12	1.402	4.0

9.2. Analytical Sensitivity

The analytical sensitivity (according to CLSI EP17-A) is defined as the apparent concentration of the analyte that can be distinguished from the zero calibrator. It is 0.07 $\mu\text{IU/mL}$.

9.3. Specificity

The following hormones were tested for cross-reactivity:

Hormone tested	Concentration [mIU/mL]	Produced intensity equivalent to TSH [$\mu\text{IU/mL}$]
LH (WHO 2 nd IS 80/552)	70	-
	175	-
	350	< 0.1
hCG (WHO 5 th IS 07/364)	1,000	-
	5,000	-
	25,000	< 0.1
FSH (WHO IS 83/575)	10	-
	20	< 0.5
	50	< 0.5

9.4. Correlation

The TSH ELISA was compared with a reference electrochemiluminescence immunoassay. 47 biological samples were used (the values ranged from approx. 0 - approx. 11 $\mu\text{IU/mL}$). The linear regression curve was calculated.

$$y = 1.06x - 0.11$$

$$r^2 = 0.97$$

9.5. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

9.6. Hook Effect

No hook effect was observed applying up to 5,000 $\mu\text{IU/mL}$ TSH.

9.7. Measurement range

The measurement range is 0.07 $\mu\text{IU/mL}$ – 20 $\mu\text{IU/mL}$.

10. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

11. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

11.1. Safety note for reagents containing hazardous substances

Reagents may contain CMIT/MIT (3:1) or MIT (refer to 3.1).

Therefore, the following hazard and precautionary statements apply.

Warning	H317	May cause an allergic skin reaction.
	P261	Avoid breathing spray.
	P280	Wear protective gloves/ protective clothing.
	P302+P352	IF ON SKIN: Wash with plenty of soap and water.
	P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
	P362+P364	Take off contaminated clothing and wash it before reuse.

Reagents may contain 5-Bromo-5-nitro-1,3-dioxane (refer to 3.1).

Therefore, the following hazard and precautionary statements apply.

Warning	H315	Causes skin irritation.
	H319	Causes serious eye irritation
	P280	Wear protective gloves/ protective clothing.
	P302+P352	IF ON SKIN: Wash with plenty of soap and water.
	P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
	P337+P313	If eye irritation persists: Get medical advice/attention.

Further information can be found in the safety data sheet.

11.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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
SYMBOLS KEY / SYMBOLSCHLÜSSEL / SIMBOLOS

	Manufactured by / Hergestellt von / Fabricado por
IVD	In Vitro Diagnostic Medical Device / In Vitro Diagnosticum / Producto para diagnóstico In vitro
LOT	Lot Number / Chargenbezeichnung / Número de lote
	Expiration Date / Verfallsdatum / Fecha de caducidad
	Storage Temperature / Lagertemperatur / Temperatura de almacenamiento
CE	CE marking / CE-Kennzeichnung / Marca CE
UDI	Unique Device Identifier / Eindeutige Produktidentifizierung / identificación única del producto
REF	Catalogue Number / Katalog Nummer / Número de Catálogo
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulte las Instrucciones de Uso
MTP	Microtiterplate / Mikrotiterplatte / Placa de Microtitulaciones
CONJ	Conjugate / Konjugat / Conjugado
CONTROL	Control / Kontrolle / Control
CAL	Standard or Calibrator A-G / Standard oder Kalibrator A-G / Estándar o Calibrador A-G
SOLN STOP	Stop Solution / Stopplösung / Solución de Parada
SUB TMB	TMB Substrate Solution / TMB-Substratlösung / Solución Substrato TMB
WASH BUF 20x	“Washing Buffer (20x concentrated)”; REF W0000 Washing Buffer 20x concentrated / Waschpuffer 20x konzentriert / Tampón de lavado concentrado x20
WASH BUF 1x	20-fold dilution of WASH BUF 20x / 20-fach Verdünnung von WASH BUF 20x / Dilución de 20 veces del WASH BUF 20x
	Contains sufficient for “n” tests / Ausreichend für “n” Tests / Contenido suficiente para “n” tests

SCHEME OF THE ASSAY

TSH

Test Preparation

Prepare reagents and samples as described.
 Establish the distribution and identification plan for all samples and standards/controls.
 Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate Blank (A1)	Standard A - G	Control	Sample
Standard A-G	-	50 µL	-	-
Control	-	-	50 µL	-
Sample	-	-	-	50 µL
Conjugate	-	100 µL	100 µL	100 µL
Cover wells with foil supplied in the kit Incubate for 60 min at 37 ± 1 °C Do not expose to direct sunlight Wash each well five times with 300 µL of WASH BUF 1x				
SUB TMB	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature (20...25°C) in the dark				
SOLN STOP	100 µL	100 µL	100 µL	100 µL
Photometric measurement at 450 nm (reference wavelength: 620 nm)				

Bioactiva Diagnostica GmbH

Louisenstrasse 137
 61348 Bad Homburg, Germany

Tel: +49 (0) 6172 17102-0
 Fax: +49 (0) 6172 17102-29
 Email: bioactiva@bioactiva.de
 Internet: www.bioactiva.de

Date of issue: 2024-04-30
 TSH1030BA_IFU_rev01_fromLot_047N