# **BORRELIA IgM ELISA RECOMBINANT ANTIGEN**

ENZYME IMMUNOASSAY FOR THE QUALITATIVE OR QUANTITATIVE DETERMINATION OF  
IgM-ANTIBODIES TO BORRELIA IN PLASMA, SERUM OR CEREBROSPINAL FLUID  
CAT. NO. BI-21042 12 X 8 TESTS

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*Additional Information about our products is available on our website.*

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**[www.bmgrp.com](http://www.bmgrp.com)**

## 1) INTRODUCTION

Lyme borreliosis is a bacterial infection caused by the spirochaete *Borrelia burgdorferi sensu stricto*, *Azelloi* or *Garinii*, and is characterised by a variety of clinical symptoms. It can be divided into 3 stages:

Stage 1, early dermatitis, Clinical: erythema migrans.

Stage 2, early disseminated infection, Clinical: lymphocytic meningoradiculitis (Bannwarth's syndrome), neuroborreliosis.

Stage 3, late disseminated infection, Clinical: chronic progressive encephalomyelitis, acrodermatitis chronica atrophicans (ACA), chronic arthritis.

To improve the diagnostic specificity, the Biomedica microwell strips are coated with the following recombinant antigens:

- p21 = OspC (outer surface protein C): *B. afzelii* (pKo)
- p21 = OspC (outer surface protein C): *B. garinii* (20047)
- p41/I = (inner part of flagellin): *B. bavariensis* (pBi)
- VlsE: fusion proteins of different *Borrelia* genospecies

## 2) CONTENTS OF THE KIT

All reagents except the wash buffer are ready for use. Store all reagents at 2-8°C.

CONT	KIT COMPONENTS	QUANTITY
PLATE	Microwell strips, breakable, colour coded, coated with recombinant antigen	12 x 8 strips
STD	Standards: positive (pos), cut-off (c.o.), negative (neg), ready to use	Each 1 x 1.5 ml
CONJ	Conjugate: rabbit anti human IgM HRPO, ready to use	1 x 14 ml
DIL	Sample diluent with green dye, ready to use	1 x 105 ml
SUB	Substrate: TMB substrate solution, ready to use	1 x 13 ml
WASHBUF	Wash buffer: 20x concentrated	1 x 50 ml
STOP	Stop solution: sulphuric acid, ready to use	1 x 10 ml

## 3) ADDITIONAL MATERIAL ADDED TO THE KIT

- 2 self-adhesive plastic films
- Quality control sheet
- Protocol sheet
- Instruction manual for use

## 4) EQUIPMENT REQUIRED BUT NOT SUPPLIED

- Precision pipettes calibrated to deliver 10-1000 µl and disposable tips
- Multichannel pipette
- ELISA reader for absorbance at 450 nm (or from 450 nm to 630 nm)
- Graph paper or software for calculation of results
- Plate washer is recommended for washing

## 5) REAGENTS AND SAMPLE PREPARATION

### Sample preparation:

The sample can be serum or plasma (EDTA, citrate, heparin, CPD), which must be separated as soon as possible from the blood clot after blood sampling to avoid a haemolysis. CSF can also be used. Avoid microbial contamination of the samples. Insoluble substances must be removed from the sample before incubation. The use of heat-inactivated, icteric, haemolytic, lipemic or turbid samples is not recommended.

If the tests are not carried out immediately, the samples can be stored for up to 2 weeks at 2-8°C. Prolonged storage of the samples is possible at -25°C or below. Repeated freezing and thawing of samples is not recommended due to the risk of inaccurate results. Avoid more than 3 freeze-thaw cycles.

### Sample dilution:

Dilute human serum or plasma 1:101 (1 + 100) with sample diluent, e.g. 10 µl of sample + 1 ml of sample diluent, mix well.

Dilute cerebrospinal fluid 1:2 (1 + 1) with sample diluent, e.g. 100 µl of sample + 100 µl of sample diluent, mix well.


For samples with results >30 BBU/ml the samples can be diluted with sample dilution buffer and retested. Obtained results need to be multiplied with respective dilution factor for final concentration.

Coniugate: The conjugate is ready for use and should be used within 2 weeks after opening. For long term storage please aliquot and freeze at -25°C unused conjugate. Frozen aliquots are stable until expiry date stated on the original vial.

Preparation WASHBUF (Wash buffer): Dilute 20x conc. wash buffer with distilled water to 1000 ml (50 ml conc. + 950 ml dist. water). Mix well. Crystals that may appear will dissolve at room temperature. Undiluted WASHBUF is stable 2-8°C until expiry date stated on the label. The diluted wash buffer is stable at 2-8°C for one month.

**STD: standards (pos, cut-off, neg) are ready to use and MUST NOT be diluted.**

## 6) PRINCIPLE OF THE ASSAY

## 7) ASSAY PROTOCOL

All reagents and samples must be at room temperature (18-26°C) before use in the assay.
Mark position for BLANK/STD/SAMPLE/ (Blank/Standard/Sample) on the protocol sheet.
Take microtiter strips out of the alu bag, take a minimum of one well as Blank. Store unused strips with desiccant at 2-8°C in the alu bag. Strips are stable until expiry date stated on the label.

- |  |
|--|
| 1. Add 100 µl STD/SAMPLE (Standards/diluted samples) into the respective wells, except blank.  |
| 2. <b>Cover strips with plastic film and incubate for 1 hour at 37°C in an incubator.</b>  |
| 3. Aspirate and wash wells 4x with at least 250 µl diluted WASHBUF (Wash buffer). Remove remaining WASHBUF by hitting plate against paper towel after the latest wash. |
| 4. Add 100 µl CONJ (anti human IgM-HRPO) into each well.   |
| 5. <b>Cover strips with plastic film and incubate for 30 min at room temperature (18-26°C).</b>  |
| 6. Aspirate and wash wells 4x with at least 250 µl diluted WASHBUF (Wash buffer). Remove remaining WASHBUF by hitting plate against paper towel after the latest wash. |
| 7. Add 100 µl SUB (TMB substrate) into each well.  |
| 8. <b>Incubate 15 min at room temperature (18-26°C), in the dark.</b>  |
| 9. Add 50 µl STOP (Stop solution) into each well. Swirl gently.  |
| 10. Measure absorbance immediately at 450 nm with reference 630 nm, if available.  |

The use of an automated system is possible, but sometimes needs modifications. If you need support, please contact your local representative.

## 8) CALCULATION OF RESULTS

Subtract the OD of the blank from the OD's of the samples and standards.

Calculate the mean OD's of all samples and standards.

The OD of the cutoff standard represents 10 BBU/ml (Biomedica Borrelia Units). This is the basis for calculation of the results.

$$\frac{\text{Ø OD of the sample}}{\text{Ø OD of the cut off control}} \times 10 = \text{BBU / ml}$$

The quality control protocol supplied with the kit shows the results of the final release QC for each kit. Data for optical density obtained by customers may differ due to various influences and/or due to the normal decrease of signal intensity during shelf life. However, this does not affect validity of results as long as data mentioned in chapter 9) ASSAY CHARACTERISTICS *part valid results* are achieved.

### 8.1. Qualitative calculation for serum/plasma

Positive result	samples with an OD of more than cut-off +10% ( $\geq 11$ BBU/ml).
Borderline result	samples between 9 - 11 BBU/ml. These patients should be retested with a new sample.
Negative result	samples with an OD of less than cut-off -10% ( $\leq 9$ BBU/ml).

### 8.2. Quantitative calculation for serum/plasma

Calculation of the BBU/ml is identical to the qualitative calculation, section 8.1.

The test is linear between 5 - 30 BBU/ml. Values outside this range are semi-quantitative.

Alternatively, samples with results >30 BBU/ml can be diluted with sample dilution buffer and retested. Obtained results need to be multiplied with respective dilution factor for final concentration.

Weak positive	Samples between 11 - 20 BBU/ml
High positive	Samples between 21 - 30 BBU/ml
Very high positive	Samples above 30 BBU/ml

### 8.3. Qualitative calculation for cerebrospinal fluid

The calculation of the BBU/ml is identical to the qualitative calculation, section 8.1.

Normal values CSF: The mean of 30 CSF samples showed 1.3 BBU/ml.

Neuroborreliosis negative	Samples with values $\leq 5$ BBU/ml
Suspect for Neuroborreliosis	Samples with values $> 5$ BBU/ml

## 9) ASSAY CHARACTERISTICS

Valid results:

Blank	The OD must be $\leq 0.250$
Negative standard	The OD must be $\leq 0.200$ (blank subtracted)
Cut-off standard	The OD must be at least twice the negative standard (blank subtracted).
Positive standard	The OD must be at least twice the OD of the cut-off standard (blank subtracted).

Crossreactivities:

	n	crossreactive
rheumafactor-positive	15	1
syphilis-positive	7	0

This ELISA uses recombinant antigens. Nevertheless, a crossreaction to autoimmune-positive samples cannot be completely excluded.

## 10) PRECISION

Intra-assay	pos (n=12)	high pos (n=12)	Inter-assay	pos (n=17)	high pos (n=20)
Mean (BBU/ml)	24	38	Mean (BBU/ml)	22	36
SD (BBU/ml)	2	3	SD (BBU/ml)	2	3
CV (%)	6.7	6.8	CV (%)	8.3	8.1

## 11) TECHNICAL HINTS

- Do not mix or substitute reagents with those from other lots or sources.
- Avoid foaming when mixing reagents.
- Do not mix stoppers and caps from different reagents or use reagents between lots.
- Do not use reagents beyond expiration date. Protect reagents from direct sunlight.
- Substrate solution should remain colourless until added to the plate.

## 12) PRECAUTIONS

All test components of human source were tested against HIV-Ab, HCV-Ab and HBsAg and were found negative. Nevertheless, they should be handled and disposed as if they were infectious. Liquid reagents contain  $\leq 0.1\%$  Proclin 950 as preservative. Avoid contact with skin and mucous membrane. Proclin 950 is not toxic in concentrations used in this kit. It may cause allergic skin reactions, avoid contact with skin or eyes.

- Do not pipette by mouth.
- Do not eat, drink, smoke or apply cosmetics where reagents are used.
- Avoid all contact with the reagents by using gloves.
- Sulfuric acid is irritating to eyes and skin. Flush with water if contact occurs. Avoid contact with skin and mucous. Irritations are possible. - Flush with water after contact!

## 13) LITERATURE

- "Use of Recombinant Antigens of *Borrelia burgdorferi* in Serologic Tests for Diagnosis of Lyme Borreliosis", MAGNAREELI L.A. et al. (1996), *Journal of Clin. Microbiol.* 34: 2 pp 237-240
- "Immunological and Molecular Heterogeneity of the Outer Surface Protein C (OspC) of *Borrelia burgdorferi* Sensu Lato", HEIPKE S.J et al., pp 11-14
- "Identification of a Protein in Several *Borrelia* Species Which Is Related to OspC of the Lyme Disease Spirochetes, MARKONI R.T. et al. (1993), *Journal of Clin Microbiol.*, pp 2577-2583
- "Antigenic Variation and Variation in Expression of Immunodominant B. *Burgdorferi* Antigens: Implications for using Recombinant Antigens for Serodiagnosis", WILSKE B. et al., pp 179-182
- Review. "Diagnosis of Lyme Borreliosis in Europe Vector-borne and Zoonotic Disease.3", No. 4, 215 -227 (2003), Wilske B.
- "Characterisation of the vls antigenic variation loci of the Lyme disease spirochaetes *Borrelia garinii* Ip90 and *Borrelia afzelii* ACA1" (2003). Eicken C. et al. *Mol. Microbiol.* 47 (5), 1407-1417

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# **BI-21042 ANTI BORRELIA IgM**

## **ASSAY PROTOCOL AND CHECKLIST**

- Bring all reagents to room temperature (18-26°C).
- Prepare reagents and samples as instructed.
- Bring unused and prepared components to the storage temperature mentioned in the package insert.
- Take microtiter strips out of the alu bag and mark positions on the protocol sheet.
  
- 1. Add 100 µl STD/SAMPLE (standard/ diluted sample) into each well, except blank.
- 2. Cover tightly and incubate for 1 hour at 37°C.**
- 3. Aspirate and wash wells with at least 250 µl WASHBUF (Wash buffer) four times. Remove remaining buffer by hitting plate against paper towel.
- 4. Add 100 µl CONJ (Conjugate) into each well.
- 5. Cover tightly and incubate for 30 minutes at room temperature (18-26°C).**
- 6. Aspirate and wash wells with at least 250 µl WASHBUF (Wash buffer) four times. Remove remaining buffer by hitting plate against paper towel.
- 7. Add 100 µl SUB (Substrate) into each well.
- 8. Incubate for 15 minutes at room temperature (18-26°C), in the dark.**
- 9. Add 50 µl STOP (Stop solution) into each well. Swirl gently.
- 10. Read Optical Density at 450 nm with reference 630 nm, if available.