- de** Gebrauchsanweisung
- en** User's Manual
- fr** Mode d'emploi
- it** Istruzioni per l'uso
- es** Instrucciones de empleo
- pt** Manual de operação
- ru** Руководство по применению

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
## Dermatoskope

QM-1-101F

Kontaktscheibe mit Skalierung für Dermatoskope Ø 25 mm  
Contact disk with scale setting for dermatoscopes Ø 25 mm  
Rondelle de contact avec cadran pour dermatoscopes Ø 25 mm  
Vetrino di contatto con scala graduata per dermatoscopi Ø 25 mm  
Arandela de contacto con escala para dermatoscopios Ø 25 mm  
Disco de contacto com escala para dermatoscópios Ø 25 mm  
Контактное стекло со шкалой Ø 25 мм

Fokuseinstellung  
Focus adjustment  
Plage de focalisation  
Messa a fuoco  
Margen de la focalización  
Ajuste do foco  
Установка фокуса

Ladebatterie (Akku)  
Rechargeable battery  
Pile rechargeable (accu)  
Batteria ricaricabile  
Pila recargable (acumulador)  
Pilha recarregável (acumulador)  
Аккумулятор

Batterien  
Batteries  
Piles  
Batterie  
Pilas  
Pilhas  
Батарейки

<b>de</b>	Dermatoskop .....	4
<b>en</b>	Dermatoscope .....	10
<b>fr</b>	Dermatoscope .....	16
<b>it</b>	Dermatoscopio .....	22
<b>es</b>	Dermatoscopio .....	28
<b>pt</b>	Dermatoscópico .....	34
<b>ru</b>	Дерматоскопы .....	40

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Dear Customer, we thank you for choosing a KaWe product. Our products are known for their high quality and longevity. This KaWe product meets EC standards 93/42/EEC (standards for medical products).

**Please read these instructions thoroughly and carefully before attempting to use this product and heed the given instructions.**

Familiarise yourself fully with this product before using it.

**Use:** The dermatoscopes are only to be used by authorised professional personnel.

**Purpose:** The KaWe dermatoscope is used for the optical examination of the skin. This tool magnifies the superficial skin structure approx. 10-times and illuminates the examined dermatomic area with a halogen bulb that is comparable to incident light microscopy (epiluminescence microscopy). This dermatoscope is used as a tool for the examination of pigmentary skin modifications, especially for the early diagnosis of malignant tumours. Moistening the skin with the special dermatoscope gel (KaWe Derma-Gel) prevents cutaneous reflexes and makes skin discolouration more recognizable and more easily seen.

**The unit has been exclusively designed to be used for the examination of the exterior surface of the skin.**

### Dermatoscopes:

- **KaWe EUROLIGHT®** D30 Dermatoscope 2.5 V
- **KaWe PICCOLIGHT®** D Dermatoscope 2.5 V
- **KaWe EUROLIGHT®** D30 Dermatoscope 3.5 V / 240 V (EU)

### Operation of of KaWe EUROLIGHT® D30 2.5 V and 3.5 V dermatoscopes:

Satisfactory operation of these instruments is guaranteed with the use of standard 2.5 V KaWe handles with KaWe-brand 2.5 V or 3.5 V rechargeable batteries and/or with 2, 1.5 V "C" batteries. Slip the instrument head onto the handle and lock it into place. Turn on the sliding switch by pushing and moving the white plastic button to the left. The light intensity (0=off 1=maximum brightness) can be easily regulated with the rheostat. Turn off the instrument by sliding the switch back to the right to its starting position. To change the batteries, unscrew the bottom locking cap. Rechargeable handles can be charged with the KaWe MedCharge® 4000 charging station by removing the plastic cap at the handle's bottom and then inserting the handle into the charging slot of the KaWe MedCharge® 4000. The KaWe EUROLIGHT® D30 3.5 V has two handle parts that are screwed together. These two parts can be separated from one another and the bottom part of the handle (the part with the plug) can be plugged into an outlet to be charged. (CAUTION: EU model has 240 V/3.5V!). Please only use 3.5 V KaWe-brand rechargeable batteries. The instrument head can be loosened and removed by pushing downward on the sleeve.

**!** **CAUTION:** Heed the additional care instructions that will prolong the life of the rechargeable batteries.

**Operation KaWe PICCOLIGHT® D dermatoscope:** A satisfactory operation is guaranteed with the use of KaWe PICCOLIGHT® handles with 2 batteries ("AA"-Mignon 1.5V) To turn on the instrument, use the sliding on/off switch. To turn the instrument on, slide the switch toward the dermatoscope head. Turn off the unit by sliding the switch in the opposite direction and back to its starting position. Remove and replace the battery and the lamp by removing the KaWe PICCOLIGHT®-head.

**Replacing the lamp in the KaWe PICCOLIGHT® Dermatoscope:** Satisfactory operation of the unit is only guaranteed if original KaWe lamps are used as only these ensure good light quality. Unscrew the dermatoscope head from the handle and extract the lamp by pulling it out gently. If necessary, clean the glass bulb of the new lamp with alcohol. The glass bulb must be properly clean and free of finger prints (free of grease). The lamp must be pushed carefully into the provided slot until it touches the lateral pin. Next, screw the instrument head onto the handle.

**Replacing the lamp in the KaWe EUROLIGHT® Dermatoscope:** Satisfactory operation of the unit is only guaranteed if original KaWe lamps are used as only these ensure good light quality. Unlick the dermatoscope head from the handle and extract the lamp by pulling it out gently. Then, proceed in the exact same manner as described above for bulb replacement of the KaWe PICCOLIGHT®. The glass bulb must be properly clean and free of finger prints (free

of grease). The lamp must be pushed carefully into the provided slot until it touches the lateral pin.

**Use:** Moisten the area of the body that is to be examined with the special dermatoscope gel KaWe Derma-Gel. Switch on the light and place the dermatoscope so that the contact disk rests lightly on the skin. This will give you an evenly intensive illumination of the skin area with an enlarged view ( $\emptyset$  25 mm). The area that is to be examined should be placed approx. in the middle of the contact disk. The eye should be as close to the ocular as possible. Turn the focussing ring with your free hand until a perfectly focused view is obtained. Try to look through the ocular with a relaxed eye. The contact disk with the scale setting measures pigmented skin modifications. In the case of elevated skin or skin that is sensitive to the touch (injured) the dermatoscope can be used without the contact disk. In order to ensure a precise view, the dermatoscope must be placed within a maximum distance of 5 mm from the skin.

#### Technical data:

- Magnifying glass with approx. 10-fold magnification
- Contact disk diameter: 25 mm
- Focussing area/focus adjustment: -6 to +3.5 dpt
- Light intensive halogen bulb 2.5 V and/or 3. V (KaWe EUROLIGHT®)
- Light intensive vacuum bulb 2.5 V (KaWe PICCOLIGHT® D)

#### Further information, maintenance, storage:

When used and stored properly, these otoscopes will serve you many years.

**Guarantee:** When used under normal circumstances and with attention to these instructions we guarantee this product for two years after the date of purchase (with the exception of lamps and rechargeable batteries). Should you need further information or should your instrument require repair, please contact your dealer.

 **Warning instructions:** Never bring the dermatoscope in contact with liquid and make sure that no liquid penetrates the housing! Since the dermatoscope is not designed for operations, simple cleaning and/or disinfection with an alcohol-based surface disinfectant is sufficient. Do not use any abrasive cleansing agents! The dermatoscope cannot be sterilized.

**Limitation of reprocessing:** Frequent reprocessing influences the instruments only minimally. The working life of the instruments usually depends on the wear and tear and damage by use.

### Instructions

**Workstation:** Remove dirt accumulation on surfaces with a disposable/paper cloth.

**Storage and transport:** No special requirements

**Cleaning preparation:** No special requirements

**Manual cleaning:** The exterior surface of the instrument can be cleaned with a damp, soft and lint-free cloth. Never allow the dermatoscope to come into contact with fluids and ensure that liquids do not penetrate the housing. As the dermatoscope is not intended for use during operational procedures, simple cleaning and or disinfection with an alcohol-based

surface disinfecting agent is sufficient. Do not use abrasive cleaning agents! The unit cannot be sterilized. The magnifying glass is to only be cleaned with glass cleaner.

**Automatic cleaning:** No automatic cleaning possible/required.

**Disinfection:** No special requirements. The disinfectant can be used in accordance with the instructions given on the label.

**Sterilisation:** The measuring scale disc must be sterilised in the autoclave at a temperature of 134°C for a minimum of 3 minutes.

**Additional information:** If many instruments are being sterilised in the autoclave at the same time, the maximum load given by the manufacturer must not be exceeded.

**Maintenance:** Defective parts must be removed and/or replaced. If the instrument does not function properly, follow these steps: **1.** Replace batteries and ensure that they are charged. **2.** Replace the lamp and ensure that it works properly. **3.** Send the unit in for repair.

**Check and performance test:** Check all movable parts such as the frames of the magnifying glasses and plugin connections to ensure that they all move smoothly. Turn on the unit. For all units: Visually check all parts for damage and wear. Keep all electrical contact points clean.

**Packaging:** Individually packed in zippered or cloth bags.

**Storage:** No special requirements.

**Disposal:** The product is to be disposed of at a separate waste disposal location in accordance with local ordinances.

### **User's Manual/Maintenance Instructions for Rechargeable NiMH/Li-Ion Batteries**

**Purpose:** The charging batteries are to be used as an energy source for the relevant KaWe products.

**Unsuited use/contraindication of the products:** Any use other than that described here is not in accordance with the intended use of the unit. The manufacturer is not liable for any resulting damages. There are no known risks. The user alone bears the risk. Use only the enclosed dermatoscope handle to supply power to the instrument.

**Further information, maintenance, storage:** The rechargeable batteries can be cleaned with a damp cloth and then be wiped off with a dry fabric. Please ensure that no short-circuit is produced.

The use of rechargeable batteries that are not explicitly approved for these units may cause explosions or leakages and the risk of fire, injuries and damage to surroundings. In the case that a rechargeable battery begins to leak and spilled substances come into contact with eyes, mouth, skin or clothing, immediately rinse the effected area with water and consult a doctor.

**Further maintenance instructions that may prolong the lifetime of the rechargeable battery:** Please heed the following instructions: The

non-observance of these instructions may result in the battery losing its full capacity within a short period of time. With NiMH technology, a lazy battery effect can occur that is analogous to the classic memory effect. It is triggered by the energy that has not been completely spent prior to recharging the battery. This effect can be nearly eliminated by completely discharging the battery from time to time and then by recharging it two or three times in quick succession. Subsequently, the accumulator will regain its full capacity. If the accumulator no longer provides the amount of energy necessary for powering the unit, immediately turn off the unit and charge the battery (in order to prevent it from losing its charge completely, which would result in the destruction of the battery).

Please keep in mind that the rechargeable batteries can spontaneously discharge. When the batteries are not used (charged) for a long period of time, a self-discharge may occur. In this case, the capacity will be lost and the batteries may become permanently defective. The manufacturer is not liable for any resulting damages. Therefore, when stored for more than one month, the batteries should be partially charged to approx. 50%.

Please always respect the given charging time in order to prevent overcharging and the resulting destruction of the battery (the KaWe charging station KaWe Med-Charge® 4000 has an integrated automatic shut-off). NiMH rechargeable batteries lose their capacity at low temperatures (as low as 0°C).

**Risks and notice of possible danger:** **1.** Avoid short-circuiting the rechargeable batteries. **2.** Never immerse the rechargeable batteries in any kind of

liquid. **3.** Keep the rechargeable batteries away from heating installations and open flame. **4.** Batteries are not to be immersed in water or salt water. **5.** Never try to disassemble or modify the rechargeable batteries or expose them to heat. **6.** Avoid dropping the batteries and protect them against being hit or crushed as this may damage the housing. **7.** Please only use the rechargeable batteries and accessory parts recommended by KaWe. **8.** When using a rechargeable battery for the first time, please make sure to charge it completely. **9.** When charging is completed and the charging unit and the compact adapter are no longer in use, please disconnect them from the unit and the socket in order to avoid the risk of fire and other dangers. **10.** Do not cover the charging unit with any kind of material (tablecloths, carpets, bedclothes etc.) while operating the instrument. If you use the adapter for a long period of time, it may become overheated, lose its original shape or catch fire. **11.** Please only charge the rechargeable batteries in the appropriate battery charging unit. The charging unit and the compact adapter have been designed especially for use with your unit. Please do not use them with other products or rechargeable batteries. All other products may cause the risk of excessive heat and distortion and thus the risk of fire and electric shock. **12.** Prior to disposing of the rechargeable battery, ensure that the contacts are masked with adhesive strips or other insulating material in order to prevent their direct contact with other materials. Contact with other metallic material in rubbish bins may cause a fire and injury.

### Charging time:

REF	Capacity [mAh]	Charging time [h]
12.80110.722	3000	ca. 14*
12.80120.742	700	ca. 12**
12.80220.722	2200	ca. 10*

\*Charging time after complete discharge of the rechargeable battery in the KaWe charging station KaWe MedCharge® 4000.

\*\*Charging time after complete discharge of the rechargeable batteries in a wall socket.

**Disposal:** Please dispose of your rechargeable batteries at a separate location for electric and electronic devices.

**Disposal:** Please dispose of expired NiMH rechargeable batteries and other batteries in accordance with legal regulations.

**Common accessories:** Further information about these products can be obtained in our Internet shop at: [www.kawemed.de](http://www.kawemed.de).

Manufacturer: **KaWe**

**Contacting the manufacturer:** Address or telephone number of the dealer or dial +49-7141-68188-0.

## Symbol key:

	Manufacturer
	Date of manufacture
	Heed the User's Manual
	Caution!
	Separate disposal of electric and electronic devices
	Dispose of expired batteries separately at an appropriate collection centre
	Complies with relevant EU guidelines
	Only use device in dry rooms
	GOST-R (Russia)
	Keep dry
	Temperature limit

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**KIRCHNER & WILHELM** GmbH + Co. KG