



MEDICAL PROFESSIONALS'
TRUSTED BRAND

UROLOGY

ABOUT

Blueneem Medical Devices Pvt. Ltd. is a leading Developer, Manufacturer and Marketer of medical devices in the fields of Urology, Gastroenterology and Interventional Radiology.

Founded in 2009 on core values and principles of Patient Wellbeing, Integrity & Fairness, Execution Excellence, Innovation and Continuous Improvement, we forge ahead to make a difference to all our stakeholders and society we operate in. Our core values and principles guide us in every day operations and strategic decision-making.

Headquartered in Harohalli, India, professionally managed Blueneem employs state-of-the-art manufacturing facilities and processes. Our products currently reach hospitals and patients in over 70 countries.

- ISO 13485 Certified Organization & CE Certified Products*
- Our Products are Manufactured in Best-in-class Infrastructure by Best-in-class People
- Paperless Manufacturing - 100% Digital Traceability
- World-class Clean Rooms - Class 7 and Class 8
- Bio Compatibility Tested and Certified Raw Materials
- Sterilization Validation and Confirmation Process

* Please check for product ranges



**BEST-IN-CLASS
INFRASTRUCTURE**



**WORLD-CLASS
CLEAN ROOMS**

PRODUCT LIST

GUIDEWIRES (5-12)

- PTFE Guidewire
- Twister Guidewire
- Hydro Twister™ Guidewire
- Tracer™ Guidewire
- Amplatz Extra Stiff Guidewire
- Amplatz Ultra Stiff™ Guidewire
- Prime Wire™

STENTS (13-23)

- P Flex Stent With String (Ureteral Stent / Double J Stent)
- P Flex Stent Set With String (Ureteral Stent Set / Double J Stent Set)
- Endopyelotomy Stent
- Urethral Stent

URETERAL ACCESS SHEATH (24-25)

- Praveza® Ureteral Access Sheath

CATHETERS (26-35)

- Ureteral Catheter
- Mono J Catheter
- PCN Catheter
- PCN Catheter with Locking Mechanism
- PCN Catheter with Needle
(Nephrostomy Drainage Set Direct Puncture)
- Malecot Catheter
- Re-entry Malecot Drainage Catheter
- Suprapubic Catheter (SPC)
- CI Catheter

DRAINAGE SETS (36-40)

- PCN Set - Pigtail
- PCN Set - Pigtail with Locking Mechanism
- PCN Set - Malecot
- Suprapubic Drainage Catheter Set (SDS)

STONE BASKETS (41-48)

- N Quadra™ - Nitinol Stone Basket With Tip
- N Quadra™ Z - Nitinol Stone Basket Zero Tip
- N Quadra™ Helical - Nitinol Stone Basket Helical With Tip
- N Quadra™ Helical Z - Nitinol Stone Basket Helical Zero Tip
- N Quadra™ Flexible - Nitinol Stone Basket Flexible
- Delta Basket
- PCNL Basket

PRODUCT LIST

DILATORS (49-63)

- Fascial Dilator Set
- Ureteral Dilator Set
- Single Step™ Ureteral Dilator
- Urethral Dilator Set
- Urethral Dilator Set - S Curved
- Urethral Self Dilator
- Amplatz Sheath
- Amplatz Renal Dilator
- Amplatz Renal Dilation Set - S
- Amplatz Renal Dilation Set - F
- Mini Perc Renal Dilation Set
- Ureteral Balloon Dilator
- Nephrostomy Balloon Dilator With Amplatz Sheath
- Screw Dilator

NEEDLES (64-66)

- Initial Puncture Needle - Two Part
- Sure Puncture™ Needle (Initial Puncture Needle - Three Part)

BIOPSY (67-68)

- Biopsy Gun - Semi Automatic

ACCESSORIES (69-75)

- Meatal Dilator - Male
- Dalela Female Self Urethral Dilator
- Bladder Evacuator
- Penile Clamp
- Irrigation Set (TUR Set)
- Cath Tap

GUIDEWIRES



- PTFE Guidewire
- Twister Guidewire
- Hydro Twister™ Guidewire
- Tracer™ Guidewire
- Amplatz Extra Stiff Guidewire
- Amplatz Ultra Stiff™ Guidewire
- Prime Wire™

UROLOGY

PTFE GUIDEWIRE

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures (URS).

Key Features

- PTFE Coating for smooth insertion
- Atraumatic floppy tip for less trauma
- Uniform outer diameter

Types


- J Tip
- Straight Tip

Size (in inch)

0.018, 0.025, 0.028, 0.032, 0.035
& 0.038

Length (in cm)

150



PTFE GUIDEWIRE - J TIP

| Product Code | Size (in inch) | Length (in cm) | Product Code | Size (in inch) | Length (in cm) |
|---------------|----------------|----------------|---------------|----------------|----------------|
| BPTS018150BJT | 0.018 | 150 | BPTS018150BST | 0.018 | 150 |
| BPTS025150BJT | 0.025 | 150 | BPTS025150BST | 0.025 | 150 |
| BPTS028150BJT | 0.028 | 150 | BPTS028150BST | 0.028 | 150 |
| BPTS032150BJT | 0.032 | 150 | BPTS032150BST | 0.032 | 150 |
| BPTS035150BJT | 0.035 | 150 | BPTS035150BST | 0.035 | 150 |
| BPTS038150BJT | 0.038 | 150 | BPTS038150BST | 0.038 | 150 |

PTFE GUIDEWIRE - STRAIGHT TIP

TWISTER GUIDEWIRE

Used to gain access to desired anatomical location and to assist in placement of medical devices (URS procedure, for ureteral stenting).

Key Features

- Is constructed with kink resistant nitinol core with appropriate stiffness for deployment of devices such as stents during a procedure
- Atraumatic floppy tip of 5 cm facilitates easy insertion and access
- Blue & White sleeve provides excellent visualization and easy maneuvering during diagnostic and interventional procedure

Type


Straight Tip

Size (in inch)

0.025, 0.032, 0.035 & 0.038

Length (in cm)

150



TWISTER GUIDEWIRE - STRAIGHT TIP

| Product Code | Size (in inch) | Length (in cm) |
|---------------|----------------|----------------|
| BTWN025150BST | 0.025 | 150 |
| BTWN032150BST | 0.032 | 150 |
| BTWN035150BST | 0.035 | 150 |
| BTWN038150BST | 0.038 | 150 |

HYDRO TWISTER™ GUIDEWIRE

Used in placement of medical devices while gaining easy access to desired anatomical location (URS & RIRS).

Key Features

- Hydrophilic coated tip and blue & white sleeve provides easy, low friction advancement with excellent visualization and maneuverability during diagnostic and interventional procedures
- Atraumatic floppy tip of 5 cm facilitates easy insertion and access
- Is constructed with kink resistant nitinol core with appropriate stiffness for deployment of devices such as stents during a procedure

Advantages

- Avoid exchange of wires
- Less procedural time
- Easy access in tortuous anatomy
- Serves two purpose - access and stenting with single wire

Type


Straight Tip

Size (in inch)

0.025, 0.032, 0.035 & 0.038

Length (in cm)

150



HYDRO TWISTER™ GUIDEWIRE - STRAIGHT TIP

| Product Code | Size (in inch) | Length (in cm) |
|---------------|----------------|----------------|
| BHTNO25150BST | 0.025 | 150 |
| BHTNO32150BST | 0.032 | 150 |
| BHTNO35150BST | 0.035 | 150 |
| BHTNO38150BST | 0.038 | 150 |

TRACER™ GUIDEWIRE (HYDROPHILIC GUIDEWIRE)

Used to gain smooth and easy access to desired anatomical location and to assist in placement of medical devices.

Key Features

- Lubricious hydrophilic coating provides easy passage in tortuous anatomy and through lesions
- Radiopaque polymer jacket enhances fluoroscopic visualization
- Kink resistant nitinol core for optimal performance

Types

- Straight Tip
- Angled Tip

Size (in inch)

0.018, 0.025, 0.032, 0.035 & 0.038

Length (in cm)

150



TRACER™ GUIDEWIRE - ANGLED TIP

| Product Code | Size (in inch) | Length (in cm) |
|---------------|----------------|----------------|
| BHWN018150BAT | 0.018 | 150 |
| BHWN025150BAT | 0.025 | 150 |
| BHWN032150BAT | 0.032 | 150 |
| BHWN035150BAT | 0.035 | 150 |
| BHWN038150BAT | 0.038 | 150 |

TRACER™ GUIDEWIRE - STRAIGHT TIP

| Product Code | Size (in inch) | Length (in cm) |
|---------------|----------------|----------------|
| BHWN018150BST | 0.018 | 150 |
| BHWN025150BST | 0.025 | 150 |
| BHWN032150BST | 0.032 | 150 |
| BHWN035150BST | 0.035 | 150 |
| BHWN038150BST | 0.038 | 150 |

AMPLATZ EXTRA STIFF GUIDEWIRE

Used to gain access to desired anatomical location and to assist in placement of medical devices and heavier instruments (PCN & PCNL).

Key Features

- PTFE coating for smooth insertion
- Atraumatic floppy tip for less trauma
- The stiffness of the wire provides support for the advancement of devices such as dilators, catheters and instruments

Type

- Straight Tip
- J Tip

Size (in inch)

0.032 & 0.035

Length (in cm)

80



AMPLATZ EXTRA STIFF GUIDEWIRE - J TIP

Product Code

Size
(in inch)

Length
(in cm)

BAES032080BJT
BAES035080BJT

0.032
0.035

80
80

AMPLATZ EXTRA STIFF GUIDEWIRE - STRAIGHT TIP

Product Code

Size
(in inch)

Length
(in cm)

BAES032080BST
BAES035080BST

0.032
0.035

80
80

AMPLATZ ULTRA STIFF™ GUIDEWIRE

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures (PCNL).

Key Features

- PTFE coated stainless steel mandril with 5 cm atraumatic floppy tip
- Smooth transition from stiff mandril to the flexible tip
- The stiffness of the mandril design provides support for the advancement of devices such as dilators, catheters and instruments

Type

J Tip

Size (in inch)

0.035

Length (in cm)

80



AMPLATZ ULTRA STIFF™ GUIDEWIRE - J TIP

| Product Code | Size (in inch) | Length (in cm) |
|---------------|----------------|----------------|
| BASSO35080BJT | 0.035 | 80 |

PRIME WIRE™

NEW LAUNCH (Gold Standard Guidewire for Urology)

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures.

Key Features

- Designed to combine the benefits of multiple specialty guidewires in one guidewire for urological procedures
- 5 cm hydrophilic floppy tip provides easy passage in tortuous anatomy
- Teflon coating provides for low friction navigation and easy advancement of devices

Advantages

- One wire for all urology procedures
- Serves multiple purposes - access and stenting with single wire
- Avoid exchange of wires
- Less procedural time
- Easy access in tortuous anatomy

Types


- Straight Tip
- Angled Tip

Size (in inch)

0.032, 0.035 & 0.038

Length (in cm)

150



PRIME WIRE™ - STRAIGHT TIP

| Product Code | Size (in inch) | Length (in cm) |
|---------------|----------------|----------------|
| BPTN032150BHS | 0.032 | 150 |
| BPTN035150BHS | 0.035 | 150 |
| BPTN038150BHS | 0.038 | 150 |

PRIME WIRE™ - ANGLED TIP

| Product Code | Size (in inch) | Length (in cm) |
|---------------|----------------|----------------|
| BPTN032150BHA | 0.032 | 150 |
| BPTN035150BHA | 0.035 | 150 |
| BPTN038150BHA | 0.038 | 150 |

STENTS

- P Flex Stent With String (Ureteral Stent / Double J Stent)
- P Flex Stent Set With String (Ureteral Stent Set / Double J Stent Set)
- Endopyelotomy Stent
- Urethral Stent

UROLOGY

P FLEX STENT WITH STRING

(URETERAL STENT / DOUBLE J STENT)

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder.

Key Features


- For indwelling up to 30 days
- Polyurethane based stent has uniform diameter and smooth surface reducing encrustation
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization
- Stent with string which helps in repositioning

Set Contains

Stent with string, pusher & clamp

Types

- Both Ends Open
- One End Closed
- Both Ends Multiloop
- With or Without Hydrophilic Coating



Standard

Size (in Fr)

4.7, 4.8, 5, 6, 7 & 8

Length (in cm)

12, 14, 16, 18, 20, 22,
24, 26 28 & 30

Multiloop

Size (in Fr)

4.7, 4.8, 5, 6, 7 & 8

Length (in cm)

18-24 & 24-30

Pediatric

Size (in Fr)

3 & 4

Length (in cm)

10, 12, 14 & 16

P FLEX STENT BOTH ENDS OPEN
WITH STRING

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUSB003010B1 | 3 | 10 |
| BUSB003012B1 | 3 | 12 |
| BUSB003014B1 | 3 | 14 |
| BUSB003016B1 | 3 | 16 |
| BUSB004010B1 | 4 | 10 |
| BUSB004012B1 | 4 | 12 |
| BUSB004014B1 | 4 | 14 |
| BUSB004016B1 | 4 | 16 |
| BUSB047012B1 | 4.7 | 12 |
| BUSB047014B1 | 4.7 | 14 |
| BUSB047016B1 | 4.7 | 16 |
| BUSB047018B1 | 4.7 | 18 |
| BUSB047020B1 | 4.7 | 20 |
| BUSB047022B1 | 4.7 | 22 |
| BUSB047024B1 | 4.7 | 24 |
| BUSB047026B1 | 4.7 | 26 |
| BUSB047028B1 | 4.7 | 28 |
| BUSB047030B1 | 4.7 | 30 |
| BUSB048012B1 | 4.8 | 12 |
| BUSB048014B1 | 4.8 | 14 |
| BUSB048016B1 | 4.8 | 16 |
| BUSB048018B1 | 4.8 | 18 |
| BUSB048020B1 | 4.8 | 20 |
| BUSB048022B1 | 4.8 | 22 |
| BUSB048024B1 | 4.8 | 24 |
| BUSB048026B1 | 4.8 | 26 |
| BUSB048028B1 | 4.8 | 28 |
| BUSB048030B1 | 4.8 | 30 |
| BUSB005012B1 | 5 | 12 |
| BUSB005014B1 | 5 | 14 |
| BUSB005016B1 | 5 | 16 |
| BUSB005018B1 | 5 | 18 |
| BUSB005020B1 | 5 | 20 |
| BUSB005022B1 | 5 | 22 |
| BUSB005024B1 | 5 | 24 |
| BUSB005026B1 | 5 | 26 |
| BUSB005028B1 | 5 | 28 |
| BUSB005030B1 | 5 | 30 |
| BUSB006012B1 | 6 | 12 |
| BUSB006014B1 | 6 | 14 |
| BUSB006016B1 | 6 | 16 |
| BUSB006018B1 | 6 | 18 |
| BUSB006020B1 | 6 | 20 |
| BUSB006022B1 | 6 | 22 |
| BUSB006024B1 | 6 | 24 |
| BUSB006026B1 | 6 | 26 |
| BUSB006028B1 | 6 | 28 |
| BUSB006030B1 | 6 | 30 |
| BUSB007012B1 | 7 | 12 |
| BUSB007014B1 | 7 | 14 |
| BUSB007016B1 | 7 | 16 |
| BUSB007018B1 | 7 | 18 |
| BUSB007020B1 | 7 | 20 |
| BUSB007022B1 | 7 | 22 |
| BUSB007024B1 | 7 | 24 |
| BUSB007026B1 | 7 | 26 |
| BUSB007028B1 | 7 | 28 |
| BUSB007030B1 | 7 | 30 |
| BUSB008012B1 | 8 | 12 |
| BUSB008014B1 | 8 | 14 |
| BUSB008016B1 | 8 | 16 |
| BUSB008018B1 | 8 | 18 |
| BUSB008020B1 | 8 | 20 |
| BUSB008022B1 | 8 | 22 |
| BUSB008024B1 | 8 | 24 |
| BUSB008026B1 | 8 | 26 |
| BUSB008028B1 | 8 | 28 |
| BUSB008030B1 | 8 | 30 |

P FLEX STENT ONE END CLOSED
WITH STRING

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUSB003010B2 | 3 | 10 |
| BUSB003012B2 | 3 | 12 |
| BUSB003014B2 | 3 | 14 |
| BUSB003016B2 | 3 | 16 |
| BUSB004010B2 | 4 | 10 |
| BUSB004012B2 | 4 | 12 |
| BUSB004014B2 | 4 | 14 |
| BUSB004016B2 | 4 | 16 |
| BUSB047012B2 | 4.7 | 12 |
| BUSB047014B2 | 4.7 | 14 |
| BUSB047016B2 | 4.7 | 16 |
| BUSB047018B2 | 4.7 | 18 |
| BUSB047020B2 | 4.7 | 20 |
| BUSB047022B2 | 4.7 | 22 |
| BUSB047024B2 | 4.7 | 24 |
| BUSB047026B2 | 4.7 | 26 |
| BUSB047028B2 | 4.7 | 28 |
| BUSB047030B2 | 4.7 | 30 |
| BUSB048012B2 | 4.8 | 12 |
| BUSB048014B2 | 4.8 | 14 |
| BUSB048016B2 | 4.8 | 16 |
| BUSB048018B2 | 4.8 | 18 |
| BUSB048020B2 | 4.8 | 20 |
| BUSB048022B2 | 4.8 | 22 |
| BUSB048024B2 | 4.8 | 24 |
| BUSB048026B2 | 4.8 | 26 |
| BUSB048028B2 | 4.8 | 28 |
| BUSB048030B2 | 4.8 | 30 |
| BUSB005012B2 | 5 | 12 |
| BUSB005014B2 | 5 | 14 |
| BUSB005016B2 | 5 | 16 |
| BUSB005018B2 | 5 | 18 |
| BUSB005020B2 | 5 | 20 |
| BUSB005022B2 | 5 | 22 |
| BUSB005024B2 | 5 | 24 |
| BUSB005026B2 | 5 | 26 |
| BUSB005028B2 | 5 | 28 |
| BUSB005030B2 | 5 | 30 |
| BUSB006012B2 | 6 | 12 |
| BUSB006014B2 | 6 | 14 |
| BUSB006016B2 | 6 | 16 |
| BUSB006018B2 | 6 | 18 |
| BUSB006020B2 | 6 | 20 |
| BUSB006022B2 | 6 | 22 |
| BUSB006024B2 | 6 | 24 |
| BUSB006026B2 | 6 | 26 |
| BUSB006028B2 | 6 | 28 |
| BUSB006030B2 | 6 | 30 |
| BUSB007012B2 | 7 | 12 |
| BUSB007014B2 | 7 | 14 |
| BUSB007016B2 | 7 | 16 |
| BUSB007018B2 | 7 | 18 |
| BUSB007020B2 | 7 | 20 |
| BUSB007022B2 | 7 | 22 |
| BUSB007024B2 | 7 | 24 |
| BUSB007026B2 | 7 | 26 |
| BUSB007028B2 | 7 | 28 |
| BUSB007030B2 | 7 | 30 |
| BUSB008012B2 | 8 | 12 |
| BUSB008014B2 | 8 | 14 |
| BUSB008016B2 | 8 | 16 |
| BUSB008018B2 | 8 | 18 |
| BUSB008020B2 | 8 | 20 |
| BUSB008022B2 | 8 | 22 |
| BUSB008024B2 | 8 | 24 |
| BUSB008026B2 | 8 | 26 |
| BUSB008028B2 | 8 | 28 |
| BUSB008030B2 | 8 | 30 |

**P FLEX STENT BOTH ENDS OPEN
WITH STRING & HYDROPHILIC COATING**

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUSB003010B3 | 3 | 10 |
| BUSB003012B3 | 3 | 12 |
| BUSB003014B3 | 3 | 14 |
| BUSB003016B3 | 3 | 16 |
| BUSB004010B3 | 4 | 10 |
| BUSB004012B3 | 4 | 12 |
| BUSB004014B3 | 4 | 14 |
| BUSB004016B3 | 4 | 16 |
| BUSB047012B3 | 4.7 | 12 |
| BUSB047014B3 | 4.7 | 14 |
| BUSB047016B3 | 4.7 | 16 |
| BUSB047018B3 | 4.7 | 18 |
| BUSB047020B3 | 4.7 | 20 |
| BUSB047022B3 | 4.7 | 22 |
| BUSB047024B3 | 4.7 | 24 |
| BUSB047026B3 | 4.7 | 26 |
| BUSB047028B3 | 4.7 | 28 |
| BUSB047030B3 | 4.7 | 30 |
| BUSB048012B3 | 4.8 | 12 |
| BUSB048014B3 | 4.8 | 14 |
| BUSB048016B3 | 4.8 | 16 |
| BUSB048018B3 | 4.8 | 18 |
| BUSB048020B3 | 4.8 | 20 |
| BUSB048022B3 | 4.8 | 22 |
| BUSB048024B3 | 4.8 | 24 |
| BUSB048026B3 | 4.8 | 26 |
| BUSB048028B3 | 4.8 | 28 |
| BUSB048030B3 | 4.8 | 30 |
| BUSB005012B3 | 5 | 12 |
| BUSB005014B3 | 5 | 14 |
| BUSB005016B3 | 5 | 16 |
| BUSB005018B3 | 5 | 18 |
| BUSB005020B3 | 5 | 20 |
| BUSB005022B3 | 5 | 22 |
| BUSB005024B3 | 5 | 24 |
| BUSB005026B3 | 5 | 26 |
| BUSB005028B3 | 5 | 28 |
| BUSB005030B3 | 5 | 30 |
| BUSB006012B3 | 6 | 12 |
| BUSB006014B3 | 6 | 14 |
| BUSB006016B3 | 6 | 16 |
| BUSB006018B3 | 6 | 18 |
| BUSB006020B3 | 6 | 20 |
| BUSB006022B3 | 6 | 22 |
| BUSB006024B3 | 6 | 24 |
| BUSB006026B3 | 6 | 26 |
| BUSB006028B3 | 6 | 28 |
| BUSB006030B3 | 6 | 30 |
| BUSB007012B3 | 7 | 12 |
| BUSB007014B3 | 7 | 14 |
| BUSB007016B3 | 7 | 16 |
| BUSB007018B3 | 7 | 18 |
| BUSB007020B3 | 7 | 20 |
| BUSB007022B3 | 7 | 22 |
| BUSB007024B3 | 7 | 24 |
| BUSB007026B3 | 7 | 26 |
| BUSB007028B3 | 7 | 28 |
| BUSB007030B3 | 7 | 30 |
| BUSB008012B3 | 8 | 12 |
| BUSB008014B3 | 8 | 14 |
| BUSB008016B3 | 8 | 16 |
| BUSB008018B3 | 8 | 18 |
| BUSB008020B3 | 8 | 20 |
| BUSB008022B3 | 8 | 22 |
| BUSB008024B3 | 8 | 24 |
| BUSB008026B3 | 8 | 26 |
| BUSB008028B3 | 8 | 28 |
| BUSB008030B3 | 8 | 30 |

**P FLEX STENT ONE END CLOSED
WITH STRING & HYDROPHILIC COATING**

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUSB003010B4 | 3 | 10 |
| BUSB003012B4 | 3 | 12 |
| BUSB003014B4 | 3 | 14 |
| BUSB003016B4 | 3 | 16 |
| BUSB004010B4 | 4 | 10 |
| BUSB004012B4 | 4 | 12 |
| BUSB004014B4 | 4 | 14 |
| BUSB004016B4 | 4 | 16 |
| BUSB047012B4 | 4.7 | 12 |
| BUSB047014B4 | 4.7 | 14 |
| BUSB047016B4 | 4.7 | 16 |
| BUSB047018B4 | 4.7 | 18 |
| BUSB047020B4 | 4.7 | 20 |
| BUSB047022B4 | 4.7 | 22 |
| BUSB047024B4 | 4.7 | 24 |
| BUSB047026B4 | 4.7 | 26 |
| BUSB047028B4 | 4.7 | 28 |
| BUSB047030B4 | 4.7 | 30 |
| BUSB048012B4 | 4.8 | 12 |
| BUSB048014B4 | 4.8 | 14 |
| BUSB048016B4 | 4.8 | 16 |
| BUSB048018B4 | 4.8 | 18 |
| BUSB048020B4 | 4.8 | 20 |
| BUSB048022B4 | 4.8 | 22 |
| BUSB048024B4 | 4.8 | 24 |
| BUSB048026B4 | 4.8 | 26 |
| BUSB048028B4 | 4.8 | 28 |
| BUSB048030B4 | 4.8 | 30 |
| BUSB005012B4 | 5 | 12 |
| BUSB005014B4 | 5 | 14 |
| BUSB005016B4 | 5 | 16 |
| BUSB005018B4 | 5 | 18 |
| BUSB005020B4 | 5 | 20 |
| BUSB005022B4 | 5 | 22 |
| BUSB005024B4 | 5 | 24 |
| BUSB005026B4 | 5 | 26 |
| BUSB005028B4 | 5 | 28 |
| BUSB005030B4 | 5 | 30 |
| BUSB006012B4 | 6 | 12 |
| BUSB006014B4 | 6 | 14 |
| BUSB006016B4 | 6 | 16 |
| BUSB006018B4 | 6 | 18 |
| BUSB006020B4 | 6 | 20 |
| BUSB006022B4 | 6 | 22 |
| BUSB006024B4 | 6 | 24 |
| BUSB006026B4 | 6 | 26 |
| BUSB006028B4 | 6 | 28 |
| BUSB006030B4 | 6 | 30 |
| BUSB007012B4 | 7 | 12 |
| BUSB007014B4 | 7 | 14 |
| BUSB007016B4 | 7 | 16 |
| BUSB007018B4 | 7 | 18 |
| BUSB007020B4 | 7 | 20 |
| BUSB007022B4 | 7 | 22 |
| BUSB007024B4 | 7 | 24 |
| BUSB007026B4 | 7 | 26 |
| BUSB007028B4 | 7 | 28 |
| BUSB007030B4 | 7 | 30 |
| BUSB008012B4 | 8 | 12 |
| BUSB008014B4 | 8 | 14 |
| BUSB008016B4 | 8 | 16 |
| BUSB008018B4 | 8 | 18 |
| BUSB008020B4 | 8 | 20 |
| BUSB008022B4 | 8 | 22 |
| BUSB008024B4 | 8 | 24 |
| BUSB008026B4 | 8 | 26 |
| BUSB008028B4 | 8 | 28 |
| BUSB008030B4 | 8 | 30 |

| P FLEX STENT BOTH ENDS OPEN BOTH ENDS MULTILOOP WITH STRING | | | P FLEX STENT ONE END CLOSED BOTH ENDS MULTILOOP WITH STRING | | |
|--|--------------|----------------|--|--------------|----------------|
| Product Code | Size (in Fr) | Length (in cm) | Product Code | Size (in Fr) | Length (in cm) |
| BUSB004022B5 | 4 | 22-32 | BUSB004022B6 | 4 | 22-32 |
| BUSB047018B5 BUSB047024B5 | 4.7 4.7 | 18 24 | BUSB047018B6 BUSB047024B6 | 4.7 4.7 | 18 24 |
| BUSB048018B5 BUSB048024B5 | 4.8 4.8 | 18 24 | BUSB048018B6 BUSB048024B6 | 4.8 4.8 | 18 24 |
| BUSB005018B5 BUSB005024B5 | 5 5 | 18 24 | BUSB005018B6 BUSB005024B6 | 5 5 | 18 24 |
| BUSB006018B5 BUSB006024B5 | 6 6 | 18 24 | BUSB006018B6 BUSB006024B6 | 6 6 | 18 24 |
| BUSB007018B5 BUSB007024B5 | 7 7 | 18 24 | BUSB007018B6 BUSB007024B6 | 7 7 | 18 24 |
| BUSB008018B5 BUSB008024B5 | 8 8 | 18 24 | BUSB008018B6 BUSB008024B6 | 8 8 | 18 24 |

| P FLEX STENT BOTH ENDS OPEN BEML WITH STRING & HYDROPHILIC COATING | | | P FLEX STENT ONE END CLOSED BEML WITH STRING & HYDROPHILIC COATING | | |
|---|--------------|----------------|---|--------------|----------------|
| Product Code | Size (in Fr) | Length (in cm) | Product Code | Size (in Fr) | Length (in cm) |
| BUSB047018B7 BUSB047024B7 | 4.7 4.7 | 18 24 | BUSB047018B8 BUSB047024B8 | 4.7 4.7 | 18 24 |
| BUSB048018B7 BUSB048024B7 | 4.8 4.8 | 18 24 | BUSB048018B8 BUSB048024B8 | 4.8 4.8 | 18 24 |
| BUSB005018B7 BUSB005024B7 | 5 5 | 18 24 | BUSB005018B8 BUSB005024B8 | 5 5 | 18 24 |
| BUSB006018B7 BUSB006024B7 | 6 6 | 18 24 | BUSB006018B8 BUSB006024B8 | 6 6 | 18 24 |
| BUSB007018B7 BUSB007024B7 | 7 7 | 18 24 | BUSB007018B8 BUSB007024B8 | 7 7 | 18 24 |
| BUSB008018B7 BUSB008024B7 | 8 8 | 18 24 | BUSB008018B8 BUSB008024B8 | 8 8 | 18 24 |

P FLEX STENT SET WITH STRING

(URETERAL STENT SET / DOUBLE J STENT SET)

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder.

Key Features


- For indwelling up to 30 days
- Polyurethane based stent has uniform diameter and smooth surface reducing encrustation
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization
- Stent with string which helps in repositioning

Set Contains

Stent with string, pusher, guidewire & clamp

Types

- Both Ends Open
- One End Closed
- Both Ends Multiloop
- With or Without Hydrophilic Coating



Standard

| Size (in Fr) | Length (in cm) |
|-----------------------|---|
| 4.7, 4.8, 5, 6, 7 & 8 | 12, 14, 16, 18, 20, 22, 24, 26, 28 & 30 |

Multiloop

| Size (in Fr) | Length (in cm) |
|-----------------------|----------------|
| 4.7, 4.8, 5, 6, 7 & 8 | 18-24 & 24-30 |

Pediatric

| Size (in Fr) | Length (in cm) |
|--------------|-----------------|
| 3 & 4 | 10, 12, 14 & 16 |

P FLEX STENT SET BOTH ENDS OPEN
WITH STRING

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUSB003010B9 | 3 | 10 |
| BUSB003012B9 | 3 | 12 |
| BUSB003014B9 | 3 | 14 |
| BUSB003016B9 | 3 | 16 |
| BUSB004010B9 | 4 | 10 |
| BUSB004012B9 | 4 | 12 |
| BUSB004014B9 | 4 | 14 |
| BUSB004016B9 | 4 | 16 |
| BUSB047012B9 | 4.7 | 12 |
| BUSB047014B9 | 4.7 | 14 |
| BUSB047016B9 | 4.7 | 16 |
| BUSB047018B9 | 4.7 | 18 |
| BUSB047020B9 | 4.7 | 20 |
| BUSB047022B9 | 4.7 | 22 |
| BUSB047024B9 | 4.7 | 24 |
| BUSB047026B9 | 4.7 | 26 |
| BUSB047028B9 | 4.7 | 28 |
| BUSB047030B9 | 4.7 | 30 |
| BUSB048012B9 | 4.8 | 12 |
| BUSB048014B9 | 4.8 | 14 |
| BUSB048016B9 | 4.8 | 16 |
| BUSB048018B9 | 4.8 | 18 |
| BUSB048020B9 | 4.8 | 20 |
| BUSB048022B9 | 4.8 | 22 |
| BUSB048024B9 | 4.8 | 24 |
| BUSB048026B9 | 4.8 | 26 |
| BUSB048028B9 | 4.8 | 28 |
| BUSB048030B9 | 4.8 | 30 |
| BUSB005012B9 | 5 | 12 |
| BUSB005014B9 | 5 | 14 |
| BUSB005016B9 | 5 | 16 |
| BUSB005018B9 | 5 | 18 |
| BUSB005020B9 | 5 | 20 |
| BUSB005022B9 | 5 | 22 |
| BUSB005024B9 | 5 | 24 |
| BUSB005026B9 | 5 | 26 |
| BUSB005028B9 | 5 | 28 |
| BUSB005030B9 | 5 | 30 |
| BUSB006012B9 | 6 | 12 |
| BUSB006014B9 | 6 | 14 |
| BUSB006016B9 | 6 | 16 |
| BUSB006018B9 | 6 | 18 |
| BUSB006020B9 | 6 | 20 |
| BUSB006022B9 | 6 | 22 |
| BUSB006024B9 | 6 | 24 |
| BUSB006026B9 | 6 | 26 |
| BUSB006028B9 | 6 | 28 |
| BUSB006030B9 | 6 | 30 |
| BUSB007012B9 | 7 | 12 |
| BUSB007014B9 | 7 | 14 |
| BUSB007016B9 | 7 | 16 |
| BUSB007018B9 | 7 | 18 |
| BUSB007020B9 | 7 | 20 |
| BUSB007022B9 | 7 | 22 |
| BUSB007024B9 | 7 | 24 |
| BUSB007026B9 | 7 | 26 |
| BUSB007028B9 | 7 | 28 |
| BUSB007030B9 | 7 | 30 |
| BUSB008012B9 | 8 | 12 |
| BUSB008014B9 | 8 | 14 |
| BUSB008016B9 | 8 | 16 |
| BUSB008018B9 | 8 | 18 |
| BUSB008020B9 | 8 | 20 |
| BUSB008022B9 | 8 | 22 |
| BUSB008024B9 | 8 | 24 |
| BUSB008026B9 | 8 | 26 |
| BUSB008028B9 | 8 | 28 |
| BUSB008030B9 | 8 | 30 |

P FLEX STENT SET ONE END CLOSED
WITH STRING

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BUSB003010B10 | 3 | 10 |
| BUSB003012B10 | 3 | 12 |
| BUSB003014B10 | 3 | 14 |
| BUSB003016B10 | 3 | 16 |
| BUSB004010B10 | 4 | 10 |
| BUSB004012B10 | 4 | 12 |
| BUSB004014B10 | 4 | 14 |
| BUSB004016B10 | 4 | 16 |
| BUSB047012B10 | 4.7 | 12 |
| BUSB047014B10 | 4.7 | 14 |
| BUSB047016B10 | 4.7 | 16 |
| BUSB047018B10 | 4.7 | 18 |
| BUSB047020B10 | 4.7 | 20 |
| BUSB047022B10 | 4.7 | 22 |
| BUSB047024B10 | 4.7 | 24 |
| BUSB047026B10 | 4.7 | 26 |
| BUSB047028B10 | 4.7 | 28 |
| BUSB047030B10 | 4.7 | 30 |
| BUSB048012B10 | 4.8 | 12 |
| BUSB048014B10 | 4.8 | 14 |
| BUSB048016B10 | 4.8 | 16 |
| BUSB048018B10 | 4.8 | 18 |
| BUSB048020B10 | 4.8 | 20 |
| BUSB048022B10 | 4.8 | 22 |
| BUSB048024B10 | 4.8 | 24 |
| BUSB048026B10 | 4.8 | 26 |
| BUSB048028B10 | 4.8 | 28 |
| BUSB048030B10 | 4.8 | 30 |
| BUSB005012B10 | 5 | 12 |
| BUSB005014B10 | 5 | 14 |
| BUSB005016B10 | 5 | 16 |
| BUSB005018B10 | 5 | 18 |
| BUSB005020B10 | 5 | 20 |
| BUSB005022B10 | 5 | 22 |
| BUSB005024B10 | 5 | 24 |
| BUSB005026B10 | 5 | 26 |
| BUSB005028B10 | 5 | 28 |
| BUSB005030B10 | 5 | 30 |
| BUSB006012B10 | 6 | 12 |
| BUSB006014B10 | 6 | 14 |
| BUSB006016B10 | 6 | 16 |
| BUSB006018B10 | 6 | 18 |
| BUSB006020B10 | 6 | 20 |
| BUSB006022B10 | 6 | 22 |
| BUSB006024B10 | 6 | 24 |
| BUSB006026B10 | 6 | 26 |
| BUSB006028B10 | 6 | 28 |
| BUSB006030B10 | 6 | 30 |
| BUSB007012B10 | 7 | 12 |
| BUSB007014B10 | 7 | 14 |
| BUSB007016B10 | 7 | 16 |
| BUSB007018B10 | 7 | 18 |
| BUSB007020B10 | 7 | 20 |
| BUSB007022B10 | 7 | 22 |
| BUSB007024B10 | 7 | 24 |
| BUSB007026B10 | 7 | 26 |
| BUSB007028B10 | 7 | 28 |
| BUSB007030B10 | 7 | 30 |
| BUSB008012B10 | 8 | 12 |
| BUSB008014B10 | 8 | 14 |
| BUSB008016B10 | 8 | 16 |
| BUSB008018B10 | 8 | 18 |
| BUSB008020B10 | 8 | 20 |
| BUSB008022B10 | 8 | 22 |
| BUSB008024B10 | 8 | 24 |
| BUSB008026B10 | 8 | 26 |
| BUSB008028B10 | 8 | 28 |
| BUSB008030B10 | 8 | 30 |

**P FLEX STENT SET BOTH ENDS OPEN
WITH STRING & HYDROPHILIC COATING**

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BUSB003010B11 | 3 | 10 |
| BUSB003012B11 | 3 | 12 |
| BUSB003014B11 | 3 | 14 |
| BUSB003016B11 | 3 | 16 |
| BUSB004010B11 | 4 | 10 |
| BUSB004012B11 | 4 | 12 |
| BUSB004014B11 | 4 | 14 |
| BUSB004016B11 | 4 | 16 |
| BUSB047012B11 | 4.7 | 12 |
| BUSB047014B11 | 4.7 | 14 |
| BUSB047016B11 | 4.7 | 16 |
| BUSB047018B11 | 4.7 | 18 |
| BUSB047020B11 | 4.7 | 20 |
| BUSB047022B11 | 4.7 | 22 |
| BUSB047024B11 | 4.7 | 24 |
| BUSB047026B11 | 4.7 | 26 |
| BUSB047028B11 | 4.7 | 28 |
| BUSB047030B11 | 4.7 | 30 |
| BUSB048012B11 | 4.8 | 12 |
| BUSB048014B11 | 4.8 | 14 |
| BUSB048016B11 | 4.8 | 16 |
| BUSB048018B11 | 4.8 | 18 |
| BUSB048020B11 | 4.8 | 20 |
| BUSB048022B11 | 4.8 | 22 |
| BUSB048024B11 | 4.8 | 24 |
| BUSB048026B11 | 4.8 | 26 |
| BUSB048028B11 | 4.8 | 28 |
| BUSB048030B11 | 4.8 | 30 |
| BUSB005012B11 | 5 | 12 |
| BUSB005014B11 | 5 | 14 |
| BUSB005016B11 | 5 | 16 |
| BUSB005018B11 | 5 | 18 |
| BUSB005020B11 | 5 | 20 |
| BUSB005022B11 | 5 | 22 |
| BUSB005024B11 | 5 | 24 |
| BUSB005026B11 | 5 | 26 |
| BUSB005028B11 | 5 | 28 |
| BUSB005030B11 | 5 | 30 |
| BUSB006012B11 | 6 | 12 |
| BUSB006014B11 | 6 | 14 |
| BUSB006016B11 | 6 | 16 |
| BUSB006018B11 | 6 | 18 |
| BUSB006020B11 | 6 | 20 |
| BUSB006022B11 | 6 | 22 |
| BUSB006024B11 | 6 | 24 |
| BUSB006026B11 | 6 | 26 |
| BUSB006028B11 | 6 | 28 |
| BUSB006030B11 | 6 | 30 |
| BUSB007012B11 | 7 | 12 |
| BUSB007014B11 | 7 | 14 |
| BUSB007016B11 | 7 | 16 |
| BUSB007018B11 | 7 | 18 |
| BUSB007020B11 | 7 | 20 |
| BUSB007022B11 | 7 | 22 |
| BUSB007024B11 | 7 | 24 |
| BUSB007026B11 | 7 | 26 |
| BUSB007028B11 | 7 | 28 |
| BUSB007030B11 | 7 | 30 |
| BUSB008012B11 | 8 | 12 |
| BUSB008014B11 | 8 | 14 |
| BUSB008016B11 | 8 | 16 |
| BUSB008018B11 | 8 | 18 |
| BUSB008020B11 | 8 | 20 |
| BUSB008022B11 | 8 | 22 |
| BUSB008024B11 | 8 | 24 |
| BUSB008026B11 | 8 | 26 |
| BUSB008028B11 | 8 | 28 |
| BUSB008030B11 | 8 | 30 |

**P FLEX STENT SET ONE END CLOSED
WITH STRING & HYDROPHILIC COATING**

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BUSB003010B12 | 3 | 10 |
| BUSB003012B12 | 3 | 12 |
| BUSB003014B12 | 3 | 14 |
| BUSB003016B12 | 3 | 16 |
| BUSB004010B12 | 4 | 10 |
| BUSB004012B12 | 4 | 12 |
| BUSB004014B12 | 4 | 14 |
| BUSB004016B12 | 4 | 16 |
| BUSB047012B12 | 4.7 | 12 |
| BUSB047014B12 | 4.7 | 14 |
| BUSB047016B12 | 4.7 | 16 |
| BUSB047018B12 | 4.7 | 18 |
| BUSB047020B12 | 4.7 | 20 |
| BUSB047022B12 | 4.7 | 22 |
| BUSB047024B12 | 4.7 | 24 |
| BUSB047026B12 | 4.7 | 26 |
| BUSB047028B12 | 4.7 | 28 |
| BUSB047030B12 | 4.7 | 30 |
| BUSB048012B12 | 4.8 | 12 |
| BUSB048014B12 | 4.8 | 14 |
| BUSB048016B12 | 4.8 | 16 |
| BUSB048018B12 | 4.8 | 18 |
| BUSB048020B12 | 4.8 | 20 |
| BUSB048022B12 | 4.8 | 22 |
| BUSB048024B12 | 4.8 | 24 |
| BUSB048026B12 | 4.8 | 26 |
| BUSB048028B12 | 4.8 | 28 |
| BUSB048030B12 | 4.8 | 30 |
| BUSB005012B12 | 5 | 12 |
| BUSB005014B12 | 5 | 14 |
| BUSB005016B12 | 5 | 16 |
| BUSB005018B12 | 5 | 18 |
| BUSB005020B12 | 5 | 20 |
| BUSB005022B12 | 5 | 22 |
| BUSB005024B12 | 5 | 24 |
| BUSB005026B12 | 5 | 26 |
| BUSB005028B12 | 5 | 28 |
| BUSB005030B12 | 5 | 30 |
| BUSB006012B12 | 6 | 12 |
| BUSB006014B12 | 6 | 14 |
| BUSB006016B12 | 6 | 16 |
| BUSB006018B12 | 6 | 18 |
| BUSB006020B12 | 6 | 20 |
| BUSB006022B12 | 6 | 22 |
| BUSB006024B12 | 6 | 24 |
| BUSB006026B12 | 6 | 26 |
| BUSB006028B12 | 6 | 28 |
| BUSB006030B12 | 6 | 30 |
| BUSB007012B12 | 7 | 12 |
| BUSB007014B12 | 7 | 14 |
| BUSB007016B12 | 7 | 16 |
| BUSB007018B12 | 7 | 18 |
| BUSB007020B12 | 7 | 20 |
| BUSB007022B12 | 7 | 22 |
| BUSB007024B12 | 7 | 24 |
| BUSB007026B12 | 7 | 26 |
| BUSB007028B12 | 7 | 28 |
| BUSB007030B12 | 7 | 30 |
| BUSB008012B12 | 8 | 12 |
| BUSB008014B12 | 8 | 14 |
| BUSB008016B12 | 8 | 16 |
| BUSB008018B12 | 8 | 18 |
| BUSB008020B12 | 8 | 20 |
| BUSB008022B12 | 8 | 22 |
| BUSB008024B12 | 8 | 24 |
| BUSB008026B12 | 8 | 26 |
| BUSB008028B12 | 8 | 28 |
| BUSB008030B12 | 8 | 30 |

| P FLEX STENT SET BOTH ENDS OPEN BEML WITH STRING | | | P FLEX STENT SET ONE END CLOSED BEML WITH STRING | | |
|---|--------------|----------------|---|--------------|----------------|
| Product Code | Size (in Fr) | Length (in cm) | Product Code | Size (in Fr) | Length (in cm) |
| BUSB004022B13 | 4 | 22-32 | BUSB004022B14 | 4 | 22-32 |
| BUSB047018B13 | 4.7 | 18 | BUSB047018B14 | 4.7 | 18 |
| BUSB047024B13 | 4.7 | 24 | BUSB047024B14 | 4.7 | 24 |
| BUSB048018B13 | 4.8 | 18 | BUSB048018B14 | 4.8 | 18 |
| BUSB048024B13 | 4.8 | 24 | BUSB048024B14 | 4.8 | 24 |
| BUSB005018B13 | 5 | 18 | BUSB005018B14 | 5 | 18 |
| BUSB005024B13 | 5 | 24 | BUSB005024B14 | 5 | 24 |
| BUSB006018B13 | 6 | 18 | BUSB006018B14 | 6 | 18 |
| BUSB006024B13 | 6 | 24 | BUSB006024B14 | 6 | 24 |
| BUSB007018B13 | 7 | 18 | BUSB007018B14 | 7 | 18 |
| BUSB007024B13 | 7 | 24 | BUSB007024B14 | 7 | 24 |
| BUSB008018B13 | 8 | 18 | BUSB008018B14 | 8 | 18 |
| BUSB008024B13 | 8 | 24 | BUSB008024B14 | 8 | 24 |

| P FLEX STENT SET BOTH ENDS OPEN BEML WITH STRING & HYDROPHILIC COATING | | | P FLEX STENT SET ONE END CLOSED BEML WITH STRING & HYDROPHILIC COATING | | |
|---|--------------|----------------|---|--------------|----------------|
| Product Code | Size (in Fr) | Length (in cm) | Product Code | Size (in Fr) | Length (in cm) |
| BUSB047018B15 | 4.7 | 18 | BUSB047018B16 | 4.7 | 18 |
| BUSB047024B15 | 4.7 | 24 | BUSB047024B16 | 4.7 | 24 |
| BUSB048018B15 | 4.8 | 18 | BUSB048018B16 | 4.8 | 18 |
| BUSB048024B15 | 4.8 | 24 | BUSB048024B16 | 4.8 | 24 |
| BUSB005018B15 | 5 | 18 | BUSB005018B16 | 5 | 18 |
| BUSB005024B15 | 5 | 24 | BUSB005024B16 | 5 | 24 |
| BUSB006018B15 | 6 | 18 | BUSB006018B16 | 6 | 18 |
| BUSB006024B15 | 6 | 24 | BUSB006024B16 | 6 | 24 |
| BUSB007018B15 | 7 | 18 | BUSB007018B16 | 7 | 18 |
| BUSB007024B15 | 7 | 24 | BUSB007024B16 | 7 | 24 |
| BUSB008018B15 | 8 | 18 | BUSB008018B16 | 8 | 18 |
| BUSB008024B15 | 8 | 24 | BUSB008024B16 | 8 | 24 |

ENDOPYELOTOMY STENT

Used for internal drainage from uretero pelvic junction to the bladder following incision of a stricture.

Key Features

- Dual diameter design
- Smooth transition from larger diameter to smaller diameter
- Maintain ureteral patency post Endopyelotomy

Also available in Set, which Includes


Stent, pusher, guidewire & clamp

Size (in Fr)

6/12 & 7/14

Length (in cm)

24 & 26



ENDOPYELOTOMY STENT

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BEPB612024B1 | 6/12 | 24 |
| BEPB612026B1 | 6/12 | 26 |
| BEPB714024B1 | 7/14 | 24 |
| BEPB714026B1 | 7/14 | 26 |

ENDOPYELOTOMY STENT SET

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BEPB612024B2 | 6/12 | 24 |
| BEPB612026B2 | 6/12 | 26 |
| BEPB714024B2 | 7/14 | 24 |
| BEPB714026B2 | 7/14 | 26 |

URETHRAL STENT

Used for stenting the urethra after Hypospedias or Epispedias repair and to provide post-operative drainage of the bladder.

Key Features

- Smooth surface for ease of introduction
- Radiopaque for better visualization
- Hydrophilic coating for easy insertion

Size (in Fr)

5, 6, 7, 8, 10 & 12

Length (in cm)

25



URETHRAL STENT OPEN END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BHSB005025B1 | 5 | 25 |
| BHSB006025B1 | 6 | 25 |
| BHSB007025B1 | 7 | 25 |
| BHSB008025B1 | 8 | 25 |
| BHSB010025B1 | 10 | 25 |
| BHSB012025B1 | 12 | 25 |

URETHRAL STENT CLOSED END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BHSB005025B2 | 5 | 25 |
| BHSB006025B2 | 6 | 25 |
| BHSB007025B2 | 7 | 25 |
| BHSB008025B2 | 8 | 25 |
| BHSB010025B2 | 10 | 25 |
| BHSB012025B2 | 12 | 25 |

URETERAL ACCESS SHEATH

- Praveza[®] Ureteral Access Sheath

UROLOGY

PRAVEZA[®] URETERAL ACCESS SHEATH

Used to establish a conduit during urological procedures facilitating the passage for Ureteroscope / Endoscope and other surgical instruments into the urinary tract (RIRS).

Key Features

- Kink resistant – reinforced with a stainless steel coil designed to provide stability and strength through tortuous anatomy
- Tapered tip dilator with hydrophilic coating provides easy access
- Transparent tip with platinum ring helps in better visualization and placement of sheath under fluoroscopy
- Sizes available for all types of flexible scopes (Storz, Wolf & Single use flexible scopes)
- Seamless transition between dilator and sheath for easy and atraumatic insertion

Size (in Fr)

| | | | | |
|----|------|----|------|----|
| ID | 9.5 | 10 | 11.5 | 12 |
| OD | 11.5 | 12 | 13.5 | 14 |

Length (in cm)

28, 35, 45 & 55



PRAVEZA[®] URETERAL ACCESS SHEATH

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUAB095028B1 | 9.5 | 28 |
| BUAB095035B1 | 9.5 | 35 |
| BUAB095045B1 | 9.5 | 45 |
| BUAB095055B1 | 9.5 | 55 |
| BUAB010028B1 | 10 | 28 |
| BUAB010035B1 | 10 | 35 |
| BUAB010045B1 | 10 | 45 |
| BUAB010055B1 | 10 | 55 |
| BUAB115028B1 | 11.5 | 28 |
| BUAB115035B1 | 11.5 | 35 |
| BUAB115045B1 | 11.5 | 45 |
| BUAB115055B1 | 11.5 | 55 |
| BUAB012028B1 | 12 | 28 |
| BUAB012035B1 | 12 | 35 |
| BUAB012045B1 | 12 | 45 |
| BUAB012055B1 | 12 | 55 |

CATHETERS

- Ureteral Catheter
- Mono J Catheter
- PCN Catheter
- PCN Catheter with Locking Mechanism
- PCN Catheter with Needle
(Nephrostomy Drainage Set Direct Puncture)
- Malecot Catheter
- Re-entry Malecot Drainage Catheter
- Suprapubic Catheter (SPC)
- CI Catheter

UROLOGY

URETERAL CATHETER

- Used for catheterization of the urinary tract during intraoperative procedures
- Used for injection of contrast media

Key Features

- Optimal radiopacity for better visualization
- With 1 cm marking which helps in catheter position, measurement of ureter length
- Detachable hub for scope insertion

Types


- Open End
- Closed End
- Cone Tip
- Angled Tip

Size (in Fr)

3, 4, 5, 6, & 7

Length (in cm)

70



URETERAL CATHETER STRAIGHT OPEN END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUCP003070B1 | 3 | 70 |
| BUCP004070B1 | 4 | 70 |
| BUCP005070B1 | 5 | 70 |
| BUCP006070B1 | 6 | 70 |
| BUCP007070B1 | 7 | 70 |

URETERAL CATHETER STRAIGHT CLOSED END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUCP003070B2 | 3 | 70 |
| BUCP004070B2 | 4 | 70 |
| BUCP005070B2 | 5 | 70 |
| BUCP006070B2 | 6 | 70 |
| BUCP007070B2 | 7 | 70 |

URETERAL CATHETER CONE TIP OPEN END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUCP003070B3 | 3 | 70 |
| BUCP004070B3 | 4 | 70 |
| BUCP005070B3 | 5 | 70 |
| BUCP006070B3 | 6 | 70 |
| BUCP007070B3 | 7 | 70 |

URETERAL CATHETER CONE TIP CLOSED END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUCP003070B4 | 3 | 70 |
| BUCP004070B4 | 4 | 70 |
| BUCP005070B4 | 5 | 70 |
| BUCP006070B4 | 6 | 70 |
| BUCP007070B4 | 7 | 70 |

URETERAL CATHETER ANGLED TIP OPEN END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUCP003070B5 | 3 | 70 |
| BUCP004070B5 | 4 | 70 |
| BUCP005070B5 | 5 | 70 |
| BUCP006070B5 | 6 | 70 |
| BUCP007070B5 | 7 | 70 |

URETERAL CATHETER ANGLED TIP CLOSED END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUCP003070B6 | 3 | 70 |
| BUCP004070B6 | 4 | 70 |
| BUCP005070B6 | 5 | 70 |
| BUCP006070B6 | 6 | 70 |
| BUCP007070B6 | 7 | 70 |

MONO J CATHETER

Used for temporary drainage and irrigation of urine from the kidney (Pigtail ureteral catheter).

Key Features

- Pigtail loop helps in retention of distal tip in the renal pelvis
- Holes in shaft helps in better drainage
- Used to inject contrast
- Used for short term external drainage to avoid stenting

Also available in Set, which Includes

Mono J catheter, guidewire & urine bag connector

Type


- Open End
- Closed End

Size (in Fr)

4.7, 5, 6, 7 & 8

Length (in cm)

70



MONO J CATHETER OPEN END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BMJBO47070B1 | 4.7 | 70 |
| BMJBO05070B1 | 5 | 70 |
| BMJBO06070B1 | 6 | 70 |
| BMJBO07070B1 | 7 | 70 |
| BMJBO08070B1 | 8 | 70 |

MONO J CATHETER CLOSED END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BMJBO47070B2 | 4.7 | 70 |
| BMJBO05070B2 | 5 | 70 |
| BMJBO06070B2 | 6 | 70 |
| BMJBO07070B2 | 7 | 70 |
| BMJBO08070B2 | 8 | 70 |

MONO J CATHETER SET OPEN END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BMJBO47070B3 | 4.7 | 70 |
| BMJBO05070B3 | 5 | 70 |
| BMJBO06070B3 | 6 | 70 |
| BMJBO07070B3 | 7 | 70 |
| BMJBO08070B3 | 8 | 70 |

MONO J CATHETER SET CLOSED END

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BMJBO47070B4 | 4.7 | 70 |
| BMJBO05070B4 | 5 | 70 |
| BMJBO06070B4 | 6 | 70 |
| BMJBO07070B4 | 7 | 70 |
| BMJBO08070B4 | 8 | 70 |

PCN CATHETER

PCN catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine.

Key Features

- The pigtail tip configuration helps in the retention of the distal tip in renal pelvis
- Eye holes for better drainage
- Hydrophilic coating for easy insertion

Package Includes

PCN catheter with metal & plastic straightener and urine bag connector

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

30



PCN CATHETER

PCN CATHETER (BALLOON TYPE)

| Product Code | Size (in Fr) | Length (in cm) | Product Code | Size (in Fr) |
|--------------|--------------|----------------|--------------|--------------|
| BPCP006030B1 | 6 | 30 | BPCP010000B8 | 10 |
| BPCP007030B1 | 7 | 30 | BPCP012000B8 | 12 |
| BPCP008030B1 | 8 | 30 | BPCP014000B8 | 14 |
| BPCP010030B1 | 10 | 30 | BPCP016000B8 | 16 |
| BPCP012030B1 | 12 | 30 | BPCP018000B8 | 18 |
| BPCP014030B1 | 14 | 30 | BPCP020000B8 | 20 |
| BPCP016030B1 | 16 | 30 | | |

PCN CATHETER WITH LOCKING MECHANISM

PCN catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine. The locking mechanism prevents accidental slippage of the catheter.

Key Features

- The pigtail tip configuration helps in the retention of the distal tip in renal pelvis
- The locking mechanism prevents accidental slippage of the catheter
- Eye holes for better drainage
- Hydrophilic coating for easy insertion

Package Includes

PCN catheter with locking mechanism, metal & plastic straightener and urine bag connector

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

35



Note: Once the catheter is deployed, pull the string, wind and lock to the hub of the catheter.

PCN CATHETER WITH LOCKING MECHANISM

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BPCP006035B2 | 6 | 35 |
| BPCP007035B2 | 7 | 35 |
| BPCP008035B2 | 8 | 35 |
| BPCP010035B2 | 10 | 35 |
| BPCP012035B2 | 12 | 35 |
| BPCP014035B2 | 14 | 35 |
| BPCP016035B2 | 16 | 35 |

PCN CATHETER WITH NEEDLE

(NEPHROSTOMY DRAINAGE SET DIRECT PUNCTURE)

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage by direct puncture.

Key Features

- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Eye holes for better drainage
- Hydrophilic coating for easy insertion
- Trocar design and catheter design prevents buckling of catheter during the procedure
- Available with or without locking mechanism

Package Includes


PCN catheter with metal straightner, metal trocar, urine bag connector and scalpel blade

PCN CATHETER WITH NEEDLE

| Size (in Fr) | Length (in cm) |
|--------------------------------|----------------|
| 7, 8, 9, 10, 12, 14 16 & 18 | 30 |

PCN CATHETER WITH NEEDLE & LOCKING MECHANISM

| Size (in Fr) | Length (in cm) |
|--------------------------|----------------|
| 7, 8, 10, 12, 14 & 16 | 35 |



PCN CATHETER WITH NEEDLE

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BPCP007030B5 | 7 | 30 |
| BPCP008030B5 | 8 | 30 |
| BPCP009030B5 | 9 | 30 |
| BPCP010030B5 | 10 | 30 |
| BPCP012030B5 | 12 | 30 |
| BPCP014030B5 | 14 | 30 |
| BPCP016030B5 | 16 | 30 |
| BPCP018030B5 | 18 | 30 |

PCN CATHETER WITH NEEDLE & LOCKING MECHANISM

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BPCP007035B6 | 7 | 35 |
| BPCP008035B6 | 8 | 35 |
| BPCP010035B6 | 10 | 35 |
| BPCP012035B6 | 12 | 35 |
| BPCP014035B6 | 14 | 35 |
| BPCP016035B6 | 16 | 35 |

MALECOT CATHETER

Malecot catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine.

Key Features

- The Malecot flower design is ideally suitable for the drainage of thick viscous fluids
- Integral flower design matches with the catheter size when closed and cause less trauma during insertion

Package Includes

Malecot catheter with plastic & metal straightner and urine bag connector

Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

Length (in cm)

30



MALECOT CATHETER

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BMCPO08030B1 | 8 | 30 |
| BMCPO10030B1 | 10 | 30 |
| BMCPO12030B1 | 12 | 30 |
| BMCPO14030B1 | 14 | 30 |
| BMCPO16030B1 | 16 | 30 |
| BMCPO18030B1 | 18 | 30 |
| BMCPO20030B1 | 20 | 30 |
| BMCPO22030B1 | 22 | 30 |

RE-ENTRY MALECOT DRAINAGE CATHETER

Used to be inserted percutaneously into the renal pelvis upto the ureter to perform external drainage of urine from the renal pelvis and internal drainage into the bladder.

Key Features

- The Malecot flower design is ideally suitable for the drainage of thick viscous fluids
- Integral flower design matches with the catheter size when closed and cause less trauma during insertion
- Used as a nephrostomy catheter and also as ureteral stent

Package Includes


Re-Entry Malecot Drainage Catheter with plastic straightener & urine bag connector

Size (in Fr)

10/6, 12/6, 14/6, 16/8 & 18/8

Length (in cm)

45



RE-ENTRY MALECOT DRAINAGE CATHETER

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BRMN610045B1 | 10/6 | 45 |
| BRMN612045B1 | 12/6 | 45 |
| BRMN614045B1 | 14/6 | 45 |
| BRMN816045B1 | 16/8 | 45 |
| BRMN818045B1 | 18/8 | 45 |

SUPRAPUBIC CATHETER (SPC)

Suprapubic catheter is intended to be inserted percutaneously into the urinary bladder and perform external drainage of urine from the urinary bladder.

Key Features

- Integral flower design matches with the catheter size when closed and cause less trauma during insertion
- Trocar and catheter designed to prevent buckling of catheter during the procedure

Package Includes

Suprapubic catheter with metal straightner, metal trocar and urine bag connector

Size (in Fr)

8, 10, 12, 14, 16 & 18

Length (in cm)

25



SUPRAPUBIC CATHETER

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BSCP008025B1 | 8 | 25 |
| BSCP010025B1 | 10 | 25 |
| BSCP012025B1 | 12 | 25 |
| BSCP014025B1 | 14 | 25 |
| BSCP016025B1 | 16 | 25 |
| BSCP018025B1 | 18 | 25 |

CI CATHETER

Used to drain urine from bladder.

Key Features

- Easy advancement minimizing the potential trauma
- Firm enough for easy insertion, yet gentle on soft urethral tissue
- Eye holes for better drainage

CI Catheter Female

| Size (in Fr) | Length (in cm) |
|---------------------|----------------|
| 10, 12, 14, 16 & 18 | 15 |

CI Catheter Male

| Size (in Fr) | Length (in cm) |
|---------------------|----------------|
| 10, 12, 14, 16 & 18 | 35 |



CI CATHETER FEMALE

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BCI0010015B0 | 10 | 15 |
| BCI0012015B0 | 12 | 15 |
| BCI0014015B0 | 14 | 15 |
| BCI0016015B0 | 16 | 15 |
| BCI0018015B0 | 18 | 15 |

CI CATHETER MALE

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BCI0010035B0 | 10 | 35 |
| BCI0012035B0 | 12 | 35 |
| BCI0014035B0 | 14 | 35 |
| BCI0016035B0 | 16 | 35 |
| BCI0018035B0 | 18 | 35 |

DRAINAGE SETS

- PCN Set - Percutaneous Nephrostomy Drainage Set "Pigtail"
- PCN Set - Pigtail with Locking Mechanism
- PCN Set - Percutaneous Nephrostomy Drainage Set "Malecot"
- Suprapubic Drainage Catheter Set

UROLOGY

PCN SET (PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL")

PCN Set is used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage.

Advantages

Kit provides entire devices for performing Nephrostomy procedure

Set Contains

IP Needle | Guidewire J Tip PTFE Coated
| Fascial Dilators | Urine Bag Connector
PCN Catheter with Hydrophilic Coating
Scalpel Blade | Catheter Fixing Disc

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

30



PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL"

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BPCP006030B3 | 6 | 30 |
| BPCP007030B3 | 7 | 30 |
| BPCP008030B3 | 8 | 30 |
| BPCP010030B3 | 10 | 30 |
| BPCP012030B3 | 12 | 30 |
| BPCP014030B3 | 14 | 30 |
| BPCP016030B3 | 16 | 30 |

PCN SET - PIGTAIL WITH LOCKING MECHANISM

(PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL")

PCN Set is used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage.

Advantages

Kit provides entire devices for performing Nephrostomy procedure

Set Contains

IP Needle | Guidewire J Tip PTFE Coated | Fascial Dilators | Urine Bag Connector
PCN Catheter with Locking Mechanism & Hydrophilic Coating | Scalpel Blade
Catheter Fixing Disc

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

35



PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL" WITH LOCKING MECHANISM

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BPCP006035B4 | 6 | 35 |
| BPCP007035B4 | 7 | 35 |
| BPCP008035B4 | 8 | 35 |
| BPCP010035B4 | 10 | 35 |
| BPCP012035B4 | 12 | 35 |
| BPCP014035B4 | 14 | 35 |
| BPCP016035B4 | 16 | 35 |

PCN SET (PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "MALECOT")

PCN Set is used for percutaneous placement of a Malecot catheter in the renal pelvis for nephrostomy drainage.

Advantages

Kit provides entire devices for performing Nephrostomy procedure

Set Contains

IP Needle | Guidewire J Tip PTFE Coated | Fascial Dilators | Urine Bag Connector
Malecot Catheter with Hydrophilic Coating | Scalpel Blade | Catheter Fixing Disc

Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

Length (in cm)

30



PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "MALECOT"

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BMCP008030B2 | 8 | 30 |
| BMCP010030B2 | 10 | 30 |
| BMCP012030B2 | 12 | 30 |
| BMCP014030B2 | 14 | 30 |
| BMCP016030B2 | 16 | 30 |
| BMCP018030B2 | 18 | 30 |
| BMCP020030B2 | 20 | 30 |
| BMCP022030B2 | 22 | 30 |

SUPRAPUBIC DRAINAGE CATHETER SET (SDS)

Suprapubic drainage catheter set is intended to create suprapubic access to the bladder and place a catheter for bladder drainage.

Advantages

- Silicone foley catheter for patient comfort and safety
- Specially designed trocar with needle tip enhances smooth puncture
- Peel away sheath forms a conduit for catheter placement and can be peeled easily to remove from the catheter
- Seamless transition between trocar and peel away sheath for easy and atraumatic puncture

Set Contains

Foley catheter, trocar with peel-away sheath and scalpel blade

Size (in Fr)

10, 12, 14, 16 & 18



SUPRAPUBIC DRAINAGE CATHETER SET - SILICONE

| Product Code | Size (in Fr) |
|--------------|--------------|
| BSD0010025B1 | 10 |
| BSD0012025B1 | 12 |
| BSD0014025B1 | 14 |
| BSD0016025B1 | 16 |
| BSD0018025B1 | 18 |

STONE BASKETS

- N Quadra™ - Nitinol Stone Basket With Tip
- N Quadra™ Z - Nitinol Stone Basket Zero Tip
- N Quadra™ Helical - Nitinol Stone Basket Helical With Tip
- N Quadra™ Helical Z - Nitinol Stone Basket Helical Zero Tip
- N Quadra™ Flexible - Nitinol Stone Basket Flexible
- Delta Basket
- PCNL Basket

UROLOGY

N QUADRA™ - NITINOL STONE BASKET WITH TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly
- Atraumatic tip

Size (in Fr)

3 & 4

Length (in cm)

70, 90 & 120



N QUADRA™ - NITINOL STONE BASKET WITH TIP

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BSBN003070BT1 | 3 | 70 |
| BSBN003090BT1 | 3 | 90 |
| BSBN003120BT1 | 3 | 120 |
| BSBN004070BT1 | 4 | 70 |
| BSBN004090BT1 | 4 | 90 |
| BSBN004120BT1 | 4 | 120 |

N QUADRA™ Z - NITINOL STONE BASKET ZERO TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Zero tip design reduces chances of perforation and mucosal damage, which also avoids stone migration during the procedure
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

3 & 4

Length (in cm)

70, 90 & 120



N QUADRA™ Z - NITINOL STONE BASKET ZERO TIP

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BSBN003070BTO | 3 | 70 |
| BSBN003090BTO | 3 | 90 |
| BSBN003120BTO | 3 | 120 |
| BSBN004070BTO | 4 | 70 |
| BSBN004090BTO | 4 | 90 |
| BSBN004120BTO | 4 | 120 |

N QUADRA™ HELICAL - NITINOL STONE BASKET HELICAL WITH TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- Helical design helps in capturing small stones and fragments effectively
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

3 & 4

Length (in cm)

70, 90 & 120



N QUADRA™ HELICAL - NITINOL STONE BASKET HELICAL WITH TIP

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BSBN003070BHI | 3 | 70 |
| BSBN003090BHI | 3 | 90 |
| BSBN003120BHI | 3 | 120 |
| BSBN004070BHI | 4 | 70 |
| BSBN004090BHI | 4 | 90 |
| BSBN004120BHI | 4 | 120 |

N QUADRA™ HELICAL Z - NITINOL STONE BASKET HELICAL ZERO TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Zero tip design reduces chances of perforation and mucosal damage, which also avoids stone migration during the procedure
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- Helical design helps in capturing small stones and fragments effectively
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

3 & 4

Length (in cm)

70, 90 & 120



N QUADRA™ HELICAL Z - NITINOL STONE BASKET HELICAL ZERO TIP

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BSBN003070BHO | 3 | 70 |
| BSBN003090BHO | 3 | 90 |
| BSBN003120BHO | 3 | 120 |
| BSBN004070BHO | 4 | 70 |
| BSBN004090BHO | 4 | 90 |
| BSBN004120BHO | 4 | 120 |

N QUADRA™ FLEXIBLE - NITINOL STONE BASKET FLEXIBLE

Flexible Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from calyx.

Key Features

- Flexible design shaft which can angulate up to 270 degrees as per the flexible scope direction during RIRS Procedure
- Kink resistant braided shaft provides durability, making removal of several fragments possible
- Tipless design helps in easy capturing of stones by allowing the basket to open at the base of calyx and encircle the stone
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

1.8, 1.9, 2, 2.2 & 2.4

Length (in cm)

120



N QUADRA™ FLEXIBLE - NITINOL STONE BASKET FLEXIBLE

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BFBB018120BTO | 1.8 | 120 |
| BFBB019120BTO | 1.9 | 120 |
| BFBB002120BTO | 2 | 120 |
| BFBB022120BTO | 2.2 | 120 |
| BFBB024120BTO | 2.4 | 120 |

DELTA BASKET

Steerable Basket for rigid nephroscope for stone extraction from calices which are inaccessible due to acute angle.

Design Features

- Can be angled inside kidney to reach inner calyx
- Can be used for straight path stone manipulation and also can be angled both directions
- Unique basket shape which prevents stone migration

Procedural Features

- Prevents injury to the infundibulum and reduce bleeding
- Avoids the need for multiple tracts
- Can be used without the assistance of flexible nephroscope
- Helps to remove big blood clots

Size (in Fr)

7

Length (in cm)

40



DELTA BASKET

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BDBR007040B0 | 7 | 40 |

PCNL BASKET

PCNL Basket is intended to be used with rigid nephroscope for stone extraction from calices.

Design Features

- Easy to handle - single handed operation
- Can be used for straight path stone manipulation
- Unique basket shape which prevents stone migration

Size (in Fr)

7 & 9

Length (in cm)

40



PCNL BASKET

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BPBR007040B0 | 7 | 40 |
| BPBR009040B0 | 9 | 40 |

DILATORS

- Fascial Dilator Set
- Ureteral Dilator Set
- Single Step™ Ureteral Dilator
- Urethral Dilator Set
- Urethral Dilator Set - S Curved
- Urethral Self Dilator
- Amplatz Sheath
- Amplatz Renal Dilator
- Amplatz Renal Dilation Set - S
- Amplatz Renal Dilation Set - F
- Mini Perc Renal Dilation Set
- Ureteral Balloon Dilator
- Nephrostomy Balloon Dilator With Amplatz Sheath
- Screw Dilator

UROLOGY

FASCIAL DILATOR SET

Used for fascial dilation, intended for stretching or enlarging the fascial tissue covering a cavity, tract or opening prior to a diagnostic or interventional procedure.

Key Features


- Gradual tapered end & smooth surface for atraumatic entry
- Also available with hydrophilic coating for ease of introduction
- Optimal radiopacity for better visualization

Size (in Fr)

6, 8, 10, 12, 14, 16, 18 & 20

Length (in cm)

20 & 22



FASCIAL DILATOR - SINGLE

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BFDN006020BS | 6 | 20 |
| BFDN008020BS | 8 | 20 |
| BFDN010020BS | 10 | 20 |
| BFDN012020BS | 12 | 20 |
| BFDN014020BS | 14 | 20 |
| BFDN016020BS | 16 | 20 |
| BFDN018020BS | 18 | 20 |
| BFDN020020BS | 20 | 20 |
| BFDN006022BS | 6 | 22 |
| BFDN008022BS | 8 | 22 |
| BFDN010022BS | 10 | 22 |
| BFDN012022BS | 12 | 22 |
| BFDN014022BS | 14 | 22 |
| BFDN016022BS | 16 | 22 |
| BFDN018022BS | 18 | 22 |
| BFDN020022BS | 20 | 22 |

FASCIAL DILATOR SET

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BFDN616020BSE | 6-16 | 20 |
| BFDN620020BSE | 6-20 | 20 |
| BFDN616022BSE | 6-16 | 22 |
| BFDN620022BSE | 6-20 | 22 |

URETERAL DILATOR SET

Used for dilation of the ureter through natural orifice prior to a diagnostic or interventional procedure.

Key Features


- Tapered end & smooth surface for atraumatic entry during dilation
- Also available with hydrophilic coating for ease of introduction
- Optimal radiopacity for better visualization

Size (in Fr)

6, 8, 10, 12, 14 & 16

Length (in cm)

70



URETERAL DILATOR - SINGLE

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BUDN006070BS | 6 | 70 |
| BUDN008070BS | 8 | 70 |
| BUDN010070BS | 10 | 70 |
| BUDN012070BS | 12 | 70 |
| BUDN014070BS | 14 | 70 |
| BUDN016070BS | 16 | 70 |

URETERAL DILATOR SET

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BUDN616070BSE | 6-16 | 70 |

SINGLE STEP[™] URETERAL DILATOR

Used for one step dilation of the ureter through natural orifice prior to a diagnostic or interventional procedure.

Key Features

- Long tapered tip avoids use of sequential dilation
- Smooth surface for atraumatic entry during dilation
- Also available with hydrophilic coating for ease of introduction
- Optimal radiopacity for better visualization

Size (in Fr)

10/6, 12/6, 14/6, 16/6

Length (in cm)

70



SINGLE STEP[™] URETERAL DILATOR

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BNDP610070BO | 10/6 | 70 |
| BNDP612070BO | 12/6 | 70 |
| BNDP614070BO | 14/6 | 70 |
| BNDP616070BO | 16/6 | 70 |

URETHRAL DILATOR SET

Used for dilation of the Urethra.

Key Features

- Tapered end for atraumatic entry during dilation
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization
- To be used for sequential dilation
- Side hole for drainage

Type


Straight

Size (in Fr)

8-24 (8, 10, 12, 14, 16, 18, 20, 22 & 24)
8-30 (8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28 & 30)

Length (in cm)

40



URETHRAL DILATOR SET

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BUD0824040BSE | 8-24 | 40 |
| BUD0830040BSE | 8-30 | 40 |

URETHRAL DILATOR SET - S CURVED

Used for dilation of the Male Urethral strictures.

Key Features

- Tapered end for atraumatic entry during dilation
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization
- To be used for sequential dilation
- Side hole for drainage
- S curve supports natural curvature of male urethra

Type


S Curved

Size (in Fr)

8-24 (8, 10, 12, 14, 16,18, 20, 22 & 24)
8-30 (8, 10, 12, 14, 16,18, 20, 22, 24, 26, 28 & 30)

Length (in cm)

40



URETHRAL DILATOR SET - S CURVED

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BUD0824040BSS | 8-24 | 40 |
| BUD0830040BSS | 8-30 | 40 |

URETHRAL SELF DILATOR

Used for self-dilation of Urethral strictures.

Key Features

- Tapered end for atraumatic entry during dilation
- Smooth surface for ease of introduction

Types

- Straight With Hole
- Curved With Hole

Size (in Fr)

6, 8, 10, 12, 14 & 16

Length (in cm)

20 & 30



URETHRAL SELF DILATOR - STRAIGHT WITH HOLE

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BSD0610030B1 | 10 | 30 |
| BSD0612030B1 | 12 | 30 |
| BSD0614030B1 | 14 | 30 |
| BSD0616030B1 | 16 | 30 |

URETHRAL SELF DILATOR - CURVED WITH HOLE

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BSD0610030B3 | 10 | 30 |
| BSD0612030B3 | 12 | 30 |
| BSD0614030B3 | 14 | 30 |
| BSD0616030B3 | 16 | 30 |

AMPLATZ SHEATH

Intended for percutaneous dilation of tissue surrounding a dilator previously placed. It is used as a working channel to maintain a previously established nephrostomy tract for introduction of instruments and catheters during urological procedures.

Key Features

- Optimal bevel tip for atraumatic entry
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization

Size (in Fr)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

Length (in cm)

17 & 20



AMPLATZ SHEATH

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BAST012017B1 | 12 | 17 |
| BAST014017B1 | 14 | 17 |
| BAST016017B1 | 16 | 17 |
| BAST018017B1 | 18 | 17 |
| BAST020017B1 | 20 | 17 |
| BAST022017B1 | 22 | 17 |
| BAST024017B1 | 24 | 17 |
| BAST026017B1 | 26 | 17 |
| BAST028017B1 | 28 | 17 |
| BAST030017B1 | 30 | 17 |
| BAST012020B1 | 12 | 20 |
| BAST014020B1 | 14 | 20 |
| BAST016020B1 | 16 | 20 |
| BAST018020B1 | 18 | 20 |
| BAST020020B1 | 20 | 20 |
| BAST022020B1 | 22 | 20 |
| BAST024020B1 | 24 | 20 |
| BAST026020B1 | 26 | 20 |
| BAST028020B1 | 28 | 20 |
| BAST030020B1 | 30 | 20 |

AMPLATZ RENAL DILATOR

(RENAL DILATOR WITH AMPLATZ SHEATH)

Used for dilation of a percutaneous tract and setup of a working channel into the kidney for PCNL procedures by placing a suitable amplatz sheath.

Key Features

- Gradual taper designed for atraumatic placement and maneuverability
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization

Dilator Size (OD of Dilator in Fr)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

Length (in cm)

30

Amplatz Sheath Size (ID of Sheath in Fr)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

Length (in cm)

17



AMPLATZ RENAL DILATOR

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BRD0012030BS | 12 | 30 |
| BRD0014030BS | 14 | 30 |
| BRD0016030BS | 16 | 30 |
| BRD0018030BS | 18 | 30 |
| BRD0020030BS | 20 | 30 |
| BRD0022030BS | 22 | 30 |
| BRD0024030BS | 24 | 30 |
| BRD0026030BS | 26 | 30 |
| BRD0028030BS | 28 | 30 |
| BRD0030030BS | 30 | 30 |

AMPLATZ RENAL DILATION SET - S

Intended for stretching or enlarging the fascial tissue prior to a nephrostomy procedure and to introduce the amplatz sheath which is intended for percutaneous dilation of tissue surrounding a dilator previously placed and used as a working channel during urological procedures.


Set Includes

- 8 Fr guiding catheter
- 10 Fr radiopaque dilation sheath
- Radiopaque fascial dilators of sizes from 6 Fr – 10 Fr of length 20 cm
- Radiopaque renal dilators of sizes from 12 Fr – 30 Fr of length 30 cm
- Teflon amplatz sheath for renal dilators of sizes 24, 26, 28 & 30 Fr of length 17 cm
- Packed in tray

AMPLATZ RENAL DILATION SET - S

Product Code

BAR0630000BSE



AMPLATZ RENAL DILATION SET - F

Intended for stretching or enlarging the fascial tissue prior to a nephrostomy procedure and to introduce the amplatz sheath which is intended for percutaneous dilation of tissue surrounding a dilator previously placed and used as a working channel during urological procedures.


Set Includes

- Scalpel blade
- 2 part 18 G x 15 cm IP needle with echogenic tip
- 8 Fr guiding catheter
- 0.035" PTFE extra stiff J tip guidewire
- 10 Fr radiopaque dilation sheath
- Radiopaque fascial dilators of sizes from 6 Fr – 10 Fr of length 20 cm
- Radiopaque renal dilators of sizes from 12 Fr – 30 Fr of length 30 cm
- Teflon amplatz sheath for renal dilators of sizes 24, 26, 28 & 30 Fr of length 17 cm
- Packed in tray

AMPLATZ RENAL DILATION SET - F

Product Code

BAR0630000B53



MINI PERC RENAL DILATION SET

Intended for stretching or enlarging the fascial tissue prior to a nephrostomy procedure and to introduce the amplatz sheath which is intended for percutaneous dilation of tissue surrounding a dilator previously placed and used as a working channel during urological procedures.


Set Includes

- Scalpel blade
- Trocar Needle with sheath 6 Fr x 15 cm
- 8 Fr guiding catheter with compatible 0.035" PTFE extra stiff J tip guidewire
- 10 Fr radiopaque dilation sheath
- Radiopaque fascial dilators of sizes from 6 Fr - 10 Fr of length 20 cm
- Radiopaque renal dilators of sizes from 12 Fr - 22 Fr of length 30 cm
- Teflon amplatz sheath for renal dilators of sizes 14, 16, 18, 20 & 22 Fr of length 17 cm
- Packed in tray

MINI PERC RENAL DILATION SET

Product Code

BAM0622000BSE



URETERAL BALLOON DILATOR


Used for dilation of narrow ureteral segments.

Product Features

- Controlled inflation pressure allows dilating force at a stricture or ureteral narrowing
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy

Product Specifications

- Shaft Size: 5 Fr
- Catheter Length: 75 cm
- Balloon Length: 4 cm
- Balloon Expansion Diameter: 5 mm
- Nominal Pressure: 7 atm
- Rate Burst Pressure: 14 atm
- Use Guidewire: 0.035 inch



URETERAL BALLOON DILATOR

URETERAL BALLOON DILATOR WITH INFLATION DEVICE

| Product Code | Size (in Fr) | Length (in cm) | Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|--------------|--------------|----------------|
| BUB0005075B0 | 5 | 75 | BUB0005075B1 | 5 | 75 |

NEPHROSTOMY BALLOON DILATOR WITH AMPLATZ SHEATH


Used for dilation of percutaneous tract for nephrostomy procedure.

Product Features

- Controlled inflation pressure allows dilating force at a nephrostomy tract
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy

Product Specifications

- Shaft Size: 5 Fr
- Catheter Length: 55 cm
- Balloon Length: 12 cm
- Balloon Expansion Diameter: 10 mm
- Nominal Pressure: 2 atm
- Rate Burst Pressure: 11 atm
- Use Guidewire: 0.035 inch



NEPHROSTOMY BALLOON DILATOR

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BNB0005055B0 | 5 | 55 |

NEPHROSTOMY BALLOON DILATOR WITH AMPLATZ SHEATH

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BNB0005055B1 | 5 | 55 |

NEPHROSTOMY BALLOON DILATOR WITH AMPLATZ SHEATH & INFLATION DEVICE

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BNB0005055B2 | 5 | 55 |

SCREW DILATOR

Used for fascial dilation in one step.

Key Features

- Radiopaque for better visualization
- Screw tip for one step dilation of fascia during mini PCNL procedure
- Also available with compatible amplatz sheath
- Avoids sequential dilatation, reducing procedural time

Size (in Fr)

12, 14 & 16

Length (in cm)

25



SCREW DILATOR

| Product Code | Size (in Fr) | Length (in cm) |
|--------------|--------------|----------------|
| BOS0012025BO | 12 | 25 |
| BOS0014025BO | 14 | 25 |
| BOS0016025BO | 16 | 25 |

SCREW DILATOR WITH AMPLATZ SHEATH

| Product Code | Size (in Fr) | Length (in cm) |
|---------------|--------------|----------------|
| BOS0012025BAS | 12 | 25 |
| BOS0014025BAS | 14 | 25 |
| BOS0016025BAS | 16 | 25 |

NEEDLES

- Initial Puncture Needle - Two Part
- Sure Puncture™ Needle (Initial Puncture Needle - Three Part)

UROLOGY

INITIAL PUNCTURE NEEDLE - TWO PART

Intended to make the initial percutaneous puncture to access the kidney and to enable the introduction of a guidewire during urological procedures (PCN, PCNL, Mini PERC & Micro PERC).

Key Features

- Siliconised needle for less trauma
- Echogenic tip for better visualization during ultrasound
- 1 cm marking on the shaft

Advantage of Locking Hub IP Needle

During the procedure the locking hub helps in holding the inner needle with outer cannula together, which helps in easy initial puncture

Type

Two Part

Size (in G)

18

Length (in cm)

15 & 20



INITIAL PUNCTURE NEEDLE - TWO PART WITH ECHOGENIC TIP

| Product Code | Size (in G) | Length (in cm) |
|---------------|-------------|----------------|
| BIP50180152PB | 18 | 15 |
| BIP50180202PB | 18 | 20 |

INITIAL PUNCTURE NEEDLE - TWO PART WITH ECHOGENIC TIP & LOCKING HUB

| Product Code | Size (in G) | Length (in cm) |
|---------------|-------------|----------------|
| BIP5018015LHB | 18 | 15 |
| BIP5018020LHB | 18 | 20 |

SURE PUNCTURE™ NEEDLE

(INITIAL PUNCTURE NEEDLE - THREE PART)

Intended to make the initial percutaneous puncture to access the kidney and to enable the introduction of a guidewire during urological procedures (PCN, PCNL, Mini PERC & Micro PERC).

Key Features

- Siliconised needle for less trauma
- Echogenic tip for better visualization during ultrasound
- 1 cm marking

Advantage of Sure Puncture™ Needle

- Combination of 21 G Bevel tip and 16 G Cannula Needle
- Inner 21 G provides minimal site puncture and contrast injection to ensure desired puncture location
- Coaxial needle provides secured access
- Can introduce working guidewire with 21 G puncture

Type

Three Part

Size (in G)

21/16

Length (in cm)

20



SURE PUNCTURE™ NEEDLE WITH ECHOGENIC TIP

| Product Code | Size (in G) | Length (in cm) |
|---------------|-------------|----------------|
| BSPSO160203PB | 21/16 | 20 |

BIOPSY

- Biopsy Gun - Semi Automatic

UROLOGY

BIOPSY GUN - SEMI AUTOMATIC

Used for obtaining histological core samples from soft tissue and organs.

Key Features

- Echogenic and atraumatic needle tip
- Ergonomic and light weight
- Cannula with depth marking
- Easy and one hand operated with two stage arming
- Coaxial available with trocar and blunt needle

Biopsy Gun - Semi Automatic With Coaxial & Blunt Needle

- 14, 16 G x 10, 15 cm - For breast biopsy
- 16, 18 G x 15, 20 cm - For liver & kidney
- 18, 20 G x 15, 20 cm - For lung



BIOPSY GUN - SEMI AUTOMATIC WITH COAXIAL & BLUNT NEEDLE

| Product Code | Size (in G) | Length (in cm) |
|---------------|-------------|----------------|
| BBG0014010BS2 | 14 | 10 |
| BBG0014015BS2 | 14 | 15 |
| BBG0016010BS2 | 16 | 10 |
| BBG0016015BS2 | 16 | 15 |
| BBG0016020BS2 | 16 | 20 |
| BBG0018015BS2 | 18 | 15 |
| BBG0018020BS2 | 18 | 20 |
| BBG0020015BS2 | 20 | 15 |
| BBG0020020BS2 | 20 | 20 |

ACCESSORIES

- Meatal Dilator - Male
- Dalela Female Self Urethral Dilator
- Bladder Evacuator
- Penile Clamp
- Irrigation Set (TUR Set)
- Cath Tap

UROLOGY

MEATAL DILATOR - MALE

Used for self-dilation of penile urethral meatus.

Key Features

- Designed specifically for dilation of meatal strictures
- Atraumatic design
- Uniform taper and smooth surface for easy insertion
- Superior handling grip

Size (in Fr)

6/30

MEATAL DILATOR - MALE

| Product Code | Size (in Fr) |
|--------------|--------------|
| BMD0030000BM | 6/30 |



DALELA FEMALE SELF URETHRAL DILATOR

PATENT PENDING

Used for self dilation of female urethra.

Key Features

- Novel design specifically for dilation of female urethra
- Made of durable, biocompatible material
- Soft tip for atraumatic and easy insertion
- Less trauma and pain during dilation
- Easy to use & extremely user friendly
- Uniform taper and smooth surface
- Can empty the bladder while device in place

Size (in Fr)

8/30

DALELA FEMALE SELF URETHRAL DILATOR

| Product Code | Size (in Fr) |
|--------------|--------------|
| BMD0030000BD | 8/30 |



BLADDER EVACUATOR

Used for irrigation and tissue collection during transurethral bladder and prostate surgery.

Key Features

- Used for TURP procedures
- Easy to fit connection to resectoscope instrument
- Soft easy to squeeze bottle design, allows fast and efficient removal of tissue during TURP

BLADDER EVACUATOR

Product Code

BBE000000B1



PENILE CLAMP

Used for temporary control of male urinary incontinence.

PENILE CLAMP

Product Code

BPC000000B0



IRRIGATION SET (TUR SET)

Used for endoscopic irrigation during transurethral resection of prostate gland.

Key Features

Non kinkable PVC tube

IRRIGATION SET - Y

Product Code

BTS000000B0



CATH TAP

Used for intermittent drainage of foley catheter.

Key Features

- Intermittent drainage of foley catheter
- Reduce infection rate
- Avoid carrying of urine bags
- Can be fixed to patient thigh
- Useful for patients with long term catheterization

CATH TAP

Product Code

BCT000000B0





BLUENEEM[®]



CE 0123

www.blueneem.com

Manufactured By:

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5, Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112, Karnataka, India

BNMD/MKG/CAT/10 | Rev.01 | 27-Sep-2021

Marketed & Distributed By:

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal Code: 61348
Bad Homburg, Germany

T: +49 6172 17102-20
E: bioactiva@bioactiva.de