

RE-ENTRY MALECOT DRAINAGE CATHETER

Used to gain access to the upper urinary tract & the ureter to perform external drainage of urine from the renal pelvis and internal drainage into the bladder.



Key Features

- Used as a nephrostomy catheter and also as ureteral stent
- The Malecot flower design is ideally suitable for the drainage of thick viscous fluids
- Integral flower is designed to match the catheter size. This reduces trauma during insertion

Package Includes

Re-entry Malecot Catheter with Plastic Straightener & Urine Bag Connector

Size (in Fr)	Length (in cm)
10/6, 12/6, 14/6, 16/8 & 18/8	45

RE-ENTRY MALECOT DRAINAGE CATHETER

Product Code	Size (in Fr)	Length (in cm)
BRMN610045B1	10/6	45
BRMN612045B1	12/6	45
BRMN614045B1	14/6	45
BRMN816045B1	16/8	45
BRMN818045B1	18/8	45

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/10 | Rev.02 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de