

SUPRAPUBIC CATHETER (SPC)

Is intended to be inserted percutaneously into the urinary bladder and perform external drainage of urine from the urinary bladder.



Key Features

- Integral flower is designed to match the catheter size. This reduces trauma during insertion
- Trocar design and catheter design prevents buckling of catheter during the procedure

Package Includes

Catheter with Metal Straightner, Trocar and Urine Bag Connector

Size (in Fr)

8, 10, 12, 14, 16 & 18

Length (in cm)

25

SUPRAPUBIC CATHETER (SPC)

Product Code	Size (in Fr)	Length (in cm)
BSCP008025B1	8	25
BSCP010025B1	10	25
BSCP012025B1	12	25
BSCP014025B1	14	25
BSCP016025B1	16	25
BSCP018025B1	18	25

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/12 | Rev.02 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de