

# URETERAL ACCESS SHEATH

Used to provide ureteral dilatation and working channel for the introduction of Ureteroscopes and devices during ureteroscopy procedures (RIRS).



## Key Features

- Kink resistant – reinforced with a stainless steel coil designed to provide stability and strength through tortuous anatomy
- Transparent Tip with Platinum Ring helps in better visualization and placement of sheath under fluoroscopy
- Tapered Tip Dilator with Hydrophilic coating provides easy access
- Sizes available for all types of flexible scopes (Storz, Wolf & Single use flexible scopes)
- Seamless transition between dilator and sheath for easy and atraumatic insertion

Size (in Fr)					Length (in cm)
ID	9.5	10	11	12	28, 35, 45 & 55
OD	11.5	12	13	14	

## URETERAL ACCESS SHEATH

Product Code	Size (in Fr)	Length (in cm)
BUAB095028B1	9.5	28
BUAB095035B1	9.5	35
BUAB095045B1	9.5	45
BUAB095055B1	9.5	55
BUAB010028B1	10	28
BUAB010035B1	10	35
BUAB010045B1	10	45
BUAB010055B1	10	55
BUAB011028B1	11	28
BUAB011035B1	11	35
BUAB011045B1	11	45
BUAB011055B1	11	55
BUAB012028B1	12	28
BUAB012035B1	12	35
BUAB012045B1	12	45
BUAB012055B1	12	55

[www.bluneem.com](http://www.bluneem.com)

### Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,  
Harohalli Industrial Area - Phase II,  
Kanakapura Taluk, Ramanagara District -562112  
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855  
E: [contact@bluneem.com](mailto:contact@bluneem.com), [marketing@bluneem.com](mailto:marketing@bluneem.com)

BNMD/MKG/LFT/13 | Rev.02 | 12-Dec-2020

Marketed & Distributed By

### Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348  
Bad Homburg, Germany

Tel.: +49 6172 17102-20  
Email: [bioactiva@bioactiva.de](mailto:bioactiva@bioactiva.de)