

INITIAL PUNCTURE NEEDLE

Used for gaining access during Percutaneous Nephrostomy procedure (PCN, PCNL, Mini PERC & Micro PERC)



Key Features

- Siliconised Needle for less trauma
- Echogenic Tip for better visualization during ultrasound
- 1 cm marking

Advantage of Locking Hub IP Needle

The locking hub firmly holds the inner needle with outer cannula. This helps in easy initial puncture.

Types	Size (in G)	Length (in cm)
• Two Part	18	15, 20
• Three Part	18	20

INITIAL PUNCTURE NEEDLE - TWO PART

Product Code	Size (in G)	Length (in cm)
BIPSO180152PB	18	15
BIPSO180202PB	18	20

INITIAL PUNCTURE NEEDLE - THREE PART

Product Code	Size (in G)	Length (in cm)
BIPSO180203PB	18	20

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/15 | Rev.02 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de