


# P FLEX URETERAL STENT SET WITH STRING

Intended to be implanted for treatment of obstructed ureter for drainage of urine from renal pelvis to urinary bladder.



## Key Features

- Biocompatibility tested for indwelling up to 30 Days.
- Polyurethane based stent has uniform diameter and smooth surface reducing encrustation.
- Radiopaque for better visualization
- Stent with String which helps in repositioning

## Set Contains

Stent with String, Pusher, Guide Wire & Clamp

## Types

- Both Ends Open
- One End Closed
- Both Ends Multiloop
- With or Without Hydrophilic coating

## Standard

**Size (in Fr):** 4.7, 4.8, 5, 6, 7 & 8

**Length (in cm):** 12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

## Multiloop

**Size (in Fr):** 4.7, 4.8, 5, 6, 7 & 8

**Length (in cm):** 18-24 & 24-30

## Pediatric

**Size (in Fr):** 3 & 4

**Length (in cm):** 10, 12, 14 & 16

**P FLEX URETERAL STENT SET BOTH ENDS OPEN WITH STRING**

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B9	3	10
BUSB003012B9	3	12
BUSB003014B9	3	14
BUSB003016B9	3	16
BUSB004010B9	4	10
BUSB004012B9	4	12
BUSB004014B9	4	14
BUSB004016B9	4	16
BUSB047012B9	4.7	12
BUSB047014B9	4.7	14
BUSB047016B9	4.7	16
BUSB047018B9	4.7	18
BUSB047020B9	4.7	20
BUSB047022B9	4.7	22
BUSB047024B9	4.7	24
BUSB047026B9	4.7	26
BUSB047028B9	4.7	28
BUSB047030B9	4.7	30
BUSB048012B9	4.8	12
BUSB048014B9	4.8	14
BUSB048016B9	4.8	16
BUSB048018B9	4.8	18
BUSB048020B9	4.8	20
BUSB048022B9	4.8	22
BUSB048024B9	4.8	24
BUSB048026B9	4.8	26
BUSB048028B9	4.8	28
BUSB048030B9	4.8	30

## P FLEX URETERAL STENT SET BOTH ENDS OPEN WITH STRING

Product Code	Size (in Fr)	Length (in cm)
BUSB005012B9	5	12
BUSB005014B9	5	14
BUSB005016B9	5	16
BUSB005018B9	5	18
BUSB005020B9	5	20
BUSB005022B9	5	22
BUSB005024B9	5	24
BUSB005026B9	5	26
BUSB005028B9	5	28
BUSB005030B9	5	30
BUSB006012B9	6	12
BUSB006014B9	6	14
BUSB006016B9	6	16
BUSB006018B9	6	18
BUSB006020B9	6	20
BUSB006022B9	6	22
BUSB006024B9	6	24
BUSB006026B9	6	26
BUSB006028B9	6	28
BUSB006030B9	6	30
BUSB007012B9	7	12
BUSB007014B9	7	14
BUSB007016B9	7	16
BUSB007018B9	7	18
BUSB007020B9	7	20
BUSB007022B9	7	22
BUSB007024B9	7	24
BUSB007026B9	7	26
BUSB007028B9	7	28
BUSB007030B9	7	30

## P FLEX URETERAL STENT SET BOTH ENDS OPEN WITH STRING

Product Code	Size (in Fr)	Length (in cm)
BUSB008012B9	8	12
BUSB008014B9	8	14
BUSB008016B9	8	16
BUSB008018B9	8	18
BUSB008020B9	8	20
BUSB008022B9	8	22
BUSB008024B9	8	24
BUSB008026B9	8	26
BUSB008028B9	8	28
BUSB008030B9	8	30

**P FLEX URETERAL STENT SET ONE END CLOSED WITH STRING**

Product Code	Size (in Fr)	Length (in cm)
BUSBO03010B10	3	10
BUSBO03012B10	3	12
BUSBO03014B10	3	14
BUSBO03016B10	3	16
BUSBO04010B10	4	10
BUSBO04012B10	4	12
BUSBO04014B10	4	14
BUSBO04016B10	4	16
BUSBO47012B10	4.7	12
BUSBO47014B10	4.7	14
BUSBO47016B10	4.7	16
BUSBO47018B10	4.7	18
BUSBO47020B10	4.7	20
BUSBO47022B10	4.7	22
BUSBO47024B10	4.7	24
BUSBO47026B10	4.7	26
BUSBO47028B10	4.7	28
BUSBO47030B10	4.7	30
BUSBO48012B10	4.8	12
BUSBO48014B10	4.8	14
BUSBO48016B10	4.8	16
BUSBO48018B10	4.8	18
BUSBO48020B10	4.8	20
BUSBO48022B10	4.8	22
BUSBO48024B10	4.8	24
BUSBO48026B10	4.8	26
BUSBO48028B10	4.8	28
BUSBO48030B10	4.8	30

**P FLEX URETERAL STENT SET ONE END CLOSED WITH STRING**

Product Code	Size (in Fr)	Length (in cm)
BUSB005012B10	5	12
BUSB005014B10	5	14
BUSB005016B10	5	16
BUSB005018B10	5	18
BUSB005020B10	5	20
BUSB005022B10	5	22
BUSB005024B10	5	24
BUSB005026B10	5	26
BUSB005028B10	5	28
BUSB005030B10	5	30
BUSB006012B10	6	12
BUSB006014B10	6	14
BUSB006016B10	6	16
BUSB006018B10	6	18
BUSB006020B10	6	20
BUSB006022B10	6	22
BUSB006024B10	6	24
BUSB006026B10	6	26
BUSB006028B10	6	28
BUSB006030B10	6	30
BUSB007012B10	7	12
BUSB007014B10	7	14
BUSB007016B10	7	16
BUSB007018B10	7	18
BUSB007020B10	7	20
BUSB007022B10	7	22
BUSB007024B10	7	24
BUSB007026B10	7	26
BUSB007028B10	7	28
BUSB007030B10	7	30

**P FLEX URETERAL STENT SET ONE END CLOSED WITH STRING**

Product Code	Size (in Fr)	Length (in cm)
BUSB008012B10	8	12
BUSB008014B10	8	14
BUSB008016B10	8	16
BUSB008018B10	8	18
BUSB008020B10	8	20
BUSB008022B10	8	22
BUSB008024B10	8	24
BUSB008026B10	8	26
BUSB008028B10	8	28
BUSB008030B10	8	30



**P FLEX URETERAL STENT SET BOTH ENDS OPEN WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B11	3	10
BUSB003012B11	3	12
BUSB003014B11	3	14
BUSB003016B11	3	16
BUSB004010B11	4	10
BUSB004012B11	4	12
BUSB004014B11	4	14
BUSB004016B11	4	16
BUSB047012B11	4.7	12
BUSB047014B11	4.7	14
BUSB047016B11	4.7	16
BUSB047018B11	4.7	18
BUSB047020B11	4.7	20
BUSB047022B11	4.7	22
BUSB047024B11	4.7	24
BUSB047026B11	4.7	26
BUSB047028B11	4.7	28
BUSB047030B11	4.7	30
BUSB048012B11	4.8	12
BUSB048014B11	4.8	14
BUSB048016B11	4.8	16
BUSB048018B11	4.8	18
BUSB048020B11	4.8	20
BUSB048022B11	4.8	22
BUSB048024B11	4.8	24
BUSB048026B11	4.8	26
BUSB048028B11	4.8	28
BUSB048030B11	4.8	30

**P FLEX URETERAL STENT SET BOTH ENDS OPEN WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB005012B11	5	12
BUSB005014B11	5	14
BUSB005016B11	5	16
BUSB005018B11	5	18
BUSB005020B11	5	20
BUSB005022B11	5	22
BUSB005024B11	5	24
BUSB005026B11	5	26
BUSB005028B11	5	28
BUSB005030B11	5	30
BUSB006012B11	6	12
BUSB006014B11	6	14
BUSB006016B11	6	16
BUSB006018B11	6	18
BUSB006020B11	6	20
BUSB006022B11	6	22
BUSB006024B11	6	24
BUSB006026B11	6	26
BUSB006028B11	6	28
BUSB006030B11	6	30
BUSB007012B11	7	12
BUSB007014B11	7	14
BUSB007016B11	7	16
BUSB007018B11	7	18
BUSB007020B11	7	20
BUSB007022B11	7	22
BUSB007024B11	7	24
BUSB007026B11	7	26
BUSB007028B11	7	28
BUSB007030B11	7	30

**P FLEX URETERAL STENT SET BOTH ENDS OPEN WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB008012B11	8	12
BUSB008014B11	8	14
BUSB008016B11	8	16
BUSB008018B11	8	18
BUSB008020B11	8	20
BUSB008022B11	8	22
BUSB008024B11	8	24
BUSB008026B11	8	26
BUSB008028B11	8	28
BUSB008030B11	8	30

**P FLEX URETERAL STENT SET ONE END CLOSED WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B12	3	10
BUSB003012B12	3	12
BUSB003014B12	3	14
BUSB003016B12	3	16
BUSB004010B12	4	10
BUSB004012B12	4	12
BUSB004014B12	4	14
BUSB004016B12	4	16
BUSB047012B12	4.7	12
BUSB047014B12	4.7	14
BUSB047016B12	4.7	16
BUSB047018B12	4.7	18
BUSB047020B12	4.7	20
BUSB047022B12	4.7	22
BUSB047024B12	4.7	24
BUSB047026B12	4.7	26
BUSB047028B12	4.7	28
BUSB047030B12	4.7	30
BUSB048012B12	4.8	12
BUSB048014B12	4.8	14
BUSB048016B12	4.8	16
BUSB048018B12	4.8	18
BUSB048020B12	4.8	20
BUSB048022B12	4.8	22
BUSB048024B12	4.8	24
BUSB048026B12	4.8	26
BUSB048028B12	4.8	28
BUSB048030B12	4.8	30

**P FLEX URETERAL STENT SET ONE END CLOSED WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB005012B12	5	12
BUSB005014B12	5	14
BUSB005016B12	5	16
BUSB005018B12	5	18
BUSB005020B12	5	20
BUSB005022B12	5	22
BUSB005024B12	5	24
BUSB005026B12	5	26
BUSB005028B12	5	28
BUSB005030B12	5	30
BUSB006012B12	6	12
BUSB006014B12	6	14
BUSB006016B12	6	16
BUSB006018B12	6	18
BUSB006020B12	6	20
BUSB006022B12	6	22
BUSB006024B12	6	24
BUSB006026B12	6	26
BUSB006028B12	6	28
BUSB006030B12	6	30
BUSB007012B12	7	12
BUSB007014B12	7	14
BUSB007016B12	7	16
BUSB007018B12	7	18
BUSB007020B12	7	20
BUSB007022B12	7	22
BUSB007024B12	7	24
BUSB007026B12	7	26
BUSB007028B12	7	28
BUSB007030B12	7	30

**P FLEX URETERAL STENT SET ONE END CLOSED WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB008012B12	8	12
BUSB008014B12	8	14
BUSB008016B12	8	16
BUSB008018B12	8	18
BUSB008020B12	8	20
BUSB008022B12	8	22
BUSB008024B12	8	24
BUSB008026B12	8	26
BUSB008028B12	8	28
BUSB008030B12	8	30

**P FLEX URETERAL STENT SET BOTH ENDS OPEN BOTH ENDS MULTILOOP WITH STRING**

Product Code	Size (in Fr)	Length (in cm)
BUSBO47018B13	4.7	18
BUSBO47024B13	4.7	24
BUSBO48018B13	4.8	18
BUSBO48024B13	4.8	24
BUSBO05018B13	5	18
BUSBO05024B13	5	24
BUSBO06018B13	6	18
BUSBO06024B13	6	24
BUSBO07018B13	7	18
BUSBO07024B13	7	24
BUSBO08018B13	8	18
BUSBO08024B13	8	24

**P FLEX URETERAL STENT SET ONE END CLOSED BOTH ENDS MULTILOOP WITH STRING**

Product Code	Size (in Fr)	Length (in cm)
BUSBO47018B14	4.7	18
BUSBO47024B14	4.7	24
BUSBO48018B14	4.8	18
BUSBO48024B14	4.8	24
BUSBO05018B14	5	18
BUSBO05024B14	5	24
BUSBO06018B14	6	18
BUSBO06024B14	6	24
BUSBO07018B14	7	18
BUSBO07024B14	7	24
BUSBO08018B14	8	18
BUSBO08024B14	8	24

**P FLEX URETERAL STENT SET BOTH ENDS OPEN BEML WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSBO47018B15	4.7	18
BUSBO47024B15	4.7	24
BUSBO48018B15	4.8	18
BUSBO48024B15	4.8	24
BUSBO05018B15	5	18
BUSBO05024B15	5	24
BUSBO06018B15	6	18
BUSBO06024B15	6	24
BUSBO07018B15	7	18
BUSBO07024B15	7	24
BUSBO08018B15	8	18
BUSBO08024B15	8	24

**P FLEX URETERAL STENT SET ONE END CLOSED BEML WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSBO47018B16	4.7	18
BUSBO47024B16	4.7	24
BUSBO48018B16	4.8	18
BUSBO48024B16	4.8	24
BUSBO05018B16	5	18
BUSBO05024B16	5	24
BUSBO06018B16	6	18
BUSBO06024B16	6	24
BUSBO07018B16	7	18
BUSBO07024B16	7	24
BUSBO08018B16	8	18
BUSBO08024B16	8	24



[www.bluneem.com](http://www.bluneem.com)

## **Blue Neem Medical Devices Pvt Ltd**

Plot Nos. 270 & 271, Road No. 5,  
Harohalli Industrial Area - Phase II,  
Kanakapura Taluk, Ramanagara District -562112  
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855  
E: [contact@bluneem.com](mailto:contact@bluneem.com), [marketing@bluneem.com](mailto:marketing@bluneem.com)

BNMD/MKG/LFT/17 | Rev.01 | 12-Dec-2020

## **Marketed & Distributed By**

## **Bioactiva Diagnostica GmbH**

Louisen Str. 137, Postal code: 61348  
Bad Homburg, Germany

Tel.: +49 6172 17102-20  
Email: [bioactiva@bioactiva.de](mailto:bioactiva@bioactiva.de)