

ENDOPYELOTOMY STENTS

Used for internal drainage from uretero pelvic junction to the bladder following incision of a stricture.



Key Features

- Dual Diameter Design
- Smooth transition from larger diameter to smaller diameter
- Maintain ureteral patency post Endopyelotomy

Also Available In Set Which Includes

Stent, Pusher, Guide Wire & Clamp

Size (in Fr)	Length (in cm)
6/12, 7/14	24 & 26

ENDOPYELOTOMY STENT

Product Code	Size (in Fr)	Length (in cm)
BEPB612024B1	6/12	24
BEPB612026B1	6/12	26
BEPB714024B1	7/14	24
BEPB714026B1	7/14	26

ENDOPYELOTOMY STENT SET

Product Code	Size (in Fr)	Length (in cm)
BEPB612024B2	6/12	24
BEPB612026B2	6/12	26
BEPB714024B2	7/14	24
BEPB714026B2	7/14	26

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/19 | Rev.01 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de