



URETERAL CATHETER

- Used for catheterization of the urinary tract during intraoperative procedures
- Used for injection of contrast media



Key Features

- Radiopaque for better visualization
- 1 cm marking helps in catheter positioning and measuring the length of ureter
- Detachable hub for scope insertion

Size (in Fr)

Length (in cm)

3, 4, 5, 6 & 7

70

Types

- Open End
- Closed End
- Cone Tip
- Angled Tip

URETERAL CATHETER STRAIGHT OPEN END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B1	3	70
BUCP004070B1	4	70
BUCP005070B1	5	70
BUCP006070B1	6	70
BUCP007070B1	7	70

URETERAL CATHETER STRAIGHT CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B2	3	70
BUCP004070B2	4	70
BUCP005070B2	5	70
BUCP006070B2	6	70
BUCP007070B2	7	70

URETERAL CATHETER CONE TIP OPEN END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B3	3	70
BUCP004070B3	4	70
BUCP005070B3	5	70
BUCP006070B3	6	70
BUCP007070B3	7	70

URETERAL CATHETER CONE TIP CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B4	3	70
BUCP004070B4	4	70
BUCP005070B4	5	70
BUCP006070B4	6	70
BUCP007070B4	7	70

URETERAL CATHETER ANGLE TIP OPEN END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B5	3	70
BUCP004070B5	4	70
BUCP005070B5	5	70
BUCP006070B5	6	70
BUCP007070B5	7	70

URETERAL CATHETER ANGLE TIP CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B6	3	70
BUCP004070B6	4	70
BUCP005070B6	5	70
BUCP006070B6	6	70
BUCP007070B6	7	70

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/01 | Rev.02 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de