

MONO J CATHETER

Used for temporary drainage and irrigation of urine from the kidney (Pigtail Ureteral Catheter).



Key Features

- Pigtail loop helps in retention of distal tip in the renal pelvis
- Holes in shaft helps in better drainage
- Used to inject contrast
- Used for short term external drainage to avoid stenting

Types

Open End & Closed End

Also available in Set, which Includes

Mono J Catheter, Guide Wire & Clamp

Size (in Fr)	Length (in cm)
3, 4, 5, 6 & 7	70

MONO J CATHETER OPEN END

Product Code	Size (in Fr)	Length (in cm)
BMJB003070B1	3	70
BMJB004070B1	4	70
BMJB005070B1	5	70
BMJB006070B1	6	70
BMJB007070B1	7	70

MONO J CATHETER CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BMJB003070B2	3	70
BMJB004070B2	4	70
BMJB005070B2	5	70
BMJB006070B2	6	70
BMJB007070B2	7	70

MONO J CATHETER SET OPEN END

Product Code	Size (in Fr)	Length (in cm)
BMJB003070B3	3	70
BMJB004070B3	4	70
BMJB005070B3	5	70
BMJB006070B3	6	70
BMJB007070B3	7	70

MONO J CATHETER SET CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BMJB003070B4	3	70
BMJB004070B4	4	70
BMJB005070B4	5	70
BMJB006070B4	6	70
BMJB007070B4	7	70

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/20 | Rev.01 | 12-Dec-2020

Marketed & Distributed By**Bioactiva Diagnostica GmbH**

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de