

HYDRO WIRE (BLACK COLOUR)

Used to gain access to desired anatomical location and to assist in placement of medical devices.



Key Features

- It is a kink resistant Nitinol core hydrophilic coated wire.
- The lubricious coating makes the wire slippery and helps in easy passage in tortuous anatomy.
- Radiopaque Polymer jacket for better fluoroscopic visibility.

Types

Straight Tip & Angled Tip

Size (in Inch)	Length (in cm)
0.018, 0.025, 0.032, 0.035 & 0.038	150

HYDRO WIRE - ANGLED TIP		
Product Code	Size (in Fr)	Length (in cm)
BHWN018150BAT	0.018	150
BHWN025150BAT	0.025	150
BHWN032150BAT	0.032	150
BHWN035150BAT	0.035	150
BHWN038150BAT	0.038	150

HYDRO WIRE - STRAIGHT TIP		
Product Code	Size (in Fr)	Length (in cm)
BHWN018150BST	0.018	150
BHWN025150BST	0.025	150
BHWN032150BST	0.032	150
BHWN035150BST	0.035	150
BHWN038150BST	0.038	150

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/23 | Rev.01 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de