

AMPLATZ PTFE SUPER STIFF GUIDEWIRE

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures. (PCNL)



Key Features

- PTFE coated stainless steel mandrel with 5 cm atraumatic floppy tip.
- Smooth transition from stiff mandrel to the flexible tip helps in easy navigation.
- The stiffness of the mandrel design provides support for the advancement of devices such as dilators, catheters and instruments.

Types

J Tip

Size (in Fr)

0.035

Length (in cm)

80

AMPLATZ PTFE SUPER STIFF GUIDEWIRE

Product Code

BASS035080BJT

Size (in Fr)

0.035

Length (in cm)

80

www.blueneem.com**Blue Neem Medical Devices Pvt Ltd**

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@blueneem.com, marketing@blueneem.com

BNMD/MKG/LFT/25 | Rev.01 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de