

PCN CATHETER

PCN Catheter is intended to be inserted percutaneously into the renal pelvis and to gain access to the upper urinary tract to perform external drainage of urine from the renal pelvis.



Key Features

- The pigtail tip configuration helps in retention of the distal tip in the renal pelvis
- Eye holes for better drainage
- Hydrophilic coating for easy insertion

Package Includes

PCN Catheter with Metal & Plastic Straightner, Urine Bag Connector

Size (in Fr)	Length (in cm)
6, 7, 8, 9, 10, 12, 14 & 16	30

PCN CATHETER

Product Code	Size (in Fr)	Length (in cm)
BPCP006030B1	6	30
BPCP007030B1	7	30
BPCP008030B1	8	30
BPCP009030B1	9	30
BPCP010030B1	10	30
BPCP012030B1	12	30
BPCP014030B1	14	30
BPCP016030B1	16	30

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/O2 | Rev.02 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de