



Chirana® 3-part disposable Catheter syringes

Chirana 3-part disposable Catheter syringes are made of polypropylene, polyethylene and lubricated synthetic rubber. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The syringes are intended for irrigation application and available in Catheter tip and with Luer adapter for wider application range.

USER CONVENIENCE: Ergonomic design, smooth flow, highly transparent barrel, tight, black extended graduation and easy to read.

COMFORT: Designed tip for best connection and scratch-free for glove manipulation.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

Specification	Pcs./Box	Pcs./Carton	Article No.
50 ml Catheter		100	CH03050C
50 ml Catheter with Luer adapter		100	CH03050CLA
100 ml Catheter with Luer adapter	25	100	CH03100C
150 ml Catheter Janette®		50	CH03150CS *

* under declared conditions suitable for autoclavable process

Shelf life 5 years

According to EN, ISO standards. Quality management certified according to EN ISO 13485:2016. CE marked.