

NEPHROSTOMY BALLOON DILATOR

Used for dilation of Nephrostomy tract.



Product Specifications

- Shaft Size: 5 Fr
- Balloon Diameter: 10 mm
- Balloon Length: 12 cm
- Catheter Length: 55 cm
- RBP: 11 atm
- Use Guide Wire: 0.035 inch

Product Features

- Controlled inflation pressure allows dilating force at a nephrostomy tract
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy

NEPHROSTOMY BALLOON DILATOR

Product Code	Size (in Fr)	Length (in cm)
BNB0005055BO	5	55

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/32 | Rev.01 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de