

SCREW DILATOR

Used for fascial dilation in one step.



Product Specifications

- Radiopaque for better visualization
- Screw Tip for one step dilation of fascia during Mini PCNL procedure
- Also available with compatible Amplatz sheath.
- Avoids sequential dilatation, reduces procedural time for practitioner

Size (in Fr)	Length (in cm)
12, 14 & 16	25

SCREW DILATOR		
Product Code	Size (in Fr)	Length (in cm)
BOS0012025BO	12	25
BOS0014025BO	14	25
BOS0016025BO	16	25

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/33 | Rev.01 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de