

AMPLATZ RENAL DILATOR

Used for dilation of a percutaneous tract and setup of a working channel into the kidney for PCNL procedures by placing a suitable amplatz sheath.



Key Features

- Gradual taper designed for atraumatic placement and maneuverability
- Smooth surface for ease of introduction
- Radiopaque for better visualization
- Dilator with suitable amplatz sheath

Dilator Size (OD of Dilator in Fr)	Length (in cm)
12, 14, 16, 18, 20, 22, 24, 26, 28 & 30	30

Amplatz Sheath Size (ID of the Sheath in Fr)	Length (in cm)
12, 14, 16, 18, 20, 22, 24, 26, 28 & 30	17

AMPLATZ RENAL DILATOR		
Product Code	Size (in Fr)	Length (in cm)
BRD0012030BS	12	30
BRD0014030BS	14	30
BRD0016030BS	16	30
BRD0018030BS	18	30
BRD0020030BS	20	30
BRD0022030BS	22	30
BRD0024030BS	24	30
BRD0026030BS	26	30
BRD0028030BS	28	30
BRD0030030BS	30	30

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/26 | Rev.01 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de