


N QUADRA HELICAL Z - NITINOL HELICAL BASKET WITHOUT TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.



Key Features

- Zero tip design reduces chances of perforation and mucosal damage, which also avoids stone migration during the procedure
- Nitinol wire construction which has high strength and high memory helps the basket to retain its shape during extreme torsion
- Helical design helps in capturing small stones and fragments effectively

Size (in Fr)

3 & 4

Length (in cm)

70, 90 & 120

N QUADRA HELICAL Z - (NITINOL HELICAL BASKET WITHOUT TIP)

Product Code	Size (in Fr)	Length (in cm)
BSBN003070BHO	3	70
BSBN003090BHO	3	90
BSBN003120BHO	3	120
BSBN004070BHO	4	70
BSBN004090BHO	4	90
BSBN004120BHO	4	120

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/37 | Rev.01 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de