


## N QUADRA HELICAL - NITINOL HELICAL WITH TIP BASKET

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.



### Key Features

- Nitinol wire construction which has high strength and high memory helps the basket to retain its shape during extreme torsion
- Helical design helps in capturing small stones and fragments effectively

### Size (in Fr)

3 & 4

### Length (in cm)

70, 90 & 120

## N QUADRA HELICAL (NITINOL STONE BASKET HELICAL WITH TIP)

Product Code	Size (in Fr)	Length (in cm)
BSBN003070BH1	3	70
BSBN003090BH1	3	90
BSBN003120BH1	3	120
BSBN004070BH1	4	70
BSBN004090BH1	4	90
BSBN004120BH1	4	120

[www.bluneem.com](http://www.bluneem.com)

### Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,  
Harohalli Industrial Area - Phase II,  
Kanakapura Taluk, Ramanagara District -562112  
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855  
E: [contact@bluneem.com](mailto:contact@bluneem.com), [marketing@bluneem.com](mailto:marketing@bluneem.com)

BNMD/MKG/LFT/38| Rev.01 | 12-Dec-2020

### Marketed & Distributed By

### Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348  
Bad Homburg, Germany

Tel.: +49 6172 17102-20  
Email: [bioactiva@bioactiva.de](mailto:bioactiva@bioactiva.de)