

DELTA BASKET

Steerable Basket for rigid nephroscope for stone extraction from calices which are inaccessible due to acute angle.



Design Features

- Can be angled inside kidney to reach inner calyx
- Can be used for straight path stone manipulation and also can be angled both directions
- Unique basket shape which prevents stone migration

Procedural Features

- Prevents injury to the infundibulum and reduce bleeding
- Avoids the need for multiple tracts
- Special advantage to centres which do not have flexible Nephroscopes
- Helps to remove big blood clots

| Size (in Fr) | Length (in cm) |
|--------------|----------------|
| 7 | 40 |

| DELTA BASKET | | |
|--------------|--------------|----------------|
| Product Code | Size (in Fr) | Length (in cm) |
| BDBRO07040BO | 7 | 40 |

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/41 | Rev.01 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de