

BIOPSY GUN - SEMI AUTOMATIC

Used for obtaining histological core samples from soft tissue and organs



Key Features

- Echogenic and atraumatic needle tip
- Ergonomic and light weight
- Cannula with depth marking
- Easy and one hand operated with two stage arming
- Coaxial available with trocar and blunt needle

Semi Automatic - with Coaxial

- 14, 16 G x 10, 15 cm - For breast biopsy
- 16, 18 G x 15, 20 cm - For liver & kidney
- 18, 20 G x 15, 20 cm - For lung

SEMI AUTOMATIC BIOPSY GUN WITH COAXIAL

Product Code	Size (in G)	Length (in cm)
BBG0014010BS1	14	10
BBG0014015BS1	14	15
BBG0016010BS1	16	10
BBG0016015BS1	16	15
BBG0016020BS1	16	20
BBG0018015BS1	18	15
BBG0018020BS1	18	20
BBG0020015BS1	20	15
BBG0020020BS1	20	20

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/42 | Rev.01 | 12-Dec-2020

Marketed & Distributed By
Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de