


BIOPSY GUN - FULLY AUTOMATIC

Used for obtaining histological core samples from soft tissue and organs



Key Features

- Echogenic and atraumatic needle tip
- Ergonomic and light weight
- Cannula with depth marking
- Easy and one hand operated with two stage arming
- Coaxial available with trocar and blunt needle

Fully Automatic - with Coaxial

- 18 G x 20, 25 & 30 cm - For prostate
- 16, 18 G x 15, 20 cm - For liver & kidney
- 14, 16 G x 10, 15 cm - For lung

FULLY AUTOMATIC BIOPSY GUN WITH COAXIAL

Product Code	Size (in G)	Length (in cm)
BBG0014010BF1	14	10
BBG0014015BF1	14	15
BBG0016010BF1	16	10
BBG0016015BF1	16	15
BBG0016020BF1	16	20
BBG0018015BF1	18	15
BBG0018020BF1	18	20
BBG0018025BF1	18	25
BBG0018030BF1	18	30

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/43 | Rev.01 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de