


DURAFLEX T™ STENT (TUMOR STENT)

Used for temporary stenting of ureter in patients with extrinsic ureteral obstructions.



Key Features

- Manufactured by using silicone blend polyurethane for enhanced patient comfort
- Reinforced stainless steel coil layer provides superior compressive strength
- Stiff shaft wall enables stent to be placed in tight strictures
- Compatible with 0.035" guide wire

Types	Size (in Fr)	Length (in cm)
Both Ends Open	6 & 7	24 & 26

Also available in Set, which Includes

Stent with or without string, pusher, guidewire & clamp

DURAFLEX T STENT BOTH ENDS OPEN

Product Code	Size (in Fr)	Length (in cm)
BTSS006024B1	6	24
BTSS006026B1	6	26
BTSS007024B1	7	24
BTSS007026B1	7	26

DURAFLEX T STENT BOTH ENDS OPEN WITH STRING

Product Code	Size (in Fr)	Length (in cm)
BTSS006024B2	6	24
BTSS006026B2	6	26
BTSS007024B2	7	24
BTSS007026B2	7	26

DURAFLEX T STENT SET BOTH ENDS OPEN

Product Code	Size (in Fr)	Length (in cm)
BTSS006024B3	6	24
BTSS006026B3	6	26
BTSS007024B3	7	24
BTSS007026B3	7	26

DURAFLEX T STENT SET BOTH ENDS OPEN WITH STRING

Product Code	Size (in Fr)	Length (in cm)
BTSS006024B4	6	24
BTSS006026B4	6	26
BTSS007024B4	7	24
BTSS007026B4	7	26

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/52

Rev.01

15-Jun-2021

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20

Email: bioactiva@bioactiva.de