


# DUAL LUMEN URETERAL CATHETER

Used for

- Catheterization of the urinary tract during intraoperative procedures
- Injection of contrast media
- Access, advancement or exchange of ancillary devices including guidewire and drainage of fluids



## Key Features

- Has double lumen, one lumen is for insertion over the guidewire and another lumen is for injection of contrast medium
- Optimal radiopacity for better visualization
- With 1 cm marking which helps in catheter position, measurement of ureter length
- Exchange of wires without removing catheter

Size (in Fr)

8, 10 & 12

Length (in cm)

50 & 70

## Advantages

- Can inject contrast while guidewire in position
- Reduced procedural time

## DUAL LUMEN URETERAL CATHETER

Product Code	Size (in Fr)	Length (in cm)
BUCP008050B19	8	50
BUCP010050B19	10	50
BUCP012050B19	12	50
BUCP008070B19	8	70
BUCP010070B19	10	70
BUCP012070B19	12	70

[www.bluneem.com](http://www.bluneem.com)

### Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,  
Harohalli Industrial Area - Phase II,  
Kanakapura Taluk, Ramanagara District -562112  
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855  
E: [contact@bluneem.com](mailto:contact@bluneem.com), [marketing@bluneem.com](mailto:marketing@bluneem.com)

BNMD/MKG/LFT/53

Rev.01

15-Jun-2021

### Marketed & Distributed By

### Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348  
Bad Homburg, Germany

Tel.: +49 6172 17102-20

Email: [bioactiva@bioactiva.de](mailto:bioactiva@bioactiva.de)