

# PCN SET WITH LOCKING MECHANISM (PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL" TYPE)

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage. The locking mechanism prevents accidental slippage of the catheter.



## Set Includes

IP Needle | Guidewire J Tip PTFE Coated | Fascial Dilators | Urine Bag Connector  
PCN Catheter with Locking Mechanism & Hydrophilic Coating | Scalpel Blade  
Catheter Fixing Disc

## Advantages

Kit provides entire set of devices to perform Nephrostomy procedure

Size (in Fr)	Length (in cm)
6, 7, 8, 10, 12, 14 & 16	35

**PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL TYPE" WITH LOCKING MECHANISM**

Product Code	Size (in Fr)	Length (in cm)
BPCP006035B4	6	35
BPCP007035B4	7	35
BPCP008035B4	8	35
BPCP010035B4	10	35
BPCP012035B4	12	35
BPCP014035B4	14	35
BPCP016035B4	16	35

[www.bluneem.com](http://www.bluneem.com)

**Blue Neem Medical Devices Pvt Ltd**

Plot Nos. 270 & 271, Road No. 5,  
Harohalli Industrial Area - Phase II,  
Kanakapura Taluk, Ramanagara District -562112  
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855  
E: [contact@bluneem.com](mailto:contact@bluneem.com), [marketing@bluneem.com](mailto:marketing@bluneem.com)

BNMD/MKG/LFT/06 | Rev.02 | 12-Dec-2020

**Marketed & Distributed By**

**Bioactiva Diagnostica GmbH**

Louisen Str. 137, Postal code: 61348  
Bad Homburg, Germany

Tel.: +49 6172 17102-20  
Email: [bioactiva@bioactiva.de](mailto:bioactiva@bioactiva.de)