

NEPHROSTOMY DRAINAGE SET DIRECT PUNCTURE

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage by direct puncture.



Key Features

- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Eye holes for better drainage
- Hydrophilic coating for easy insertion
- Trocar design and catheter design prevents buckling of catheter during the procedure

Package Includes

PCN Catheter with Metal Straightner, Trocar & Urine Bag Connector

Size (in Fr)

7, 8, 9, 10, 12, 14, 16 & 18

Length (in cm)

30

NEPHROSTOMY DRAINAGE SET DIRECT PUNCTURE

Product Code	Size (in Fr)	Length (in cm)
BPCP007030B5	7	30
BPCP008030B5	8	30
BPCP009030B5	9	30
BPCP010030B5	10	30
BPCP012030B5	12	30
BPCP014030B5	14	30
BPCP016030B5	16	30
BPCP018030B5	18	30

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/07 | Rev.02 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de