

# NEPHROSTOMY DRAINAGE SET DIRECT PUNCTURE WITH LOCKING MECHANISM

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage by direct puncture. The locking mechanism prevents accidental slippage of the catheter.



## Key Features

- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Eye holes for better drainage
- Hydrophilic coating for easy insertion
- Trocar and catheter design prevents buckling of catheter during the procedure
- The locking mechanism prevents accidental slippage of the catheter

### Package Includes

PCN Catheter with Locking Mechanism & Metal Straightner, Trocar & Urine bag Connector

Size (in Fr)	Length (in cm)
7, 8, 10, 12, 14 & 16	35

### NEPHROSTOMY DRAINAGE SET DIRECT PUNCTURE WITH LOCKING MECHANISM

Product Code	Size (in Fr)	Length (in cm)
BPCP007035B6	7	35
BPCP008035B6	8	35
BPCP010035B6	10	35
BPCP012035B6	12	35
BPCP014035B6	14	35
BPCP016035B6	16	35

[www.bluneem.com](http://www.bluneem.com)

#### Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,  
Harohalli Industrial Area - Phase II,  
Kanakapura Taluk, Ramanagara District -562112  
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855  
E: [contact@bluneem.com](mailto:contact@bluneem.com), [marketing@bluneem.com](mailto:marketing@bluneem.com)

BNMD/MKG/LFT/08 | Rev.02 | 12-Dec-2020

#### Marketed & Distributed By

#### Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348  
Bad Homburg, Germany

Tel.: +49 6172 17102-20  
Email: [bioactiva@bioactiva.de](mailto:bioactiva@bioactiva.de)