



Vitrex End-to-End Pipettes

Vitrex End-to-End Pipettes are intended for easy and accurate capillary tube volume measurement. The pipettes are in conformity with ISO 7550.

The pipettes are packed in glass vials with a dispenser top, that makes it easy to shake out one tube at a time and thereby ensuring hygienic handling.

The dispenser tops are colour-coded to indicate the tube volume in conformity with ISO 1769. Ten vials are packaged in a sturdy carton in order to protect the tubes from damage during handling or transport.

Vitrex tolerances

Vol. (µl)	Accuracy	CV
0.5 - 10	± 0.50%	0.75%
20 - 100	± 0.50%	0.60%





Technical Specifications

Colour Code	Vol. (µl)	L (mm) +/- 1mm	Pcs.	Plain	Natrium Hep.	Ammon. Hep.	K2EDTA
○	0.5	15	10x100	177613	*	*	*
●	1	29	10x100	165013	167513	169013	*
●	2	29	10x100	165113	167613	169113	*
●	3	29	10x100	165213	167713	169213	*
●	4	29	10x100	165313	167813	169313	*
○	5	29	10x100	165413	167913	169413	170113
●	10	29	10x100	174313	177113	179213	174113
●	20	29	10x100	174513	177213	179313	174413
○	25	29	10x100	174613	177313	179413	178213
●	40	29	10x100	174813	177413	179513	*
●	44.7	29	10x50	175013	177713	178613	*
●	50	29	10x50	174913	177513	179613	171113

* Available upon request
For lithium hep. please contact us

Shelf life 4 years



Accessories

Product	Article No.
Pipette holder	140413

