

# AMPLATZ SHEATH

Used as a working channel to maintain a previously established nephrostomy tract for introduction of instruments and catheters during urological PCNL procedures.



## Key Features

- Bevel tip for easy insertion.
- Smooth surface for ease of introduction
- Radiopaque for better visualization

## Size (in Fr)

12, 14, 16, 18, 20, 22, 24, 26,  
28 & 30

## Length (in cm)

17 & 20

### AMPLATZ SHEATH - 17 cm

Product Code	Size (in Fr)	Length (in cm)
BAST012017B1	12	17
BAST014017B1	14	17
BAST016017B1	16	17
BAST018017B1	18	17
BAST020017B1	20	17
BAST022017B1	22	17
BAST024017B1	24	17
BAST026017B1	26	17
BAST028017B1	28	17
BAST030017B1	30	17

### AMPLATZ SHEATH - 20 cm

Product Code	Size (in Fr)	Length (in cm)
BAST012020B1	12	20
BAST014020B1	14	20
BAST016020B1	16	20
BAST018020B1	18	20
BAST020020B1	20	20
BAST022020B1	22	20
BAST024020B1	24	20
BAST026020B1	26	20
BAST028020B1	28	20
BAST030020B1	30	20

[www.bluneem.com](http://www.bluneem.com)

#### Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,  
Harohalli Industrial Area - Phase II,  
Kanakapura Taluk, Ramanagara District -562112  
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855  
E: [contact@bluneem.com](mailto:contact@bluneem.com), [marketing@bluneem.com](mailto:marketing@bluneem.com)

BNMD/MKG/LFT/16 | Rev.02 | 12-Dec-2020

#### Marketed & Distributed By

#### Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348  
Bad Homburg, Germany

Tel.: +49 6172 17102-20  
Email: [bioactiva@bioactiva.de](mailto:bioactiva@bioactiva.de)