

MALECOT CATHETER

Malecot Catheter is intended to be inserted percutaneously into the renal pelvis, and to gain access to the upper urinary tract to perform external drainage of urine from the renal pelvis.



Key Features

- The Malecot flower design is ideally suitable for the drainage of thick viscous fluids
- Integral flower is designed to match the catheter size. This reduces trauma during insertion

Package Includes

Malecot Catheter with Plastic or Metal Straightner & Urine Bag Connector

Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

Length (in cm)

30

MALECOT CATHETER

Product Code	Size (in Fr)	Length (in cm)
BMCPO08030B1	8	30
BMCPO10030B1	10	30
BMCPO12030B1	12	30
BMCPO14030B1	14	30
BMCPO16030B1	16	30
BMCPO18030B1	18	30
BMCPO20030B1	20	30
BMCPO22030B1	22	30

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/09 | Rev.02 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de