

PCN SET - PIGTAIL TYPE

(PERCUTANEOUS NEPHROSTOMY DRAINAGE SET)

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage.



Set Includes

IP Needle | Guidewire J Tip PTFE Coated | Fascial Dilators | Urine Bag Connector | PCN Catheter with Hydrophilic Coating | Scalpel Blade | Catheter Fixing Disc

Advantages

Kit provides entire set of devices to perform Nephrostomy procedure

Size (in Fr)	Length (in cm)
6, 7, 8, 9, 10, 12, 14, 16 & 18	30

PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL TYPE"

Product Code	Size (in Fr)	Length (in cm)
BPCP006030B3	6	30
BPCP007030B3	7	30
BPCP008030B3	8	30
BPCP009030B3	9	30
BPCP010030B3	10	30
BPCP012030B3	12	30
BPCP014030B3	14	30
BPCP016030B3	16	30
BPCP018030B3	18	30

www.bluneem.com

Blue Neem Medical Devices Pvt Ltd

Plot Nos. 270 & 271, Road No. 5,
Harohalli Industrial Area - Phase II,
Kanakapura Taluk, Ramanagara District -562112
Karnataka, India

T: +91 80 2976 1335/36, M: +91 97399 72854, 97399 72855
E: contact@bluneem.com, marketing@bluneem.com

BNMD/MKG/LFT/04 | Rev.02 | 12-Dec-2020

Marketed & Distributed By

Bioactiva Diagnostica GmbH

Louisen Str. 137, Postal code: 61348
Bad Homburg, Germany

Tel.: +49 6172 17102-20
Email: bioactiva@bioactiva.de